
# Varolast<sup>®</sup> Plus

Zinkleimbinde, längselastisch  
Zinc paste bandage, lengthways elastic  
Bande à l'oxyde de zinc, élastique en longueur  
Zinklijmzwachtel, elastisch in de lengte  
Benda all'ossido di zinco, elasticità longitudinale  
Ligadura de óxido de zinc, elasticidade longitudinal  
Επίδεσμος με κόλλα ψευδαργύρου, ελαστικός κατά μήκος  
Zinkoklihoné obinadlo, podléne elastický  
Zinkoglejový obvaz, pozdĺžne elastický  
Opaska z pastą cynkową, podłużnie elastyczna  
Cinkenyn-rólya, hosszában nyújtható  
Бинт с цинковой массой, продольно эластичный  
Бинт с цинкова паста, еластичен по дължина  
Zinklimbind, elastisk i længden  
縱向彈性鋅膏繃帶  
رباط بمعجون الزنك قابل للتمط بالطول

HARTMANN

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## EN Instructions for use

### General description of the device

Varolast Plus is a ready-to-use zinc paste bandage which is lengthways elastic. The bandage has a carrier material made of viscose and polyamide, and uses zinc paste as impregnation mass, which contains zinc oxide, glycerin and gelatin. Due to technical reasons, zinc oxide is added in order to achieve a semi-rigid consistency of the bandage; an ancillary pharmacological action is not intended. The bandage is white in colour. The single-use, non-active and non-sterile device is available in various widths and lengths.

### Intended purpose

Varolast Plus is a single-use, non-sterile zinc paste bandage that is used for semi-rigid support and as a compression dressing. It can be applied for long-term treatment in compression therapy of chronic venous diseases, thromboembolic venous diseases, and oedemata. This zinc paste bandage is indicated with an ABPI of 0.8-1.3 only. The product is to be applied on intact skin only. This zinc paste bandage is indicated for humans (not for children under 3 years of age) and can be applied by healthcare professionals in clinical and home environments.

### Composition

69% viscose, 31% zinc paste (packed in an aluminum foil pouch).

### Indications

The device can be used for:

- Chronic venous diseases, prevention and treatment of venous oedemata, prevention and treatment of venous skin changes, eczema and pigmentation, lipodermatosclerosis and atrophie blanche, treatment of venous leg ulcers, prevention of recurrent venous leg ulcers, varicose veins, initial phase after treatment of varicose veins, functional venous insufficiency (associated with obesity, jobs that involve sitting or standing), venous malformations.
- Thromboembolic venous diseases (superficial venous thrombosis, deep arm and leg venous thrombosis, post-thrombotic syndrome, thrombosis prophylaxis in mobile patients).
- Oedemata (lymphoedema, oedema in pregnancy, post-traumatic oedema, postoperative oedema, postoperative reperfusion oedema, idiopathic cyclic oedema, lipoeidema from stage II, congestive states as a result of immobility (arthrogenic congestive syndrome, paresis and partial paresis of extremities), drug-induced oedema if not possible to switch).

### Contraindications

The device must not be used in cases of advanced peripheral arterial occlusive disease (if one of the following parameters applies: ABPI < 0.5, ankle blood pressure < 60 mmHg, toe blood pressure < 30 mmHg or TcPO<sub>2</sub> < 20 mmHg instep), decompensated cardiac insufficiency, phlegmasia coerulea dolens, septic phlebitis, severe exuding dermatosis, intolerance to compression material or zinc paste ingredients, severe sensory disturbances of the extremities, advanced

peripheral neuropathy (e.g. associated with diabetes mellitus), primary chronic polyarthritits.

### Precautions

Special precautionary measures are required in the event of severe weeping dermatosis, severe sensory disturbances of the extremities, advanced peripheral neuropathy (e.g. associated with diabetes mellitus), primary chronic rheumatoid arthritis and compression material or zinc paste ingredients. The medical device should not be applied on children under 3 years of age.

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

### Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

### Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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PAUL HARTMANN AG