

Pütter®binde

Kräftige Kurzzugbinden
Strong short-stretch bandages
Bandes résistantes à allongement court
Sterke zwachtels met korte rek
Bende a corta trazione forte
Vendas resistentes de corta elasticidad
Ligaduras resistentes de curta tração
Ανεκτικοί επίδεσμοί μικρής τάσης
Mocne opaski o niskiej rozciągliwości
Erős, rövid megnyúlású pólyák
Прочные бинты короткой растяжимости
Bandaje rezistente cu elasticitate redusă
Starka kortsträcksbindor
強力短拉力繃帶
أربطة متينة قصيرة المدى

1

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EN Instructions for use

General description of the device

Pütterbinde* is a brown-coloured, especially strong short-stretch bandage with textile elasticity and an extensibility of approximately 85 % for very high compression.

Intended purpose

Pütterbinde can be used for long-term compression therapy:

- treatment of chronic venous insufficiency (CVI) with and without venous leg ulcers (VLU)
- thrombosis/thrombophlebitis
- chronic oedemata of other origins (ankle brachial pressure index – ABPI 0.8–1.3)
- lymph oedemata

The product can be used on humans by healthcare professionals in clinical and home environment in combination with padding bandages and inner lining material. Pütterbinde is reusable.

Properties

Extensibility: approx. 85%.
Washable: Wash the bandage in a regular wash cycle at a temperature of 95°C.

Indication

The product can be used for the treatment of varicose veins, thrombosis/thrombophlebitis, chronic venous insufficiency (CVI), primary and secondary prevention of leg ulcer venous/venous insufficiency and oedemata of other origins (ankle brachial pressure index - ABPI 0.8 - 1.3).

Intended users

Pütterbinde is designed to be used by healthcare professionals who have received appropriate training. Pütterbinde can be used for the treatment of: varicose veins, thrombosis/thrombophlebitis, post-thrombotic syndrome, thromboprophylaxis, chronic venous insufficiency (CVI), venous insufficiency, angiodysplasia and other oedemata.

Composition

100 % cotton

Contraindications

- Pütterbinde must not be used in case of:
- advanced peripheral arterial occlusive disease (PAOD)
 - decompensated heart failure (NYHA II +IV)
 - phlegmasia coerulea dolens

- in case of ankle-brachial pressure index (ABPI) > 1.3 and ≤ 0.5
- septic phlebitis

Precautions

Special precautionary measures are required in the event of:

- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral neuropathy)
- pronounced weeping dermatoses
- incompatibility with compression bandage material
- primary chronic polyarthritis

The application technique may have to be adapted to the patient's anatomical and pathophysiological characteristics (e.g. in regard to stretch, padding, etc.). If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. If an inelastic system is created by close to full stretch application, high ankle pressures should be closely monitored by a physician.

Application (see illustrations)

Modified bandaging technique according to Pütter, if not otherwise prescribed by the doctor. (1)

Further remarks

Pütterbinde is indicated for use by healthcare professionals who have received appropriate training. Before applying Pütterbinde, clean the wound and cover it with a suitable wound dressing. If required, excess padding bandage material may be used for additional padding of bony prominences and tissue recesses (e.g. shinbone, Bisgaard's backdrop). If additional fixation is required, use tape but no bandage clips. After application, the compression bandage must be checked (e.g. circulation in the patient's toes, pain reported by the patient) in order to be certain that the pressure of the bandage is not too high. If at any time during therapy with Pütterbinde the patient develops severe pain, numbness, tingling or noticeable discolouration of the toes, a healthcare professional must be contacted immediately and the bandage must be removed. During therapy with Pütterbinde, the patient should be encouraged to move.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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