

# Pütter-haft®

Kohäsive Kurzzugbinde  
Cohesive short-stretch bandage  
Bande cohésive à allongement court  
Cohesieve zwachtel met korte rek  
Benda coesiva a corta trazione  
Venda cohesiva de corta elasticidad  
Ligadura coesiva de curta tração  
Συνεκτικός επίδεσμος μικρής τάνσεως  
Kohezyjna opaska o niskiej rozciągliwości  
Kohéziós, rövid megnyúlású pólya  
Когезионный бинт короткой растяжимости  
Bandaj coeziv cu elasticitate redusă  
內聚短拉力繃帶  
رباط لاصق قصير المدى

1

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## EN Instructions for use

### General description of the device

Pütter-haft® is a brown-coloured, cohesive short-stretch bandage with textile elasticity and an extensibility of approximately 65 % for very high compression.

### Intended purpose

Pütter-haft can be used for long-term compression therapy:

- treatment of chronic venous insufficiency (CVI) with and without venous leg ulcers (VLU)
- thrombosis/thrombophlebitis
- chronic oedemata of other origins (ankle brachial pressure index – ABPI 0.8–1.3)
- lymph oedemata

The product can be used on humans by healthcare professionals in clinical and home environment in combination with padding bandages and inner lining material. If re-used, the functionality of the bandage will be reduced or lost.

### Properties

Extensibility: approx. 65%.

### Indication

The product can be used for the treatment of varicose veins, thrombosis / thrombophlebitis, chronic venous insufficiency (CVI), primary and secondary prevention of leg ulcer venous / venous insufficiency and oedemata of other origins (ankle brachial pressure index – ABPI 0.8 – 1.3).

### Intended users

Pütter-haft is designed to be used by healthcare professionals who have received appropriate training. Pütter-haft can be used for the treatment of: varicose veins, thrombosis / thrombophlebitis, post-thrombotic syndrome, thromboprophylaxis, chronic venous insufficiency (CVI), venous insufficiency, angiodysplasia and other oedemata.

### Composition

100 % cotton

#### Contraindications

Pütter-haft must not be used in case of:

- advanced peripheral arterial occlusive disease (PAOD)
- decompensated heart failure (NYHA II +IV)
- phlegmasia coerulea dolens
- in case of ankle-brachial pressure index (ABPI) > 1.3 and ≤ 0.5
- septic phlebitis

#### Precautions

Special precautionary measures are required in the event of:

- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral polyneuropathy)
- pronounced weeping dermatoses
- incompatibility with compression bandage material
- primary chronic polyarthritis

The application technique may have to be adapted to the patient's anatomical and pathophysiological characteristics (e.g. in regard to stretch, padding, etc.). If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. If an inelastic system is created by close to full stretch application, high ankle pressures should be closely monitored by a physician.

#### Application (see illustrations)

Modified bandaging technique according to Pütter, if not otherwise prescribed by the doctor. (1)

#### Further remarks

Pütter-haft is indicated for use by healthcare professionals who have received appropriate training. Before applying Pütter-haft, clean the wound and cover it with a suitable wound dressing. If required, excess padding bandage material may be used for additional padding of bony prominences and tissue recesses (e.g. shinbone, Bisgaard's backdrop). If additional fixation is required, use tape but no bandage clips. After application, the compression bandage must be checked (e.g. circulation in the patient's toes, pain reported by the patient) in order to be certain that the pressure of the bandage is not too high. If at any time during therapy with Pütter-haft the patient develops severe pain, numbness, tingling or noticeable discolouration of the toes, a healthcare professional must be contacted immediately and the bandage must be removed. During therapy with Pütter-haft, the patient should be encouraged to move.

#### Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

#### Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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