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Produktbeschreibung

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Instructions for use

Product description
Zetivut Plus Silicone Border is a sterile self-adhesive superabsorbent dressing with a silicone interface for the treatment of moderately to heavily exuding, chronic and acute wounds.

The dressing consists of a silicone adhesive top layer, a green nonwoven (100% polypropylene) to ensure water resistance and a soft, white, hydrophilic nonwoven (viscose and polyamide) facing away from the wound.

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Composition

Zetivut Plus Silicone Border comprises a semi-permeable polyurethane backing film, a perforated silicone film towards the side facing the wound and an absorbent pad with hydrophilic cellulose fibers and superabsorbent polyacrylate to absorb and retain the exudate. The backing film is permeable to air but waterproof, which allows the patient to shower. The perforated silicone film facing the wound side facilitates the application of the dressing

and promotes non-adherence to the wound. The absorbent pad comprises cellulose and superabsorbent polyacrylate for the absorption and retention of exudate and is wrapped in a hydrophilic cellulose tissue. It is covered with a hydrophilic nonwoven (viscose and polyamide).

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Intended purpose

Zetivut Plus Silicone Border is a single use sterile self-adhesive superabsorbent dressing with a silicone interface for long-term treatment of injured skin, acute and chronic wounds, with moderate to high levels of exudate. It is used on adults only, by health professionals with adequate training in the use of medical devices. The dressing can be used in conjunction with secondary dressings.

Indications

Zetivut Plus Silicone Border is suitable for the treatment of wounds with moderate to high levels of exudate: acute wounds (burns, post-op, trauma wounds) and chronic wounds (venous ulcers, pressure ulcers, diabetic foot ulcers, tumours).

Mode of application
Clean the wound using a physiological solution or as instructed by a clinician. Choose a suitable dressing size for the wound, ensuring that the central pad exceeds the edges of the wound by approximately 1 to 2 cm. Do not cut the absorbent green pad. The transparent border can be cut in order to adapt the dressing e.g. to a concave form.

Before applying the dressing, after rinsing the wound, make sure that the surrounding skin is dry by using wicks or gauze. The dressing should be applied to a clean, dry wound.

Remove the two protective sheets. Place the adherent side of the dressing onto the wound, so that the green side faces away from the wound.

Press down on the adhesive edges firmly and avoid wrinkles in the border to ensure the dressing stays securely in place.

No additional fixation is required, if indicated, compression bandages can be applied to the dressing to secure it in place.

The dressing can remain on the wound for up to 5 days depending on the condition of the wound and the established clinical protocol.

The dressing must be changed if clinically indicated or when exudate reaches the back of the absorbent pad.

Contraindications

Do not use Zetivut Plus Silicone Border in case of hypersensitivity to any of its components.

Special precautions

Do not cut the absorbent green pad. The transparent border can be cut in order to adapt the dressing to the patient's anatomy. Do not use the dressing on patients with the following signs and symptoms:

- Before treating a wound with healing impairment, a clinician has to be consulted to establish a clinical diagnosis. Treatment with Zetivut Plus Silicone Border cannot replace a clinical diagnosis. The alimhodilla absorbente el tratamiento de heridas con signos de infección. The treated area must be controlled with appropriate treatment before this dressing can be used.
- In all cases, follow the established clinical protocol.
- The dressing is available in different sizes and is suitable for dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary, on these population groups this dressing should be used with caution and following the established clinical protocol.
- It is not recommended to stand on the dressing, e.g. in the case of wounds on the feet.

Contraindications

Use of this dressing which are packed in paper-paper pack packages. The sealed seam of this pad contains natural rubber latex, which may cause allergic reactions!

Resuing a single-use medical device is dangerous. Reprocessing devices in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Side effects

Excessive exudate, maceration, increased pain, adherence of the dressing to the wound, itching, burning, oedema, allergic ulceration, skin torn, blisters, erythema. If any of these side effects occur, consult a clinician.

Incident reporting

For a patient-related incident in the European Union and in countries with identical regulatory requirements (Regulation 2017/745/EU on Medical Devices), if during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Produkt disposal

The risk of potential infection hazards, or environmental pollution, of disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device using household waste.

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