Veroval®

Netzgerät

Deutsch

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AC/DC Adaptor for HARTMANN Veroval® Blood Pressure Monitors

Product description

The mains part plug converts the alternative voltage of 100-240 V to a safety extra-low voltage. The device is tested according to the medical standard IEC 60601-1 and is intended for use with compatible Veroval® blood pressure monitors.

Safety and precautionary measures

- Read the operating instructions before initial use.
- Observe the Veroval® blood pressure monitor instructions for use.
- Only connect the power supply unit to power sources that are suitable for this purpose. Ensure the polarity is correct.
- The power supply unit must never be exposed to rain, moisture or high heat.
- Only operate power supply unit indoors.
- Do not use the power supply unit when damaged.
- Do not open or modify the power supply unit. Repairs may only be carried out by authorised technical personnel.
- Do not cover the power supply unit. Ensure sufficient cooling. Do not expose to direct sunlight.
- Keep power supply unit and cables away from flammable items.
- Only clean with a damp cloth. Do not use any cleaning products. Do not sterilise device.

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these Instructions for Use, including domestic environments.
- It may only be possible to use the device to a limited extent in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly beside other devices or stacked on top of other devices, as this could result in faulty operation. However, if it is necessary to use the device in the manner described above, this device and the other devices must be monitored to ensure that they are working properly.
- The use of applied parts other than those specified or provided by the manufacturer of this device may result in increased electromagnetic emissions or a decrease in the device's electromagnetic immunity and lead to faulty operation.
- Failure to comply with the above can impair the performance characteristics of the device.

Initial use

The power supply unit is ready for use once it is plugged into a mains socket. The power supply unit must not be used outside of the specified input values as this may lead to the device being destroyed.

Technical data

Type:	LXCP12-006060BEH
Item number:	925 391
Input voltage:	100 – 240V AC, 50/60Hz
Max input current:	500mA
Output voltage:	6V DC, 600mA, only in combination with HARTMANN Veroval® blood pressure monitors
Max output current:	600mA
⚡⚡⚡	Polarity of the DC voltage connection
Output current protection:	Protection against short circuit
No-load losses:	< 0.075 watts
Approvals:	Medical standard IEC 60601-1 (in compliance with CISPR 11, IEC 61000-3-2, IEC 61000-3-3)
Conformity:	CE compliant as per European guidelines
Operating conditions:	Ambient temperature: 10 °C to 40 °C Relative humidity: 20 – 90 %, non-condensing
Storage / transport conditions:	Ambient temperature: -25 °C to +70 °C Relative humidity: 20 – 90 %, non-condensing
Dimensions L x W x H in mm:	62 x 39 x 26 (without mains plug)

Signs and symbols:

Caution

Follow instructions for use

Do not use outdoors (indoor use only)

Protective insulation/protection class 2

Symbol for the marking of electrical and electronic equipment

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