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EN	Instructions for use
Product description	RespoSorb Super is a sterile superabsorbent dressing used on superficial, moderately to severely exuding wounds. The absorbent core absorbs and binds exudate.
Composition	RespoSorb Super is a combined superabsorbent dressing which consists of four layers of different materials. On the wound side, the product features a soft, white, hydrophilic nonwoven (viscose and polyamide). The inner dressing core consists of soft cellulose fluff blended with liquid-retaining polyacrylate polymers. This absorbent core is enclosed in a thin, nonwoven fabric that evenly distributes the liquid. On the side facing away from the wound, the product features a green layer of polypropylene nonwoven, which is water-repellent but permeable to air and allows for gas exchange.
Properties and mode of action	RespoSorb Super rapidly absorbs exudate and binds it within the absorbent core. Exudate removal eliminates inhibitory factors from the wound, e.g. proteases. The increased absorption capacity of RespoSorb Super reduces the required frequency of dressing changes. This promotes wound rest and provides additional protection against contamination. Apart from its absorbent quality, RespoSorb Super also has a padding effect.
Intended purpose	Single-use sterile superabsorbent dressing for long-term treatment of injured skin, acute and chronic wounds, with moderate to high levels of exudate. It is used on adults only in healthcare professionals in clinical or homecare environments and can be combined with local antiseptics, primary and secondary dressings.
Indications	RespoSorb Super is suitable for the treatment of superficial, moderately to severely exuding wounds: acute wounds (traumatic wounds, post-operative wounds, lymphatic wounds, arterial leg ulcers) and chronic wounds (decubitus/pressure sores, venous or mixed leg ulcers, tumor wounds). RespoSorb Super can be used under compression bandages.
Mode of application	• Select RespoSorb Super to match the wound size so that the dressing extends at least 1–2 centimeters beyond the wound margins. <ul style="list-style-type: none"> • The wound dressing may not be cut to size. • Place the white side and local antiseptics away from the wound so that the green special nonwoven layer faces away from the wound. • Secure the dressing e.g. with adhesive tape, conforming bandages, or, if necessary, compression bandages. • One dressing can remain on the wound for up to 5 days, depending on the condition of the wound. Change the dressing if clinically indicated or when exudate reaches the edges of the dressing or is noticeable through the green nonwoven top layer. • In case of infected wounds, a combination with e.g. Atramaun Ag is possible. • With decreasing amounts of exudate, the use of a suitable hydroactive wound dressing (e.g. HydroTac) is recommended.
Contraindications	Do not use RespoSorb Super on dry wounds or on exposed bones, muscles or tendons. Do not use RespoSorb Super in case of hypersensitivity to any of its components.
Precautions	• Do not cut the dressing. <ul style="list-style-type: none"> • Medical assessment of the wound condition and the causes of wound-healing impairment is necessary before treating wounds with an impaired healing tendency. Treatment with RespoSorb Super cannot replace a causal treatment of the wound-healing impairment. If there are clinical signs of infection, the infection needs to be controlled with appropriate treatment before this dressing can be used. • In all cases, follow the established clinical protocol. • In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary, on these population groups this dressing should be used with caution and following a clinician's recommendation. • For dressings which are packed in paper-paper peel packs: The sealed seam of the peel pack contains natural rubber latex, which may cause allergic reactions!
Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.	
Incident reporting	For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/45/EU on Medical Devices), if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.
Product disposal	In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to the applicable local rules, regulations and infection prevention standards. Dispose of the medical device with household waste.
Date of revision of the text: 2020-03-11	
<p> GB — PAUL HARTMANN Pty. Ltd. · Macquarie Park, NSW 2113 GB — PAUL HARTMANN LTD. · Heywood/Greater Manchester O10 21T ZA — HARTMANN South Africa · Northriding, 2169 </p>	

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