## Schürfwunden Abrasions Égratignures Escoriazioni Schaafwonden

 PAUL HARTMANN AG  
Paul-Hartmann-Straße 12  
89522 HEIDENHEIM, GERMANY

 Distributor  
AT – PAUL HARTMANN Ges.m.b.H.  
2355 Wiener Neudorf  
BE – N.V. PAUL HARTMANN S.A.  
1480 Saintes/Sint-Renelde  
CH – IVF HARTMANN AG  
8212 Neuhausen  
FR – Lab. PAUL HARTMANN S.à.r.l.  
Châtenois · 67607 Sélestat CEDEX  
GB – PAUL HARTMANN Ltd.  
Heywood/Greater Manchester OL10 2TT  
IT – PAUL HARTMANN S.p.A.  
37139 Verona  
NL – PAUL HARTMANN B.V.  
6546 BC Nijmegen  
ZA – HARTMANN South Africa  
Northriding, 2162

www.hartmann.info

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#### Produktbeschreibung

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### GB Instructions for use

#### Product description

The hydrocolloid plaster is a non-sterile, single-use, self-adhesive plaster with a hydrocolloid formulation and occlusive properties (moist wound treatment). The product is only meant for the treatment of humans by lay users.

#### Composition

The plaster consists of one layer of adhesive hydrocolloid material (containing carboxymethyl-cellulose) and another, larger layer of PU film, acting as the dressing backing material. The PU film is coated by an acrylic adhesive layer. The silicone paper protects the coated side of the plaster.

#### Properties

Whilst traditional plasters work via dry wound healing, which allows the covered tissue to dehydrate and dry out, the hydrocolloid plaster is based on the mode of action of moist wound healing (occlusive properties). Self-adhesive hydrocolloid dressings contain a layer which absorbs exudate and expands into a gel-like mass, thus providing cushioning for the treated damaged area, and helping to relieve pressure and friction.

#### Intended use

The hydrocolloid plaster is a single-use self-adhesive plaster for protection and treatment of abrasions, and for the absorption of exudate.

#### Contraindications

The hydrocolloid plaster must not be used on untreated infected abrasions, major injuries, highly bleeding or exudative wounds, for pressure ulcers or any other use apart from those mentioned in the intended use section.

#### Side effects

Adverse reactions may occur in patients who are allergic or hypersensitive to any of the product's components.

#### Product disposal

The disposable components and used plasters can be disposed of normally with household waste.

#### Mode of application

- Wash hands and ensure the affected area is clean, dry and oil free before applying the plaster.
- Open the peel pack and extract the plaster.
- Warm the plaster with the hands before applying.
- Remove the first protective paper and apply the exposed part of the plaster, to the affected skin area.
- Then gently pull the second protective paper and apply the rest of the plaster. Gently mould the plaster edges to ensure optimum adhesion.
- The application maximum of the hydrocolloid is up to 4 days per individual product. Should you wish to remove the plaster early or at the maximum time of application, first wet it briefly with warm water.

#### Special precautions

- For single use only. Not reusable.
- Only apply the plaster if the hydrocolloid pad completely covers the abrasion.
- Consult a doctor before applying the plaster on a potentially infected abrasion or shows signs of infection while plasters are in place (e.g. redness, overheating, swelling or pain).

- In cases of diabetes mellitus always consult a doctor before using plasters.
- Do not apply the product for more than 4 days.
- Store plasters in a cool (not in the fridge), dry place, out of the reach of children.

### Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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