# Vivano® Med Foam Kit

Wundverbandset für die Unterdruck-Wundtherapie  
Wound dressing set for negative pressure wound therapy  
Kit de pansements pour le traitement de plaies par pression négative  
Wondverbandset voor wondtherapie met negatieve druk  
Kit di medicazione per il trattamento di ferite tramite pressione negativa  
Kit de apósitos para el tratamiento de heridas mediante presión negativa  
Kit de pensos para a terapia por pressão negativa  
Conjunto de curativos para tratamento de feridas por pressão negativa tópica  
Σετ επιθέματος τραυμάτων για τη θεραπεία τραυμάτων με αρνητική πίεση  
Set krytí ran pro podtlakovou terapii  
Súprava krytia na rany na podtlakovú terapiu  
Zestaw opatrunkowy do podciśnieniowej terapii ran  
Kötserkészlet negatívnyomás-terápiás sebkezeléshez  
Набор перевязочный для терапии ран отрицательным давлением  
Превързочен комплект за рани за терапия на рани с негативно налягане  
Komplet potrošnog materijala za terapiju rana negativnim tlakom  
Set pansament pentru tratamentul plăgilor prin presiune negativă  
Komplet obvezilnega materiala za zdravljenje ran z negativnim tlakom  
用於負壓傷口治療的傷口敷料套裝

PAUL HARTMANN AG  
Paul-Hartmann-Straße 12  
89522 Heidenheim, Germany  
www.hartmann.info

[www.vivanosystem.info](http://www.vivanosystem.info)

0 1 2 3

060 252/6 (190620)

Vivano®  
Safety. And Simplicity.

# Deutsch

<b>DE</b>	Deutsch	3
<b>EN</b>	English	12
<b>FR</b>	Français	20
<b>NL</b>	Nederlands	29
<b>IT</b>	Italiano	38
<b>ES</b>	Español	47
<b>PT</b>	Português (PT)	55
<b>BR</b>	Português (BR)	64
<b>GR</b>	Ελληνικά	72
<b>CZ</b>	Česky	81
<b>SK</b>	Slovensky	89
<b>PL</b>	Polski	97
<b>HU</b>	Magyar	106
<b>RU</b>	Русский	114
<b>BG</b>	Български	123
<b>HR</b>	Hrvatski	132
<b>RO</b>	Română	140
<b>SI</b>	Slovenščina	148
<b>HK</b>	中文	156

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CH – IVF HARTMANN AG  
8212 Neuhausen

## Intended use

Foam dressing kit for negative pressure wound therapy intended to transfer negative pressure to the wound bed and to transport exudates from the wound.

The VivanoMed Foam wound dressing kits are only intended for use in conjunction with the Vivano System by PAUL HARTMANN AG or the ATMOS S042 NPWT negative pressure therapy system.

The VivanoMed Foam Kit is meant only for use in humans.

The VivanoMed Foam Kit may be used in hospitals, as well as in the homecare environment.

**IMPORTANT:** VivanoMed Foam Kits may only be used by a physician or a qualified person, as per the law of your country, in accordance with physician's instructions.

## Kit components identification

VivanoMed Foam

Hydrofilm transparent film dressing

VivanoTec Port

## Indications

The VivanoMed Foam is used in wounds with injured tissue to support healing by secondary intention. VivanoMed Foam may be used on intact skin and wounds healing by primary intention when direct contact with the underlying structure is prevented by a suitable wound contact layer.

Treatment with VivanoMed Foam in conjunction with the Vivano System is appropriate in the following types of wounds:

- Chronic wounds
- Acute wounds
- Traumatic wounds
- Dehiscent wounds
- Infected wounds
- Ulcers (e.g. pressure or diabetic ulcers)
- Grafts
- Partial thickness burns
- Surgical incisions (only with a suitable wound contact layer e.g. Atrauman Silicone)

## Contraindications

Contraindications for the use of the Vivano System:

- Malignant tumor wounds
- Non-enteric/unexplored fistulas
- Untreated osteomyelitis
- Necrotic tissue

**IMPORTANT:** The VivanoMed Foam must not be applied directly over exposed nerves, anastomosis points, blood vessels or organs.

Using the VivanoMed Foam Kit is not authorized in patients diagnosed with cerebrospinal fluid leakage, metastatic disease and bleeding diathesis.

**NOTE:** For more information on a particular contraindication, please refer to the **Warnings** and **Precautions** sections of this document.

## Warnings

Please pay attention to the following warnings related to the use of the VivanoMed Foam Kit:

## Bleeding

**NOTE:** The Vivano System was not developed for the prevention or stopping of bleeding.

**IMPORTANT:** Switch the negative pressure therapy unit off immediately, undertake haemostatic measures and inform the attending physician.

**NOTE:** Regardless of using the negative pressure wound therapy, certain medical conditions favor occurrence of bleeding complications.

The following circumstances increase the risk of possibly fatal bleeding, if not controlled with appropriate care:

- Surgical sutures and/or anastomoses
- Trauma
- Irradiation
- Inadequate haemostasis
- Wound infection
- Treatment with anticoagulants or coagulation inhibitors
- Protruding bone fragments or sharp edges

Patients with increased risk of bleeding complications should be monitored with additional level of care, under the responsibility of the supervising physician.

**IMPORTANT:** In patients with diagnosed acute bleeding, coagulation disorders or being treated with anticoagulants, the **800 ml canister** should not be used for the exudate collection. Instead, a 300 ml canister should be utilized. Such practice enables more frequent monitoring of the patient by the health professionals, thereby reducing the potential risk of excessive blood loss.

## Malignant tumor wounds

Negative pressure wound therapy for malignant tumor wounds is contraindicated, as it is linked to the risk of enhanced tumor formation by proliferation support effect. However, it is considered as legitimate in a palliative context. For patients in the end-of-life stage in whom a complete cure is no longer the aim, the improvement of their quality of life by controlling the three most disabling elements: the odor, exudate,

and pain associated with changing the dressings, outweighs the risk of accelerating the spread of tumors.

## Non-enteric/unexplored fistulas

Application of the wound dressing on non-enteric or unexplored fistulas is contraindicated, as it may damage the intestinal structures and/or organs.

## Untreated osteomyelitis

Application of the wound dressing in wounds with untreated osteomyelitis is contraindicated, as it may result in the spreading of the infection.

## Necrotic tissue

Application of the wound dressing on necrotic tissue is contraindicated, as it may lead to local spreading of the infection.

## Application of the VivanoMed Foam on nerves, anastomosis points, blood vessels or organs

The VivanoMed Foam must not be applied directly over exposed nerves, anastomosis points, blood vessels or organs as it may result in deterioration of the underlying structures.

## Special precautions

Please pay attention to the following precautions:

### Infected wounds

Wound dressings should be changed at regular intervals, according to instructions listed in the **Dressing change** section of this document.

Infected wounds should be monitored more frequently, and may require wound dressings to be changed more often.

**NOTE:** For more information on wound monitoring in the context of negative pressure wound therapy, please refer to the **Monitoring** section of this document.

Typical signs of wound infection are redness, swelling, itching, increased warmth of the wound itself or in the neighbouring area, foul odour etc.

Infected wounds may trigger systemic infection, manifesting by high fever, headache, dizziness, nausea, vomiting, diarrhea, disorientation, erythroderma etc.

The consequences of systemic infection can be fatal.

**IMPORTANT:** If there is any suspicion of either local or systemic infection, contact the supervising physician and consult if negative pressure wound therapy should be stopped, or an alternative therapy should be considered.

### Blood vessels and organs

Blood vessels and organs should be adequately protected by means of fascias, tissues or other types of protective layers placed above them.

**IMPORTANT:** Special precautions need to be taken, when dealing with infected, weakened, irradiated or sutured blood vessels or body organs.

### Bone fragments or sharp edges

Protruding bone fragments and sharp edges should be removed, or covered accordingly, before using the VivanoMed Foam, as they may damage blood vessels or body organs and cause bleeding.

**NOTE:** For more information related to bleeding in the context of negative pressure wound therapy, please refer to the **Bleeding** section of this document.

### Surgical incisions

Application of the VivanoMed Foam on surgical incisions may be performed only with a suitable wound contact layer e.g. Atrauman Silicone.

### Enteric fistulas

In case of treatment of wounds containing explored enteric fistulas, additional level of precautions need to be implemented, if negative pressure wound therapy is to be applied. Presence of enteric fistula in close proximity to the wound increases the risk of wound contamination and/or infection. In order to mitigate the risk associated with the potential contact

of intestinal content with the wound, enteric fistula needs to be surgically separated, following local guidelines or established surgical practices.

### Spinal cord injuries with the development of autonomic hyperreflexia

Discontinue the negative pressure wound therapy if the patient has spinal cord injuries with the development of autonomic hyperreflexia.

### Magnetic resonance imaging

The VivanoTec / VivanoTec Pro unit is not considered as MRI safe and must not be used in close proximity to an MRI unit.

### Defibrillation

The VivanoTec / VivanoTec Pro unit must be disconnected if resuscitation of the patient using a defibrillator is necessary.

### Hyperbaric oxygen therapy (HBO)

The VivanoTec / VivanoTec Pro unit must be disconnected for patients undergoing hyperbaric oxygen therapy, as its use is considered a potential fire hazard.

## General precautions

Please pay attention to the following precautions:

### Damaged, expired or contaminated product

Do not use any of the components of the kit, in case of damage, expiration or any suspicion of contamination. It may cause overall decline in the therapeutic efficiency, wound contamination and/or infection.

### Single use only

All disposable components of the VivanoMed Foam Kit are meant for single use only. Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

### Resterilization

Components of VivanoMed Foam Kit are provided

sterile. Do not resterilize any of the components as it may cause overall decline in the therapeutic efficiency of the kit and potentially lead to wound contamination and/or infection.

### Product disposal

To minimize the risk of potential infection hazards or environmental pollution, disposable components of VivanoMed Foam Kit should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

### Safety measures for infection prevention

Implement and apply adequate personal protection and institutional infection control measures, when handling the components of the VivanoMed Foam Kit (e.g. the use of sterile gloves, masks, gowns etc.)

**IMPORTANT:** Before and after using the stopper on the connector, it must be cleaned and disinfected.

### Patients population

No general restrictions for the usage of VivanoMed Foam Kits in different patient populations (e.g. adults and/or children) is given. However, the VivanoTec / VivanoTec Pro unit has not been evaluated for use in paediatrics.

**IMPORTANT:** Before prescribing its use on a child, weight and height along with the overall health condition must be medically evaluated first.

### Patient's health status

Weight and general condition of the patient should be taken into consideration during any application of negative pressure wound therapy.

### Compatibility with other systems

The VivanoMed wound dressings are compatible with Atrauman Silicone.

### Dressing size

The size of the dressing must be adapted to the size of the wound being treated with the negative pressure wound therapy.

Improper wound dressing size may either cause maceration and disintegration of the periwound tissue, or wound margins drying out and inefficient transfer of exudate.

**NOTE:** For more information on complications related to excessive coverage of the intact skin, please refer to the **Dressing of the intact skin** section of this document.

**IMPORTANT:** In order to provide optimal conditions for negative pressure wound therapy, the film dressing should cover about 5 cm of intact skin around the wound.

### Dressing placement

Use only dressings coming directly from sterile packages.

Do not force the placement of the foam, as it may lead to direct tissue damage or subsequent wound healing delay or even local necrosis, due to an elevated level of compression.

**IMPORTANT:** Always record the number of foams used for each wound.

The number of layers of film in the dressing can be adapted to each medical condition. Placement of multiple layers of film increases the risk of tissue maceration and consequently tissue irritation.

**IMPORTANT:** In case of tissue irritation, due to the use of multiple layers of the film, discontinue the Vivano negative pressure wound therapy.

### Dressing removal

**IMPORTANT:** Always record the number of foams removed from the wound in order to guarantee the removal of all foams introduced.

Foam left in the wound for a longer period of time than indicated in the **Dressing change** section may cause growth of granulation tissue into the foam. This may increase the difficulty of dressing change and may promote wound infection, among other medical complications.

# English

Dressing changes may possibly lead to disruption of the new granulation tissue, which may result in bleeding.

**IMPORTANT:** Implement additional protective measures, when changing the dressing of patients with identified increased risk of bleeding.

**NOTE:** For more information related to bleeding in the context of negative pressure wound therapy, please refer to the **Bleeding** section of this document.

## Disconnection from the VivanoTec / VivanoTec Pro unit

The decision on how long the patient can be disconnected from the VivanoTec / VivanoTec Pro unit is a clinical assessment that has to be made by the attending physician.

The time interval for safe interruption of the therapy is strongly dependent on the overall patient and wound status, as well as the exudate composition and amount of exudate extracted per time unit.

Long interruption might lead to exudate retention and local maceration effects, as well as to a blocked wound dressing due to coagulation effects within the foam matrix. Lack of the effective barrier between the wound and the non-sterile environment increases the risk of infection.

**IMPORTANT:** Do not leave the dressing with the VivanoTec / VivanoTec Pro unit being switched off for prolonged time periods. In case of leaving the dressing for a longer time, it is recommended for a physician to perform an evaluation of the wound condition along with the overall health status of the patient. According to the physician's evaluation, either rinsing of the wound along with a dressing change, or switching to an alternative therapy is recommended.

## Intermittent pressure mode

Intermittent pressure, as compared to continuous pressure, may be used for enhancing local perfusion

and granulation formation, if tolerable for the patient, the patients' health and wound condition. However, the continuous therapy is generally recommended for the treatment of patients with an increased risk of bleeding, acute enteric fistulas, highly exuding wounds, or when wound bed stabilization is required.

## Pressure settings

**PRECAUTION:** Pressure settings below 50 mm Hg may potentially lead to exudate retention and decreased therapeutic efficiency.

**PRECAUTION:** High pressure settings may increase the risk for micro-trauma, hematoma and bleeding, local hyperfusion, tissue damage or fistula formation.

Correct pressure setting for the Vivano negative pressure wound therapy must be decided by the supervising physician and should be based upon the exudate output, overall patient status as well as recommendations from therapeutic guidelines.

## Dressing of the intact skin

Dressing of the intact skin should be about 5 cm margin around the wound. Prolonged or repeated dressing of bigger areas might result in tissue irritation.

**IMPORTANT:** In case of tissue irritation, discontinue the Vivano negative pressure wound therapy.

Application of the wound dressing on intact skin may create wrinkles on the dressing surface. Formation of wrinkles significantly increases the risk of dressing leakage, and as a consequence occurrence of infection.

**IMPORTANT:** Special care should be taken when applying the wound dressing to the peri-wound area's fragile skin.

## Dressing of wounds susceptible to irritation

For wounds susceptible to constant irritation (close proximity to limbs), continuous (rather than intermittent) therapy is indicated.

## Circumferential dressings

Circumferential dressings should be performed under medical supervision. Lack of adequate protective measures might cause local hypoperfusion.

## Dressings in the vicinity of the vagus nerve

Dressings in the vicinity of the vagus nerve should be performed under medical supervision, as its stimulation might cause bradycardia.

## Allergies

Application of the Vivano negative pressure wound therapy is not recommended if the patient is allergic to any component of the VivanoMed wound dressing.


**IMPORTANT:** Hydrofilm has an acrylic adhesive coating which might present a risk of adverse reactions in patients who are allergic or hypersensitive to acrylic adhesives.

## Dressing application

Adequate personal protection and institutional infection control measures must be implemented, when applying the dressing. Always apply the components of the VivanoMed Foam Kit in the following order:

1. Before applying the dressing for the first time, and after each dressing change, the wound must be thoroughly cleaned and conditioned in accordance with the physician's instructions.


2. Cut the VivanoMed Foam with sterile scissors or scalpel, so that it fits the size and shape of the wound to be treated.

**NOTE:** In case of using pre-perforated foam, follow the same preparation procedure. The perforations may be used as a guidance to determine size and shape.

**IMPORTANT:** In order to minimize the risk of foam pieces falling into the wound, do not cut the foam in the wound vicinity.

3. Introduce the foam material into the wound.

**NOTE:** If necessary, multiple foam pieces may be applied on top of one another.

**IMPORTANT:** Make sure that no small foam pieces detach and contaminate the wound, when foam is being introduced into the wound.


**IMPORTANT:** Always record the number of foams used for each wound.

4. In order to seal the wound filled with the VivanoMed Foam, the Hydrofilm film dressing should cover about 5 cm of intact skin around the wound. If necessary, the film dressing can be cut to shape.

**IMPORTANT:** Make sure that all release liners are removed from the film dressing, as leaving them might result in lower breathability of the film dressing.


**IMPORTANT:** Do not use any moisturizing/lipid transferring skin disinfection solution, before application of the film, as it may impair its adhesion to skin.

5. Apply the film dressing in accordance with the printed instructions and mould it to shape in order to achieve an airtight cover.


**IMPORTANT:** Do not stretch the film dressing during application as it may cause skin blisters and/or shear damage in attached skin areas.

6. Select a point in the center of the film dressing, where the VivanoTec Port is to be applied and cut a hole of approx. 2 x 4 cm. In doing so, make sure that the foam material is not perforated.


7. Apply the VivanoTec Port so that the center of the port is precisely above the prepared hole in the film dressing.

**IMPORTANT:** Make sure that the port holes are perfectly aligned with the prepared hole in the film dressing to guarantee optimal functionality of the system.

8. Connect the drainage tube to the canister tube of the Vivano System and commence the therapy in accordance with the user manual.

**IMPORTANT:** Improper connection of the canister to the Vivano System may result in impaired negative pressure performance, which may prolong the overall wound healing process.

**IMPORTANT:** Check that dressing collapses upon application of negative pressure.

## Monitoring

**IMPORTANT:** The monitoring frequency must be adapted according to the overall patient's health status and the treated wound condition, evaluated by the supervising physician.

When treating the wound with negative pressure using VivanoMed Foam Kit it should be periodically verified, if the nominal exudate flow/volume

corresponds with the actual flow/volume. Additionally, the patient's health status and the wound condition should be closely monitored. The wound condition should be monitored via checking of the wound dressing for signs of leakage, and the wound margins, as well as the exudate for signs of infection.

**IMPORTANT:** If there are any signs of infection, the attending physician must be notified immediately.

## Dressing change

Dressing change must occur at regular time intervals, usually taking place between 48 to 72 hours. In case of infected wounds, the dressing must be changed at intervals of 12 to 24 hours (or even more frequently if necessary).

**IMPORTANT:** It is a prerogative of the supervising physician to decide, what frequency of wound monitoring and dressing change should be applied in a particular case. It is based upon the medical evaluation of the wound site, as well as the patient's health status.

**NOTE:** For more information regarding monitoring of infected wounds, please refer to the **Infected wounds** section of this document.

Dressing change should be performed in the following order:

**IMPORTANT:** Ensure any bleeding has been stopped prior to dressing change. The attending physician must decide, if negative pressure wound therapy may be safely continued or an alternative therapy should be applied.

1. Each time the dressing is changed the wound must be thoroughly cleansed and conditioned in accordance with the physician's instructions and applicable hygiene standards.
2. Inspect the wound and make sure that all pieces of the VivanoMed Foam have been removed.

3. Apply the new dressing according to the instructions listed in the **Dressing application** section of this document.

**NOTE:** In case of adherence of the dressing to the wound, saline solution can be added to the foam material to enable dressing removal. After 15–30 minutes, the dressing should be carefully removed from the wound.

In case of wounds with a known tendency to adhesion, the use of a non-adhering wound contact layer should be considered. If the patient complains of pain while the dressing is being changed, premedication, the use of a non-adhering wound contact layer or the administration of a local anesthetic should be considered.

## Special notes

Keep out of the reach of children.

## Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

-  Medical device
-  Manufacturer
-  Date of manufacture
-  Use-by date
-  Batch code
-  Catalogue number
-  Consult instructions for use

-  Do not re-use
-  Do not re-sterilize
-  Do not use if package is damaged
-  Keep dry
-  Keep away from sunlight
-  Caution
-  Sterilized using ethylene oxide
-  Single sterile barrier system
-  Double sterile barrier system
-  Unique Device Identifier

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HK – PAUL HARTMANN Asia-Pacific Ltd.