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Produktbeschreibung

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EN Instructions for use

Product description

Hydrocoll is a non-active, self-adhesive hydrocolloid dressing consisting of an adhesive hydrocolloid layer covered with a polyurethane film. The side of the dressing in contact with the wound is covered with a release liner composed of siliconised paper. It does not contain any medicinal substances, tissue or cells of human or animal origin or nanomaterials. It is sterilized using irradiation.

Composition

Hydrocoll consists of a gas permeable polyurethane cover. The hydrocolloid compound itself mainly contains sodium carboxymethylcellulose, a synthetic rubber mixture based on polystyrene polyisoprene copolymers and polyisobutylene.

Properties and mode of action

The product absorbs wound exudate. After the exudate has been absorbed, the hydrocolloid layer transforms into a gel which provides a moist wound healing environment. The ingress of germs is prevented by the polyurethane film. The film is water-repellent and allows moisture vapor and gas to pass through. Since the gel is produced during the application, the dressing does not adhere to the wound and the healing process is not impaired by the desing change. Newly formed tissue is not damaged.

Due to the thin flexible material, Hydrocoll molds easily to body contours. The surface of Hydrocoll is dirt-repellent and water-repellent.

Intended purpose

Hydrocoll is a sterile, single-use, non-invasive hydrocolloid dressing for the use on intact and injured human skin by healthcare professionals.

Hydrocoll (normal), Concave, Sacral can be used for the treatment of medium to slightly exuding acute and chronic wounds, burns up to grade II or wounds at skin removal sites and for prevention of facial pressure injuries.

Hydrocoll Thin can be used for the treatment of slightly exuding wounds.

Hydrocoll Concave can be used on the heel and elbow, Hydrocoll Sacral on the sacral area.

The product can remain in place up to 7 days. The accumulated use can exceed 30 days.

Modo de aplicacón Seleccione el correct size and shape of Hydrocoll. The dressing should exceed the wound by at least one centimetre to ensure that it stays firmly in place. If needed, cut the dressing to the size of the wound with sterile scissors.

Hydrocoll / Hydrocoll Thin: remove the protective product cover, place Hydrocoll on the wound, adjust the edges to fit and press on.

Hydrocoll Sacral: remove the protective product cover and press the dressing into the anal cleft. The tip of the dressing can point upwards or downwards depending on the position and depth of the ulcer. Press the edges on firmly.

Hydrocoll Concave: remove the largest part of the protective product cover, place the middle part of the Hydrocoll dressing over the wound, remove the rest of the protective product cover, then overlap the wings to fit.

Due to the absorption of exudate a whitish gel is formed that is visible on the top as a discoloration and a blister. Hydrocoll should be changed if the blister has almost reached the edge of the dressing. Should exudate leak from the dressing, it must be changed immediately. To change the dressing, lift Hydrocoll at the edges and remove the whole dressing carefully.

It may happen that a layer of hydrocoll gel remains on the wound. It will not impair the healing of the wound and can be rinsed off if necessary.

The wear time of the dressing depends on the condition of the wound and the level of exudate but can remain on the wound for up to 7 days.

Contraindications

Hydrocoll is contraindicated in any patients who may be allergic to any of its ingredients.

Hydrocoll must not be used on wounds extending into muscle, tendon or bone. It must not be applied to clinically infected wounds or third-degree burns.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Special precautions

In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary on these population groups, this dressing should be used with caution following a clinician's recommendation.

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Before opening the product packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/45 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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DE - From Information: 2024-11-21 PT – PAUL HARTMANN S.p.A. - 37139 Verona