
Vivano[®]Med

White Foam

PVA Schaumverband für die Unterdruck-Wundbehandlung
PVA foam dressing for negative pressure wound treatment
Pansement mousse en APV pour le traitement de plaies par pression négative
PVA schuimverband voor de vacuümwondbehandeling
Medicazione in schiuma di PVA per il trattamento di ferite tramite pressione negativa
Apósito de espuma de PVA para el tratamiento de heridas mediante presión negativa
Penso de espuma de PVA para o tratamento de feridas por pressão negativa
Curativo de espuma APV para tratamento de feridas por pressão negativa
Αφρώδης επίδεσμος PVA για τη θεραπεία τραυμάτων με αρνητική πίεση
Pénové krytí z PVA pro podtlakovou terapii
Penové krytie z PVA pre podtlakovú terapiu
Opatrunek piankowy z polialkoholu winylowego (PVA)
do podciśnieniowej terapii leczenia ran
PVA habkótszer negatívnyomás-terápiás sebkezeléshez
Полимерная губка на основе ПВС (PVA) для вакуумной обработки ран
ПВА пенообразна превръзка за лечение на рани с вакуум терапия
PVA spruža za terapiju rana negativnim tlakom
Pansament spuma din PVA pentru tratamentul plagilor prin presiune negativa
PVA obloga iz penastega materiala za zdravljenje ran z negativnim tlakom
PVA skumförband för sårbehandling med negativt tryck
Polyvinylalkohol-skumforbinding til sårbehandling med negativt tryk
Haavan alipainemuhoidon PVA-vaahtomuovisidos
PVA-skumforbinding for undertrykksbehandling av sår
Negatif Basıncı yardimli yara Tedavisi İçin PVA Köpük yara Örtüsü
ضمادة زبدية مصنوعة من بولي فينيل الكحول لمعالجة الجروح بالضغط السالب
用於負壓傷口治療的聚乙烯醇 (PVA) 發泡塑膠敷料

Vivano[®]
Safety. And Simplicity.

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CH – IVF HARTMANN AG
8212 Neuhausen

Instructions for use

VivanoMed White Foam wound dressings may only be used by a physician or in accordance with medical instructions.

The VivanoMed White Foam wound dressings are intended for use in conjunction with the VivanoTec negative pressure therapy system by PAUL HARTMANN AG or with the ATMOS S042 NPWT system.

Properties and mode of action

The VivanoMed White Foam wound dressings consist of pre-moistened hydrophilic polyvinyl alcohol (PVA). To apply a complete dressing for negative pressure therapy (in conjunction with a VivanoTec unit) at least the following further components are required:

- Hydrofilm film dressing
- VivanoTec Port

In conjunction with the VivanoTec negative pressure therapy system, the VivanoMed White Foam wound dressing ensures reliable removal of wound exudate from the wound. In doing so, the applied negative pressure is evenly distributed on the wound base. The Hydrofilm film dressing ensures a reliable airtight seal of the wound. Thereby, in conjunction with the VivanoTec Port system, it ensures optimum use of the negative pressure and transport of the exudate.

Indications

VivanoMed White Foam can be used in all skin defects, where primary closure is not possible or not desired.

VivanoMed White Foam is intended for sensitive acute and chronic wounds. The dense structure prevents ingrowth of granulation tissue, and delicate structures are protected. VivanoMed White Foam is especially suitable to treat sensitive wounds in combination with the Vivano Negative Pressure Wound Therapy System.



Application notes for use in combination with the VivanoTec NPWT system

VivanoMed White Foam is indicated in patients who can benefit from the application of a controlled negative pressure and the wound healing that it

promotes. In particular this includes patients for whom drainage through negative pressure and the removal of infectious material through continuous or intermittent negative pressure may be of advantage. VivanoMed White Foam protects sensitive structures. Treatment with VivanoMed White Foam in conjunction with the VivanoTec negative pressure therapy system is appropriate in the following types of wound: Chronic, acute, traumatic, dehiscent wounds, ulcers (e.g. pressure or diabetic ulcers), grafts and burns, where granulation ingrowth is not desired. Additionally, VivanoMed White Foam can be used in explored undermined and tunnelling wounds.

Application of the dressing

- 1) Before applying the dressing for the first time, and after each dressing change, the wound must be thoroughly cleansed in accordance with the physician's instructions.
- 2) Cut the VivanoMed White Foam with sterile scissors so that it fits the size and shape of the wound to be treated. Then introduce the foam material into the wound. When cutting the foam material to shape make sure that no small pieces of foam fall into the wound.



English

EN


- 3) Seal the wound with the VivanoMed White Foam (and additionally with VivanoMed Foam if necessary) and then with the Hydrofilm film dressing. It must be ensured that the film dressing covers at least 5 cm of intact skin all around the wound. If necessary, the film dressing can be cut to shape. Apply the film dressing in accordance with the printed instructions (numbering) and mould it to shape in order to achieve an airtight covering.


- 4) Then select a point on the film dressing where the VivanoTec Port is to be applied and cut a hole of approx. 2 cm in diameter into the film. In doing so make sure that the foam material is not perforated.

- 5) Apply the VivanoTec Port so that the centre of the port is precisely above the prepared hole in the film dressing.

- 6) Then connect the drainage tube to the canister tube of the VivanoTec negative pressure therapy system and commence the therapy in accordance with the user manual.

Monitoring the wound

In order to guarantee that the therapy is safe, the wound dressing has to be checked frequently. The wound dressing must be checked for imperviousness and the wound margins as well as the exudate for signs of infection. If there are any signs of infection, the attending physician must be notified immediately.

Dressing change

Each time the dressing is changed the wound must be thoroughly cleansed in accordance with the physician's instructions and applicable hygiene standards. If the patient complains of pain while the dressing is being changed, premedication, the use of a non-adhering wound contact layer or the administration of a local anaesthetic should be considered.

The dressing must be changed at regular intervals of 48 to 72 hours. In the case of infected wounds, the dressing should be changed at intervals of 12 to 24 hours. On removing the foam material, it must be ensured that all the pieces of foam that have been introduced are removed.

Contraindications

Contraindications are present in patients with malignant tumour wounds, non-enteric and unexamined fistulas, osteomyelitis or necrotic tissue.

Special precautions

Special precautions are indicated in patients with acute bleeding, coagulation disorders and patients being treated with anticoagulants. The 800 ml canister should not be used here. In the event of sudden or intensified bleeding, the treatment must be discontinued immediately. In these cases, the attending physician must be informed and appropriate haemostatic measures taken.

Special precautions under medical supervision are also indicated in the case of:

- Circumferential dressings
- Use in the vicinity of the vagus nerve
- Patients who are to undergo any of the following therapies: hyperbaric oxygen therapy, defibrillation.

- Vicinity of MRI units: The VivanoTec negative pressure unit is not compatible with MRI (magnetic resonance imaging) units. The VivanoTec negative pressure unit should therefore not be operated in the vicinity of an MRI unit.
- Enteric fistulas

In general, the weight and general condition of the patient should be taken into consideration when carrying out treatment with negative pressure therapy.

Special precautions must be taken in the case of infected, weakened, irradiated or sutured blood vessels or organs. Do not use VivanoMed White Foam directly on organs.

Bone fragments or sharp edges can damage vessels and organs and must therefore be covered or removed before using VivanoMed White Foam.

Wounds with **intestinal fistulas** must only be treated with VivanoMed White Foam in conjunction with special precautions.

In the case of **movable parts** of the body, which cannot be kept still, continuous (rather than intermittent) therapy is indicated in order to minimise irritations.

If the patient has **spinal cord injuries** with the development of **autonomic hyperreflexia**, the therapy must be discontinued.

The decision on how long the patient can be disconnected from the VivanoTec negative pressure unit is a clinical assessment that has to be made by the attending physician.

Warning

The wound must not be too firmly or forcefully filled with VivanoMed White Foam as this can increase the risk of harm to health.

Special notes

Store in a dry place and at a temperature not below 5 °C and not above 30 °C.

Do not sterilise. Keep out of the reach of children.

Contents are sterile unless package is damaged

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HK – PAUL HARTMANN Asia-Pacific Ltd.
CN – PAUL HARTMANN (Shanghai) Trade Co., Ltd.