Vivano® Tec

Y-Connector

Y-Connector für die Unterdruck-Wundbehandlung
Y-Connector for negative pressure wound therapy
Connecteur en Y pour le traitement de plaies par pression négative
Y-connector voor vacuümwondtherapie
Connettore a Y per il trattamento di ferite tramite pressione negativa
Conector en «Y» para el tratamiento de heridas mediante presión negativa
Conector em Y para o tratamento de feridas por pressão negativa
Conector em Y para tratamento de feridas por pressão negativa tópica (PNT)
Y-Connector για τη θεραπεία τραυμάτων με αρνητική πίεση
Y konektor pro podtlakovú terapii
Y konektor pre podtlakovú terapiu
Złącze Y do podciśnieniowego leczenia ran
Y csatlakozó negatívnyomás-terápiás sebkezeléshez
Разъем Y-образный для вакуумной терапии ран
Y-конектор за лечение на рани с вакуум терапия
Y konektor za terapiju rana negativnim tlakom
Conector în Y pentru tratamentul plăgilor prin presiune negativă
Y priključek za zdravljenje ran z negativnim tlakom
Negatif Basıncı Yardımlı Yara Tedavisi İçin Y-Konnektör
用於負壓傷口治療的 Y-Connector
موصل ثلاثي لمعالجة الجروح بالضغط السلبي

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CH – IVF HARTMANN AG
8212 Neuhausen

Intended use

The VivanoTec Y-Connector is used to connect two VivanoTec Ports to one VivanoTec / VivanoTec Pro on one patient.

The VivanoTec Y-Connector is applied in particular in patients where two wounds should be treated with one negative pressure unit or in patients with large wounds for which one port is not enough.

The VivanoTec Y-Connector is only intended for use with the VivanoTec Port and in conjunction with the Vivano System by PAUL HARTMANN AG or the ATMOS S042 NPWT negative pressure wound therapy system.

NOTE: To apply a complete dressing at least the following further components are required:

- VivanoMed Foam
- Hydrofilm film dressing
- VivanoTec Port

The VivanoTec Y-Connector is meant only for use in humans.

The VivanoTec Y-Connector may be used in hospitals, as well as in the homecare environment.

IMPORTANT: The VivanoTec Y-Connector may only be used by a physician or a qualified person, as

per the law of your country, in accordance with physician's instructions.

Indications

Treatment with VivanoTec Y-Connector in conjunction with the Vivano System is appropriate in the following types of wounds:

- Chronic wounds
- Acute wounds
- Traumatic wounds
- Dehiscent wounds
- Infected wounds
- Ulcers (e.g. pressure or diabetic ulcers)
- Grafts
- Partial thickness burns
- Surgical incisions (only with a suitable wound contact layer e.g. Atrauman Silicone)

Contraindications

Contraindications for the use of the Vivano System:

- Malignant tumor wounds
- Non-enteric/unexplored fistulas
- Untreated osteomyelitis
- Necrotic tissue

Two different wounds in the same patient may be connected by means of the Y connector only if:

- both wounds have the same aetiology, and
- neither of the two wounds is infected.

IMPORTANT: The VivanoMed Foam must not be applied directly over exposed nerves, anastomosis points, blood vessels or organs.

Using the VivanoMed Foam Kit is not authorized in patients diagnosed with cerebrospinal fluid leakage, metastatic disease and bleeding diathesis.

NOTE: For more information on a particular contraindication, please refer to the **Warnings** and **Precautions** sections of the instructions of use of the other Vivano System components.

Special precautions

Please pay attention to the following precautions:

If two non-contiguous wounds are managed by means of a VivanoTec Y-Connector, the VivanoTec / VivanoTec Pro can detect possible blockades and pressure deviations only to a limited degree. For this reason, the tube system must be checked regularly for kinks.

The marking on the VivanoTec Y-Connector indicates the VivanoTec Port at which pressure deviations and possible blockades cannot be detected.

For special precautions, please also observe the instructions for use of the other complementary Vivano products.

General precautions

Please pay attention to the following precautions:

Damaged, expired or contaminated product

Do not use in case of damage, expiration or any suspicion of contamination. It may cause overall decline in the therapeutic efficiency, wound contamination and/or infection.

Single use only

The VivanoTec Y-Connector is meant for single use only. Re-use may result in wound contamination and/or infection, as well as general failure of the wound healing process.

Re-sterilization

The VivanoTec Y-Connector is provided sterile. Do not re-sterilize as it may cause overall decline in the therapeutic efficiency and potentially lead to wound contamination and/or infection.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, the VivanoTec Y-Connector should be disposed of according to local regulations and infection prevention standards.

Safety measures for infection prevention

Implement and apply adequate personal protection and institutional infection control measures, when handling the VivanoTec Y-Connector (e.g. the use of sterile gloves, masks, gowns etc.)

IMPORTANT: Before and after using the stopper on the connector, it must be cleaned and disinfected.

Patients population

No general restrictions for the usage of the VivanoTec Y-Connector in different patient populations (e.g. adults and/or children) is given. However, the VivanoTec / VivanoTec Pro has not been evaluated for use in paediatrics.


IMPORTANT: Before prescribing its use on a child, weight and height along with the overall health condition must be medically evaluated first.

Patient's health status


Weight and general condition of the patient should be taken into consideration during any application of negative pressure wound therapy.

Application notes


Application of the dressing:

1. Apply the dressings and the VivanoTec Ports in accordance with the instructions for use of the VivanoMed Foam Kits.

2. Connection of the two VivanoTec Ports with the VivanoTec Y-Connector.

3. Subsequently connect the VivanoTec Y-Connector with the canister tube of the Vivano System and start the therapy in accordance with the instructions for use.

Monitoring

IMPORTANT: The monitoring frequency must be adapted according to the overall patient's health status and the treated wound condition, evaluated by the supervising physician.

When treating the wound with negative pressure using the VivanoTec Y-Connector it should be periodically verified, if the nominal exudate flow/volume corresponds with the actual flow/volume. Additionally, the patient's health status and the wound condition should be closely monitored. The wound condition should be monitored via checking of the wound dressing for signs of leakage, and the wound margins, as well as the exudate for signs of infection.

IMPORTANT: If there are any signs of infection, the attending physician must be notified immediately.

Special notes

Keep out of the reach of children.

If VivanoTec Y-Connectors are applied, please observe the instructions for use of the other complementary Vivano products.

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