
Vivano[®]Tec

Exudate Canister

Exsudatbehälter für die Unterdruck-Wundbehandlung
Exudate canister for negative pressure wound therapy
Réservoir à exsudats pour le traitement de plaies par pression négative
Exsudaatcanister voor vacuümtherapie
Contentitore dell'essudato per il trattamento di ferite tramite pressione negativa
Recipiente para exudado para el tratamiento de heridas mediante presión negativa
Recipiente para exsudado para o tratamento de feridas por pressão negativ
Recipiente para exsudato para tratamento de feridas por pressão negativa tópica
Δοχείο συλλογής εκκρίσεων για τη θεραπεία τραυμάτων με αρνητική πίεση
Sběrná nádoba na exsudát pro podtlakovou terapii
Zberná nádoba na exsudát pre podtlakovú terapiu
Zbiornik na wydzielinę do podciśnieniowego leczenia ran
Váladéktároló tartály negatívnyomás-terápiás sebkezeléshez
Контейнер для сбора экссудата для вакуумной терапии ран
Канистър за екссудат за лечение на рани с вакуум терапия
Posuda za eksudat za terapiju rana negativnim tlakom
Canistră pentru exsudat pentru tratamentul plăgilor prin presiune negativă
Zbirnik za eksudat za zdravljenje ran z negativnim tlakom
Negatif basınç yardımlı yara tedavisi için eksüda toplama kabı
用於負壓傷口治療的滲出物收集器
حاوية الإفرازات للعلاج بالضغط السالب

Vivano[®]
Safety. And Simplicity.

Deutsch

Vorsicht!

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Caution!

Important safety instructions

- Before use, check that the exudate canister and the connection tubing are not damaged.
- The exudate canister is sterile and can be used in the sterile operating environment.
- The exudate canister may only be changed by specialist medical personnel during the therapy.
- Components marked with ② must not be reused. In the event of reuse, these components lose their function resulting in a high risk of infection.
- In order to guarantee uninterrupted operation of the negative pressure unit, the filling level of the exudate canister should be regularly checked. Once the maximum filling level has been reached, it should be replaced with a fresh exudate canister.
- Avoid the risk of the tube becoming displaced or blocked by regularly checking the tubing system and its connections for leaks and kinks.
- The VivanoTec Exudate Canisters are only intended for use in conjunction with the Vivano System by PAUL HARTMANN AG or the ATMOS S042 NPWT negative pressure wound therapy system.
- Regularly monitor the patient, the unit and the wound dressing. Watch out for wound exudate, maceration, infection and loss of vacuum. In order to ensure a safe therapy, it is necessary to check the wound dressing frequently. In doing so, be sure to check the wound dressing for imperviousness, and also for negative pressure, check the wound edges for maceration and check the wound edges and the exudate for signs of infection. In the event of signs of an infection, the attending physician must be informed immediately.
- The exudate canister has to be changed after each patient, but after one week at the latest.
- The negative pressure unit incl. exudate canister must always remain in an upright position during operation.
- The negative pressure unit incl. exudate canister must not be placed on a patient's bed.

Contraindications

Contraindications for the use of the Vivano System:

- Malignant tumor wounds
- Non-enteric / unexplored fistulas
- Untreated osteomyelitis
- Necrotic tissue

IMPORTANT: The Vivano System must not be applied directly over exposed nerves, anastomosis points, blood vessels or organs.

Special precautions

- Special precautionary measures are indicated in patients with acute bleeding, coagulation disorders and patients being treated with anticoagulants. The 800 ml exudate canister should not be used in these cases. In the event of sudden or intensified bleeding, the therapy must be discontinued immediately. In such cases, the attending physician must be informed and appropriate haemostatic measures taken.
- Blood vessels and organs must be adequately protected by means of overlying fascias, tissues or other types of protective layers placed above them. Special precautionary measures must be undertaken in the case of infected, weakened, irradiated or sutured blood vessels or organs.
- Bone fragments or sharp edges can damage vessels or organs and must therefore be covered or removed before the application of negative pressure therapy.
- Wounds with intestinal fistulas may only be treated with the Vivano System in the presence of special precautionary measures.
- In the case of movable parts of the body, which cannot be kept still, continuous (rather than intermittent) therapy is indicated in order to minimise irritations.
- If the patient has spinal cord injuries, with the development of autonomic hyperreflexia, the therapy must be discontinued.
- The VivanoTec / VivanoTec Pro unit must not be used in the presence of strong magnetic fields (such as an induction stove) and not used near application of HF surgical equipment.

English

- The VivanoTec / VivanoTec Pro unit has not been evaluated for use in paediatrics. Before prescribing its use on a child, the weight and height of the patient must first be checked.
- The VivanoTec / VivanoTec Pro unit is not compatible with magnetic resonance imaging (MRI). The VivanoTec / VivanoTec Pro unit should therefore not be operated in the vicinity of an MRI unit.
- The VivanoTec / VivanoTec Pro unit must be removed from the wound dressing and from the patient during any hyperbaric oxygen therapy or defibrillation procedure.
- Circumferential dressings should be performed under medical supervision. Lack of adequate protective measures might cause local hypoperfusion.
- Dressing application in the vicinity of the vagus nerve should be performed under medical supervision, as its stimulation might cause bradycardia.

Intended use

VivanoTec Exudate Canister serves the collection and safe storage of secretion and serves the protection from contamination of patient, user and environment. The filter system of the Exudate Canister protects the unit from contamination and filters the air for the flushing.

Indications for use

The Vivano System is used in wounds with injured tissue to support healing by secondary intention. The Vivano System may be used on intact skin and wounds healing by primary intention when direct contact with the underlying structure is prevented by a suitable wound contact layer.

Treatment with VivanoTec Exudate Canister in conjunction with the Vivano system is appropriate in the following types of wounds:

- Chronic wounds
- Acute wounds
- Traumatic wounds
- Dehiscent wounds
- Infected wounds
- Ulcers (e.g. pressure or diabetic ulcers)
- Grafts
- Partial thickness burns
- Surgical incisions (only with a suitable wound contact layer e.g. Atrauman Silicone)

The VivanoTec Exudate Canister is meant only for use in humans.

The VivanoTec Exudate Canister may be used in hospitals, as well as in the homecare environment.

Users

Users are specialist medical personnel. The exudate canister in combination with the negative pressure unit may only be used by a physician or by qualified and authorised personnel acting on the instructions of a physician. Users must have knowledge and understanding of the special medical application of the negative pressure unit.

Function

The negative pressure unit generates a vacuum in the wound dressing, in the tubing system and in the exudate canister, by means of which the exudate can be aspirated off. The exudate is collected in the exudate canister. The filter technology prevents fluids, bacteria and odours from penetrating into the interior of the device.

General precautions

Please pay attention to the following precautions:

Wound conditioning

Before applying the dressing for the first time, and after each dressing change, the wound must be thoroughly cleaned and conditioned in accordance with the physician's instructions.

Damaged, expired or contaminated product

Do not use in case of damage, expiration or any suspicion of contamination. It may cause overall decline in the therapeutic efficiency, wound contamination and / or infection.

Single use only

The VivanoTec Exudate Canister is meant for single use only. Re-use may result in wound contamination and / or infection, as well as general failure of the wound healing process.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, the VivanoTec Exudate Canister should be disposed of according to local regulations and infection prevention standards.

Safety measures for infection prevention

Implement and apply adequate personal protection and institutional infection control measures, when handling the VivanoTec Exudate Canister (e.g. the use of sterile gloves, masks, gowns etc.)

IMPORTANT: Before and after using the stopper on the connector, it must be cleaned and disinfected.

Patients population

No general restrictions for the usage of VivanoTec Exudate Canister in different patient populations (e.g. adults and / or children) is given. However, the Vivano System has not been evaluated for use in paediatrics.

IMPORTANT: Before prescribing its use on a child, weight and height along with the overall health condition must be medically evaluated first.

Patient's health status

Weight and general condition of the patient should be taken into consideration during any application of negative pressure wound therapy.

Disconnection from the VivanoTec / VivanoTec Pro unit

The decision on how long the patient can be disconnected from the negative pressure unit incl. exudate canister is a clinical assessment that has to be made by the attending physician.

English

Application

The activities in this chapter may only be carried out by specialist medical personnel.

IMPORTANT: The exudate canisters of the negative pressure unit are sterile components and can therefore be used in a sterile surgical environment.

Inserting the exudate canister

- Carefully remove the exudate canister from the sterile peel package.

IMPORTANT: The attached tubing piece must not be allowed to fall onto a non-sterile surface.

- Insert the exudate canister at a slight angle into the guidance on the negative pressure unit **1**.
- Tilt the exudate canister in the direction of the negative pressure unit until it fully engages with the blue unlock **2** key.
- Pull gently on the exudate canister to make sure that it is securely attached to the negative pressure unit.

Removing the exudate canister

- Removing the exudate canister.

See "Connecting / removing the negative pressure unit to / from the wound dressing set".

- Press the blue unlock key **2** on the negative pressure unit.
- Slightly tilt the exudate canister and remove it.
- Properly dispose of the exudate canister. Observe local regulations.

Connecting and removing the negative pressure unit to and from the wound dressing


Caution!

- Make sure that the tube connectors are always correctly connected to each other in order to prevent malfunction.
- For application of the wound dressing set, follow the instructions for use of the wound dressing set.


Connecting

Attach the connectors (tube ends) of the exudate canister to the connectors of the wound dressing set **1**.


Removing

Press the unlocking device at the side of the connector and keep it pressed **2**.

Pull the two ends apart **3**.

Disposal

Before disposal, break off the stopper from the connector and plug it into the exudate lumen **4**.

This ensures that no exudate can escape from the canister.

Ensure correct disposal of the exudate canister, observing the locally applicable regulations.

English

Special notes

Keep out of the reach of children

Date of revision of the text: 2018-06

Symbols

Warning, observe carefully

This device is CE certified

Do not use if package is damaged or open

Do not re-use

Keep dry

Keep away from sunlight

Dispose of packaging in an environmentally responsible way

Dispose of packaging in an environmentally responsible way

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