# Thermoval® standard

Predalček za baterijo  
Pretinac za bateriju  
Pil yuvası  
電池格  
Отделение за батерии  
Battery compartment

Toplotno tipalo  
Mjerni senzor  
Ölçü sensörü  
探頭  
Сензор  
Thermometric sensor

Gumb za vklop/izklop  
Uključno/isključni gumb  
Аçma-Kapama tuşu  
開/關鍵  
Бутон ON/OFF  
ON/OFF button


Konica termometra  
Vrh toplomjera  
Termometre ucu  
探熱針尖  
Врѣхче на термометъра  
Thermometer tip

## Normal temperature

36.6 – 38.0 °C

35.5 – 37.5 °C

34.7 – 37.3 °C

### SI Razlaga piktogramov

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#### ■ English

Congratulations on your purchase of this HARTMANN product. Please read the instructions carefully before using the thermometer for the first time, and keep them in a safe place. This product is intended for the measurement of human body temperature.

**Explanation of symbols**  
▼: battery is empty  
 Err: for errors, misreadings  
 Lo °C: temperature under 32 °C  
 Hi °C: temperature over 42.9 °C

<span></span>	Observation of the operating instructions
<span></span>	Protection against electric shock (type B)

<span></span>	Disposal note
<b>Methods of measuring temperature</b> <b>In the rectum (rectal)</b>	
This is by far the most accurate method from a medical point of view, because it comes closest to the core body temperature. The thermometer tip is inserted carefully into the rectum to a maximum of 2 cm. The usual measurement time is approximately 40 to 60 seconds.	

**Under the arm (axillary)**  
 Placing the thermometer in the armpit provides a measurement of surface temperature that can fluctuate by around 0.5 °C to 1.5 °C from rectal temperature readings in adults. The usual measurement time for this method is approximately 90 seconds. It should be noted, however, that an exact reading cannot be obtained if, for example, the armpits have been allowed to cool. If this is the case, we recommend extending the measurement time to 5 minutes in order to obtain the most precise possible reading.

**In the mouth (oral)**  
 There are different heat zones in the mouth. As a general rule, the oral temperature is 0.3 °C to 0.8 °C lower than the rectal temperature. To ensure that the reading is as accurate as possible, place the thermometer tip under the tongue to the left or right of the root of the tongue. The sensor must have constant contact with the tissue during the reading and be placed under the tongue in one of the two heat pockets at the back. Keep the mouth closed during the reading and breathe evenly through the nose. Do not eat or drink anything before the measurement. The usual measurement time is approximately 60 seconds.

**Note:**  
 However, please note that Thermoval standard is a maximum thermometer, i.e. the displayed temperature can increase slightly if measurement continues after the beeps. We strongly recommend the rectal method as the most accurate method for identifying the basal temperature, and advise you to extend the measuring time by 3 minutes after the beeps have sounded.

**Cleaning and disinfection**  
 The best way to clean the thermometer tip is by applying disinfectant (e.g. 70 % ethyl alcohol) with a damp cloth.

- Safety precautions**
- Do not allow the device to come into contact with hot water.
  - Do not expose to high temperatures or direct sunlight.
  - Do not drop the thermometer. It is neither shock-proof nor impact-resistant.
  - Do not bend or open the device (except the battery compartment).
  - Do not clean with thinners, petrol or benzene! Only clean with water or disinfectants.
  - Do not soak the thermometer in liquid.
  - The thermometer contains small parts (battery, battery cover) that can be swallowed by children. Do not therefore leave the thermometer unattended when children are present.
  - If the ambient temperature is above 35 °C, immerse the thermometer tip in cold water for approximately 5 to 10 seconds before measuring the temperature.
  - A persistently high temperature, particularly in children, must be treated medically – please contact your doctor!

**Battery replacement**  
 The battery is empty and needs replacing when the symbol ▼ appears on the bottom right of the digital display. Remove the battery cover and replace it with a battery of the same type. Please note: the + sign on the battery must point downwards. Close the battery compartment, making sure that the catch clicks into place.

**Product disposal**  
 Please respect the environment. Thermometers and batteries should not be disposed of in household waste. Hand them into collection centres or municipal recycling centres for special waste.

**Technical data**  
 Type: Maximum thermometer  
 Measurement range: 32 °C to 42.9 °C

Accuracy:	+/- 0.1 °C in the range 35.5 °C to 42 °C +/- 0.2 °C for the remaining measurement range
Display:	3-digit LCD display, (display unit 0.1 °C)
Storage/ transportation temperature:	–25 °C to + 55 °C; relative humidity of ≤ 95 <span> </span> %

Ambient temperature during use:	10 °C to 35 °C; relative humidity of ≤ 80 <span> </span> %
Battery type:	Alkaline-manganese battery, type LR 41, 1.55 V, battery life ~3,000 measurements Approx. 11 g incl. battery
Weight:	PAUL HARTMANN AG, Paul-Hartmann-Straße, 89522 Heidenheim, Germany
<span></span>	<span></span>

**Legal requirements and guidelines**  
 This thermometer complies with the Medical Devices Directive 93/42/EEC, the European Standard for Medical Thermometers EN 12470-3: 2000 + A1: 2009, and carries the CE mark.

The CE certification confirms that this is a medical device with a measuring function in the sense of the Medical Devices Act (Medizinproduktegesetz), which has undergone a conformity assessment procedure. A "Notified Body" (TÜV Produkt Service GmbH) confirms that this product fulfils all the appropriate statutory regulations.

ELECTRICAL MEDICAL DEVICES require special precautionary measures with regard to electromagnetic compatibility (EMC).

The use of this device in the vicinity of portable mobile high-frequency communication devices or other devices with strong electromagnetic fields can adversely affect its operation.

**Calibration check**  
 We recommend a calibration check at intervals of two years in the case of professionally used devices, for example in pharmacies, medical practices or hospitals. You should also observe the national regulations determined by the legislator, such as, in Germany, the "Medizinprodukte-Betreiberverordnung" (Medical Device Operating Regulation).

**Warranty**  
 This high quality product is covered by a 2-year warranty from the date of purchase, subject to the conditions listed below. Claims must be presented with proof of purchase within the warranty period. This product has been manufactured with the greatest care in accordance with international quality guidelines. Should you still have cause for complaint, please send the instrument, together with proof of purchase, to the service address below. Batteries are excluded from the warranty. The warranty does not cover damage resulting from improper use or unauthorised operation. Claims for compensation are limited to the value of the merchandise; compensation for subsequent damage or injury is not covered.

AE – PAUL HARTMANN	
Middle East	
FZE Dubai	
AU – PAUL HARTMANN Pty. Ltd.	
Level 6, 5 Ryder Boulevard	
Rhodes, NSW 2138	
Australia	
GB – PAUL HARTMANN Ltd.	
Heywood/Lancashire OL10 2TT	
HK – PAUL HARTMANN	
Asia-Pacific Ltd.	
Hong Kong	
ZA – HARTMANN South Africa	
2194 Johannesburg	

Date of revision of the text: 2013-12