
HARTMANN

Vivano[®]Tec

Y-Connector

Y-Connector für die Unterdruck-Wundtherapie
Y-Connector for negative pressure wound therapy
Connecteur en Y pour la thérapie par pression négative
Y-Connector voor wondtherapie met negatieve druk
Y-Connector per terapia con pressione negativa
Y-Connector para la terapia de presión negativa en heridas
Conector em Y para a terapia de feridas por pressão negativa
Y-Connector para a terapia de feridas por pressão negativa
Y-Connector για τη θεραπεία τραυμάτων με αρνητική πίεση
Konektor Y-Connector pro podtlakovou terapii ran
Konektor Y-Connector na podtlakovú terapiu
Złącze Y do podciśnieniowej terapii ran
Y-Connector negatívnyomás-terápiás sebkezeléshez
Разъем Y-Connector для терапии ран отрицательным давлением
Y-конектор за терапия на рани с негативно налягане
Y-konektor za terapiju rana negativnim tlakom
Conector in Y pentru terapia plăgilor cu presiune negativă
Y-priključek za terapijo z negativnim tlakom
用於負壓傷口治療的 Y-Connector

DE Deutsch	3
EN English	7
FR Français.....	11
NL Nederlands.....	15
IT Italiano	19
ES Español	23
PT Português (PT).....	27
PT Português (BR).....	31
GR Ελληνικά.....	35
CZ Česky.....	39
SK Slovensky.....	43
PL Polski	47
HU Magyar.....	51
RU Русский.....	55
BG Български.....	59
HR Hrvatski	63
RO Română.....	67
SI Slovenščina	71
HK 中文.....	75

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Stand der Information: 2022-01-25

Intended use

VivanoTec Y-Connector is used to connect two VivanoTec Ports to one VivanoTec / VivanoTec Pro unit on one patient.

VivanoTec Y-Connector is applied in particular in patients where two wounds should be treated with one negative pressure unit or in patients with large wounds for which one port is not enough.

VivanoTec Y-Connector is only intended for use in conjunction with the Vivano System by PAUL HARTMANN AG or the ATMOS S042 NPWT negative pressure therapy system.

NOTE: To establish a functional negative pressure wound therapy dressing, at least the following further components are required:

- VivanoMed Foam
- Hydrofilm transparent film dressing
- VivanoTec Port
- VivanoTec Exudate Canister
- NPWT control unit (e.g. VivanoTec Pro)

VivanoTec Y-Connector is meant only for use in humans.

VivanoTec Y-Connector may be used in hospital, residential care and home-care settings.

IMPORTANT: VivanoTec Y-Connector may only be used by a physician or a qualified person, as per the law of your country, in accordance with physician's instructions.

Indications

Treatment with VivanoTec Y-Connector in conjunction with the Vivano System is appropriate in the following types of wounds:

- Chronic wounds
- Acute wounds
- Traumatic wounds
- Dehiscent wounds
- Infected wounds
- Ulcers (e.g. pressure or diabetic ulcers)
- Grafts
- Partial thickness burns
- Surgical incisions (only with a suitable wound contact layer e.g. Atrauman Silicone)

Contraindications

Contraindications for the use of the Vivano System:

- Malignant tumor wounds
- Non-enteric/unexplored fistulas
- Untreated osteomyelitis
- Necrotic tissue

Two different wounds in the same patient may be connected by means of the Y-Connector only if:

- both wounds have the same aetiology, and
- neither of the two wounds is infected.

IMPORTANT: VivanoMed Foam must not be applied directly over exposed nerves, anastomosis points, blood vessels or organs.

Using VivanoMed Foam Kit is not authorized in patients diagnosed with cerebrospinal fluid leakage, metastatic disease and bleeding diathesis.

NOTE: For more information on a particular contraindication, please refer to the **Warnings** and **Precautions** sections of the instructions for use of the other Vivano System components.

Special precautions

Please pay attention to the following precautions:

If two non-contiguous wounds are managed by means of a VivanoTec Y-Connector, the VivanoTec / VivanoTec Pro unit can detect possible blockages and pressure deviations only to a limited degree. For this reason, the tube system must be checked regularly for kinks.

The marking **1** on VivanoTec Y-Connector indicates the VivanoTec Port at which pressure deviations and possible blockages cannot be detected.

For special precautions, please also observe the instructions for use of the other Vivano System components.

General precautions

Please pay attention to the following precautions:

Damaged, expired or contaminated product

Do not use in case of damage, expiration or any suspicion of contamination. It may cause overall decline in the therapeutic efficiency, wound contamination and/or infection.

Before opening the product packaging, inspect it carefully for signs of damage. In particular, examine the integrity of the packaging surfaces and the entire seal area completely under suitable lighting conditions (white light or daylight). Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Single use only

VivanoTec Y-Connector is meant for single use only. Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Resterilization

VivanoTec Y-Connector is provided sterile. Do not resterilize as it may cause overall decline in the therapeutic efficiency and potentially lead to wound contamination and/or infection.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Safety measures for infection prevention

Implement and apply adequate personal protection and institutional infection control measures, when handling VivanoTec Y-Connector (e.g. the use of sterile gloves, masks, gowns etc.)

IMPORTANT: Before and after using the stopper on the connector, it must be cleaned and disinfected.

Patients population

No general restrictions for the usage of VivanoTec Y-Connector in different patient populations (e.g. adults and/or children) are given. However, the VivanoTec / VivanoTec Pro unit has not been evaluated for use in paediatrics.

IMPORTANT: Before prescribing its use on a child, weight and height along with the overall health condition must be medically evaluated first.


Patient's health status

Weight and general condition of the patient should be taken into consideration during any application of negative pressure wound therapy.


Application notes

Application of the dressing:


1. Apply the dressings and the VivanoTec Ports in accordance with the instructions for use of the VivanoMed Foam Kits.

2. Connection of the two VivanoTec Ports with VivanoTec Y-Connector.

3. Subsequently connect VivanoTec Y-Connector with the canister tube of the Vivano System and start the therapy in accordance with the instructions for use.

Monitoring

IMPORTANT: The monitoring frequency must be adapted according to the overall patient's health status and the treated wound condition, evaluated by the supervising physician.

When treating the wound with negative pressure using VivanoTec Y-Connector it should be periodically verified, if the nominal exudate flow/volume corresponds with the actual flow/volume. Additionally, the patient's health status and the wound condition should be closely monitored. The wound condition should be monitored via checking of the wound dressing for signs of leakage, and the wound margins, as well as the exudate for signs of infection.

IMPORTANT: If there are any signs of infection, the attending physician must be notified immediately.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Special notes

Keep out of the reach of children.

If VivanoTec Y-Connectors are applied, please observe the instructions for use of the other Vivano System components.

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