

Pütter® binde E

Kräftige Kurzzugbinden
Strong short-stretch bandages
Bandes résistantes à allongement court
Sterke zwachtels met korte rek
Bende a corta trazione forte
Vendas resistentes de curta elasticidad
Ligaduras resistentes de curta tração
Ανεκτικοί επιδέσμοι μικρής τόνωσης
Mocne opaski o niskiej rozciągliwości
Bandaje rezistente cu elasticitate redusă
Starka kortsträcksbindor
Kraftigt kortstræksbind
Vahvat vähäelastiset siteet
Sterke kortstrekking

1

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EN Instructions for use

General description of the device

Pütterbinde E* is a brown-coloured, strong short-stretch bandage with textile elasticity and an extensibility of approximately 85 % provided with an end finder. The non-active device is available in different widths and lengths.

Intended purpose

The compression bandages are non-sterile medical devices used by healthcare professionals on intact skin of humans (above 6 months of age) in clinical and home environment. They can be used for long-term compression therapy (ABPI 0.9-1.3) (decongestion phase) of:

- varicose veins
- chronic venous insufficiency (CVI) with and without venous leg ulcers (VLU)
- thrombosis/thrombophlebitis
- chronic oedemata of other origin
- lymph oedemata

In addition, they can be used as support and relief bandages in musculoskeletal injuries. The compression bandages can also be used in combination with primary and secondary dressings.

Composition

100 % cotton

Properties

Extensibility: approx. 85 %.

Washable: Wash the bandage in a regular wash cycle at a temperature of 95 °C.

Indications

The compression bandages can be used for compression therapy (ABPI 0.9-1.3) (decongestion phase) for the treatment of:

- varicose veins
- chronic venous insufficiency (CVI) with or without venous leg ulcers (VLU)
- thrombosis/thrombophlebitis (e.g. thrombophlebitis as well as condition after healed thrombophlebitis, deep leg and arm vein thrombosis, condition after deep vein thrombosis, post-thrombotic syndrome, thromboprophylaxis)

- chronic oedemata of other origin
- lymph oedemata

The devices can also be used as support and relief bandage in musculoskeletal injuries.

Intended users

Pütterbinde E* is designed to be used by healthcare professionals who have received appropriate training.

Contraindications

The short-stretch compression bandages must not be used in case of:

- advanced peripheral arterial occlusive disease (PAOD)
- decompensated heart failure (NYHA III + IV)
- phlegmasia coerulea dolens
- ankle-brachial pressure index (ABPI) > 1.3 and ≤ 0.5
- septic phlebitis
- florid erysipelas

Precautions

Special precautionary measures are required in the event of:

- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral polyneuropathy)
- pronounced weeping dermatoses
- incompatibility with bandage material
- primary chronic polyarthritis

The application technique, e.g. with regards to stretch and/or padding, may have to be adapted to the patient's anatomical and pathophysiological characteristics, especially in case of fragile or elderly skin.

If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. The right measurement method must be considered for the ABPI. During the use of the compression bandages, the ABPI must be controlled regularly. The therapy must be adapted accordingly.

If an inelastic system is created by close to full stretch application, high ankle pressures should be closely monitored by a physician.

The devices should not be used on infants below the age of 6 months.

When applying the compression bandages on children (from 6 months on), the pressure should be reduced appropriately.

If at any time during therapy the patient develops severe pain, numbness, tingling or noticeable discoloration of the toes, a healthcare professional must be contacted immediately, and the compression bandages must be removed.

Application (see pictograms)

Modified bandaging technique according to Pütter, if not otherwise prescribed by the doctor. (1)

Further remarks

Pütterbinde E* is indicated for use by healthcare professionals who have received appropriate training. Before applying Pütterbinde E, clean the wound and cover it with a suitable wound dressing. If required, excess padding bandage material may be used for additional padding of bony prominences and tissue recesses (e.g. shinbone, Bisgaard's backdrop). If additional fixation is required, use tape but no bandage clips. After application, the compression bandage must be checked (e.g. circulation in the patient's toes, pain reported by the patient) in order to be certain that the pressure of the bandage is not too high.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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* © MEDICE-Iserlohn
PAUL HARTMANN AG