HARTMANN

# Vivano® Med

## Gel Strip

Gelstreifen

Gel strip

Bande gel

Gelstrip

Gel in striscia

Banda de gel

Tira de gel

Gel em tira

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Stand der Information: 2022-01-27

## Intended use

VivanoMed Gel strip is a sterile, adhesive, cuttable, single-use hydrogel slab intended for use with the Vivano negative pressure wound therapy (NPWT) system from PAUL HARTMANN AG. It enables the creation of an airtight seal around difficult-to-dress wounds, and it secures at or near the wound site certain fixtures (e.g., drains, tubes) required for the therapy.

VivanoMed Gel Strip is intended for application on intact skin or on already applied wound dressings.

VivanoMed Gel Strip is suitable for use on the general human population (adults and children; however, the Vivano System has not been evaluated for use in paediatrics).

VivanoMed Gel Strip may be used only by or under the supervision/guidance of a physician, but may be applied in both hospital and homecare settings.

## Indications

VivanoMed Gel Strip is indicated for use during NPWT to:

- seal dressings around difficult-to-dress wound sites
- attach fixtures at/near a wound site
- create effective seals around drains

## Contraindications

Contraindications for the use of the Vivano System (incl. VivanoMed Gel Strip):

- wounds in which there is a malignancy in the wound bed or margin
- open or exposed joints and untreated osteomyelitis
- open peritoneal spaces, pleural cavity or non-enteric/unexplored fistulas
- exposed organs, arteries or veins, or anastomosis points
- exposed nerves
- wounds covered with necrotic tissue, eschar or more than 30% slough

VivanoMed Gel Strip should not be used on patients who may be allergic to any of its ingredients.

In case of an allergic reaction, remove the product immediately and clean the affected area.

**IMPORTANT:** VivanoMed Gel Strip must not be applied directly over exposed nerves, anastomosis points, organs, arteries or veins.

## Warnings and precautions

The product is intended for single use only.

Do not use after expiry date.

Discontinue use if excessive discomfort or irritation develops.

All offcuts of the product must be discarded. Using the product together with creams or lotions can reduce its adhesion.

Ensure that both plastic film backings have been removed before use.

Keep out of reach of children.

Special care should be taken when applying the wound dressing to the peri-wound area's fragile skin.

If VivanoMed Gel Strip is difficult to remove, irrigate the product with a sterile saline solution or any other suitable sterile rinsing solution and then remove.

Implement and apply adequate personal protection, institutional infection control measures and if applicable aseptic working conditions, when handling VivanoMed Gel Strip (e.g. the use of sterile gloves, masks, gowns etc.).

## Product disposal

To minimize the risk of potential infection hazards or environmental pollution, disposable components of VivanoMed Gel Strip should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

## Application

1. Ensure that the skin on which the gel strip is to be applied is clean and dry.
2. Remove the sterile gel strip from the peel pouch.

**IMPORTANT:** Before opening the product packaging, inspect it carefully for signs of damage.

In particular, examine the integrity of the packaging surfaces and the entire seal area completely under suitable lighting conditions (white light or daylight).

Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

3. If necessary, the gel strip can be cut to the required size with sterile scissors.
4. Remove the 1<sup>st</sup> plastic film backing before applying the gel strip.
5. Apply the gel strip on the respective area of the skin.
6. Remove the 2<sup>nd</sup> plastic film backing.
7. Apply the transparent film dressing (Hydrofilm) for the sealing case, if necessary.
8. Check the sealing and/or fixing function at the start as well as during therapy.

## Monitoring

The wound condition should be monitored via checking of the wound dressing for signs of leakage, and the wound margins, as well as the exudate for signs of infection. The wound dressing must be checked for imperviousness and the wound margins as well as the exudate for signs of infection. If there are any signs of infection, the attending physician must be notified immediately.

## Dressing change

Dressing change must occur at regular time intervals, usually taking place between 48 to 72 hours. In case of infected wounds, the dressing must be changed at intervals of 12 to 24 hours (or even more frequently if necessary).

**IMPORTANT:** It is a prerogative of the supervising physician to decide, what frequency of wound monitoring and dressing change should be applied in a particular case. It is based upon the medical evaluation of the wound site, as well as the patient's health status.

## Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Consult instructions for use
	Caution
	Unique Device Identifier
	Sterilized using irradiation
	Single sterile barrier system
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Keep dry
	Keep away from sunlight

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