HARTMANN

Vivano[®]Tec

Exudate Canister

Exsudatbehälter für die Unterdruck-Wundtherapie
Exudate canister for negative pressure wound therapy
Réservoir à exsudats pour la thérapie par pression négative
Exsudaatcanister voor wondtherapie met negatieve druk
Contentore dell'essudato per terapia con pressione negativa
Recipiente de exudado para la terapia de presión negativa en heridas
Recipiente para exudado para a terapia de feridas por pressão negativa
Coletor de exsudato para terapia de feridas por pressão negativa
Δοχείο συλλογής εκκρίσεων για τη θεραπεία τραυμάτων με αρνητική πίεση
Sběrná nádoba na exsudát pro podtlakovou terapii ran
Zberná nádoba na exsudát na podtlakovú terapiu
Zbiornik na wydzielinę do podciśnieniowej terapii ran
Váladéktároló tartály negatívnyomás-terápiás sebkezeléshez
Контейнер для экссудата для терапии ран отрицательным давлением
Канистър за ексудат за терапия на рани с негативно налягане
Posuda za eksudat za terapiju rana negativnim tlakom
Canistră pentru exsudat pentru terapia plăgilor cu presiune negativă
Zbiralnik za eksudat za terapijo z negativnim tlakom
用於負壓傷口治療的滲出液桶

DE Deutsch	3
EN English	9
FR Français	15
NL Nederlands	21
IT Italiano.....	27
ES Español.....	33
PT Português (PT).....	39
BR Português (BR).....	45
GR Ελληνικά.....	51
CZ Česky.....	57
SK Slovensky	63
PL Polski.....	69
HU Magyar	75
RU Русский	81
BG Български.....	88
HR Hrvatski.....	94
RO Română.....	100
SI Slovenščina.....	106
HK 中文.....	112

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Stand der Information: 2022-04-22

Intended use

The exudate canister is intended to connect VivanoTec Port to the VivanoTec/VivanoTec Pro unit and to collect and store wound exudates during the treatment with negative pressure wound therapy (NPWT).

VivanoTec Exudate Canister is only intended for use in conjunction with the Vivano System by PAUL HARTMANN AG or the ATMOS S042 NPWT negative pressure wound therapy system.

NOTE: To establish a functional negative pressure wound therapy dressing at least the following further components are required:

- VivanoMed Foam
- Hydrofilm transparent film dressing
- VivanoTec Port
- NPWT control unit (e.g. VivanoTec Pro)

The Vivano System is meant only for use in humans.

No general restrictions for the usage of the Vivano System in different patient populations (e.g. adults and/or children) are given. However, the Vivano System has not been evaluated for use in paediatrics.

VivanoTec Exudate Canister may be used in hospital, residential care and home-care settings.

IMPORTANT: The Vivano System may only be used by a physician or a qualified person, as per law of your country, in accordance with the physician's instructions.

Indications

The Vivano System is used in wounds with injured tissue to support healing by secondary intention. VivanoMed Foam may be used on intact skin and wounds healing by primary intention when direct contact with the underlying structure is prevented by a suitable wound contact layer. VivanoTec Exudate Canister is used to connect the port to the negative pressure unit.

Contraindications

Contraindications for the use of the Vivano System:

- malignant tumour wounds
- non-enteric/unexplored fistulas
- untreated osteomyelitis
- necrotic tissue

NOTE: For more information on a particular contraindication, please refer to the **Warnings** and **Precautions** sections of this document.

Warnings

Please pay attention to the following warnings related to the use of VivanoTec Exudate Canister.

Bleeding

NOTE: The Vivano System was not developed for the prevention or stopping of bleeding.

IMPORTANT: In the case of blood occurring suddenly or more frequently on the dressing, in the tubes or in the exudate container switch the negative pressure wound therapy unit off immediately, undertake haemostatic measures and inform the physician.

NOTE: Regardless of using the negative pressure wound therapy, certain medical conditions favour occurrence of bleeding complications.

The following circumstances increase the risk of possibly fatal bleeding, if not controlled with appropriate care:

- surgical sutures and/or anastomoses
- non-sutured haemostatic agents, e.g. spray wound sealant or bone wax
- trauma
- irradiation
- inadequate haemostasis
- wound infection
- treatment with anticoagulants or coagulation inhibitors
- protruding bone fragments or sharp edges

Patients with increased risk of bleeding complications should be monitored with additional level of care, under the responsibility of the physician.

IMPORTANT: In patients with diagnosed acute bleeding, coagulation disorders or being treated with anticoagulants, the **800 ml canister** should not be used for the exudate collection. Instead, a 300 ml canister should be utilized. Such a practice enables more frequent monitoring of the patient by the healthcare professionals, thereby reducing the potential risk of excessive blood loss.

IMPORTANT: When using non-sutured haemostatic agents, additional protective measures should be implemented, to prevent them from accidental displacement. Suitability for negative pressure wound therapy should be assessed by a physician on an individual basis.

Malignant tumour wounds

Negative pressure wound therapy for malignant tumour wounds is contraindicated, as it is linked to the risk of enhanced tumour formation by proliferation support effect. However, it is considered as legitimate in a palliative context. For patients in the end-of-life stage for whom a complete cure is no longer the aim, the improvement of their quality of life by controlling the three most disabling elements: the odour, exudate, and pain associated with changing the dressings, outweighs the risk of accelerating the spread of tumours.

Non-enteric/unexplored fistulas

Application of the wound dressing on non-enteric or unexplored fistulas is contraindicated, as it may damage the intestinal structures and/or organs.

Untreated osteomyelitis

Application of the wound dressing in wounds with untreated osteomyelitis is contraindicated, as it may result in the spreading of the infection.

Necrotic tissue

Application of the wound dressing on necrotic tissue is contraindicated, as it may lead to local spreading of the infection.

Application of VivanoMed Foam on nerves, anastomosis points, blood vessels or organs

VivanoMed Foam must not be applied directly over exposed nerves, anastomosis points, blood vessels or organs as it may result in deterioration of the underlying structures.

Special precautions

Please pay attention to the following precautions.

Infected wounds

Wound dressings should be changed at regular intervals, according to the corresponding instructions for the materials that you use. Infected wounds should be monitored more frequently, and may require wound dressings to be changed more often.

NOTE: For more information on wound monitoring in the context of negative pressure wound therapy, please refer to the **Monitoring** section of the corresponding instructions for the materials that you use.

Typical signs of wound infection are redness, swelling, itching, increased warmth of the wound itself or in the neighbouring area, foul odour etc.

Infected wounds may trigger systemic infection, manifesting by high fever, headache, dizziness, nausea, vomiting, diarrhoea, disorientation, erythroderma etc. The consequences of systemic infection can be fatal.

IMPORTANT: If there is any suspicion of either local or systemic infection, contact the physician and consult if negative pressure wound therapy should be stopped, or an alternative therapy should be considered.

Blood vessels and organs

Blood vessels and organs should be adequately protected by means of fascias, tissues or other types of protective layers placed above them.

IMPORTANT: Special precautions need to be taken, when dealing with infected, weakened, irradiated or sutured blood vessels or body organs.

Bone fragments or sharp edges

Protruding bone fragments and sharp edges should be removed, or covered accordingly, before using VivanoMed Foam, as they may damage blood vessels or body organs and cause bleeding.

NOTE: For more information related to bleeding in the context of negative pressure wound therapy, please refer to the **Bleeding** section of this document.

Surgical incisions

Application of VivanoMed Foam on surgical incisions may be performed only with a suitable wound contact layer e.g. Atrauman Silicone.

Enteric fistulas

In case of treatment of wounds containing explored enteric fistulas, additional level of precautions need to be implemented, if negative pressure wound therapy is to be applied. Presence of enteric fistula in close proximity to the wound increases the risk of wound contamination and/or infection. In order to mitigate the risk associated with the potential contact of intestinal content with the wound, enteric fistula needs to be surgically separated, following local guidelines or established surgical practices.

Spinal cord injuries with the development of autonomic hyperreflexia

Discontinue the negative pressure wound therapy if the patient has spinal cord injuries with the development of autonomic hyperreflexia.

Magnetic Resonance Imaging (MRI)

This device is not considered as MRI safe and must not be used in close proximity to a MRI unit.

Defibrillation

The VivanoTec/VivanoTec Pro unit must be disconnected if resuscitation of the patient using a defibrillator is necessary.

Hyperbaric oxygen therapy (HBO)

The VivanoTec/VivanoTec Pro unit must be disconnected for patients undergoing hyperbaric oxygen therapy, as its use is considered a potential fire hazard.

General precautions

Please pay attention to the following precautions.

Damaged, expired or contaminated product

Do not use any component of the Vivano System in case of damage, expiration or any suspicion of contamination. It may cause overall decline in the therapeutic efficiency, wound contamination and/or infection.

IMPORTANT: Before opening the product packaging, inspect it carefully for signs of damage. In particular, examine the integrity of the packaging surfaces and the entire seal area completely under suitable lighting conditions (white light or daylight). Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Single use only

All disposable components of the Vivano System are meant for single use only. Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Resterilization

The components of the Vivano System that are provided sterile are meant for single use. Do not resterilize any of these components as it may cause overall decline in the therapeutic efficiency of the kit and potentially lead to wound contamination and/or infection.

Product disposal

To minimize the risk of potential infection hazards or environmental pollution, disposable components of the Vivano System should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Safety measures for infection prevention

Implement and apply adequate personal protection and institutional infection control measures, when handling the components of the Vivano System (e.g. the use of sterile gloves, masks, gowns etc.)

IMPORTANT: Before and after using the stopper on the connector, it must be cleaned and disinfected.

Patient population

No general restrictions for the usage of the Vivano System in a different patient population (e.g. adults and/or children) are given. However, the Vivano System has not been evaluated for use in paediatrics.

IMPORTANT: Before prescribing its use on a child, weight and height along with the overall health condition must be medically evaluated first.

Patient's health status

Weight and general condition of the patient should be taken into consideration during any application of negative pressure wound therapy.

Dressing size

The size of the dressing must be adapted to the size of the wound being treated with the negative pressure wound therapy.

Improper wound dressing size may either cause maceration and disintegration of the peri-wound tissue, or wound margins drying out and inefficient transfer of exudate.

NOTE: For more information on complications related to excessive coverage of the intact skin, please refer to the **Dressing on intact skin** section of this document.

IMPORTANT: In order to provide optimal conditions for negative pressure wound therapy, the film dressing should cover about 5 cm of intact skin around the wound.

Dressing placement

Use only dressings coming directly from sterile packages.

Do not force the placement of the foam, as it may lead to direct tissue damage or subsequent wound healing delay or even local necrosis, due to an elevated level of compression.

IMPORTANT: Always record the number of foams used for each wound.

The number of layers of film in the dressing can be adapted to each medical condition. Placement of multiple layers of film increases the risk of tissue maceration and consequently tissue irritation.

IMPORTANT: In case of tissue irritation, due to the use of multiple layers of the film, discontinue the Vivano negative pressure wound therapy.

Dressing removal

IMPORTANT: Always record the number of foams removed from the wound in order to guarantee the removal of all foams introduced.

Foam left in the wound for a longer period of time than indicated in the **Dressing change** section may cause growth of granulation tissue into the foam. This may increase the difficulty of dressing change and may promote wound infection, among other medical complications.

Dressing changes may possibly lead to disruption of the new granulation tissue, which may result in bleeding.

IMPORTANT: Implement additional protective measures, when changing the dressing of patients with identified increased risk of bleeding.

NOTE: For more information related to bleeding in the context of negative pressure wound therapy, please refer to the **Bleeding** section of this document.

Disconnection from the VivanoTec/VivanoTec Pro unit

The decision on how long the patient can be disconnected from the VivanoTec/VivanoTec Pro unit is a clinical assessment that has to be made by the physician.

The time interval for safe interruption of the therapy is strongly dependent on the overall patient and wound status, as well as the exudate composition and amount of exudate extracted per time unit.

Long interruption might lead to exudate retention and local maceration effects, as well as to a blocked wound dressing due to coagulation effects within the foam matrix. Lack of the effective barrier between the wound and the non-sterile environment increases the risk of infection.

IMPORTANT: Do not leave the dressing with the VivanoTec/VivanoTec Pro unit being switched off for prolonged time periods. In case of leaving the dressing for a longer time period, it is recommended for a physician to perform an evaluation of the wound condition along with the overall health status of the patient. According to the physician's evaluation, either rinsing of the wound along with a dressing change, or switching to an alternative therapy is recommended.

Intermittent pressure mode

Intermittent pressure, as compared to continuous pressure, may be used for enhancing local perfusion and granulation formation, if tolerable for the patient, the patient's health and wound condition. However, the continuous therapy is generally recommended for the treatment of patients with an increased risk of bleeding, acute enteric fistulas, highly exuding wounds, or when wound bed stabilization is required.

Pressure settings

PRECAUTION: Pressure settings below 50 mmHg may potentially lead to exudate retention and decreased therapeutic efficiency.

PRECAUTION: High pressure settings may increase the risk for micro-trauma, haematoma and bleeding, local hyperfusion, tissue damage or fistula formation.

Correct pressure setting for the Vivano negative pressure wound therapy must be decided by the physician and should be based upon the exudate output, overall patient status as well as recommendations from therapeutic guidelines.

Dressing on intact skin

Dressing on intact skin should cover about 5 cm around the wound. Prolonged or repeated dressing on bigger areas might result in tissue irritation.

IMPORTANT: In case of tissue irritation, discontinue the Vivano negative pressure wound therapy. Application of the wound dressing on intact skin may create wrinkles on the dressing surface. Formation of wrinkles significantly increases the risk of dressing leakage, and as a consequence occurrence of infection.

IMPORTANT: Special care should be taken when applying the wound dressing to the peri-wound area's fragile skin.

Dressing of wounds susceptible to irritation

For wounds susceptible to constant irritation (close proximity to limbs), continuous (rather than intermittent) therapy is indicated.

Circumferential dressings

Circumferential dressings should be used under medical supervision. Lack of adequate protective measures might cause local hypoperfusion.

Dressings in the vicinity of the vagus nerve

Dressing application in the vicinity of the vagus nerve should be performed under medical supervision, as its stimulation might cause bradycardia.

Allergies

Application of the Vivano negative pressure wound therapy is not recommended if the patient is allergic to any component of the Vivano System.

Electromagnetic fields

The VivanoTec Pro negative pressure unit must not be used in the presence of strong magnetic fields (such as an induction stove) and not used near application of HF surgical equipment.

Exudate canister application

The activities in this chapter may only be carried out by a physician or a qualified person in accordance with the instructions of a physician:

Before use, check that the exudate canister and the connection tubing are not damaged.

IMPORTANT: The exudate canisters of the negative pressure unit are sterile components and can therefore be used in a sterile surgical environment.

Connecting the exudate canister to the negative pressure unit

- Carefully remove the exudate canister from the sterile peel package.

IMPORTANT: The attached tubing piece must not be allowed to fall onto a non-sterile surface.

- Insert the exudate canister at a slight angle into the guidance on the negative pressure unit **1**.
- Tilt the exudate canister in the direction of the negative pressure unit until it fully engages with the blue unlock **2** key.
- Pull gently on the exudate canister to make sure that it is securely attached to the negative pressure unit.


IMPORTANT: The negative pressure unit including exudate canister must always remain in an upright position during operation.

IMPORTANT: The negative pressure unit including exudate canister must not be placed on a patient's bed.

Connecting the exudate canister to the wound dressing

- Make sure that the tube connectors are always correctly connected to each other in order to prevent malfunction.
- For application of the wound dressing kit, follow the instructions for use of the wound dressing kit.

Attach the connectors (tube ends) of the exudate canister to the connectors of the wound dressing kit **1**.

IMPORTANT: Check that the dressing collapses upon application of negative pressure.

Monitoring

IMPORTANT: The monitoring frequency must be adapted according to the overall health status of the patient and the treated wound condition, evaluated by the physician.

Regularly monitor the patient, the unit and the wound dressing. Watch out for wound exudate, maceration, infection and loss of vacuum. In order to ensure a safe therapy, it is necessary to check the wound dressing frequently. In doing so, be sure to check the wound dressing for imperviousness, and also for negative pressure, check the wound edges for maceration and check the wound edges and the exudate for signs of infection. In the event of signs of an infection, the physician must be informed immediately.

Avoid the risk of the tube becoming displaced or blocked by regularly checking the tubing system and its connections for leaks and kinks.

In order to guarantee uninterrupted operation of the negative pressure unit, the filling level of the exudate canister should be regularly checked. Once the maximum filling level has been reached, it should be replaced with a fresh exudate canister.

IMPORTANT: If there are any signs of infection, the physician must be notified immediately.


Exudate canister change

Removing the exudate canister


IMPORTANT: The exudate canister has to be changed after each patient, but at least one time per week.

Removing the exudate canister from the wound dressing


Press the unlocking device at the side of the connector and keep it pressed **2**.

Pull the two ends apart **3**.

Break off the stopper from the connector and plug it into the exudate lumen **4**. This ensures that no exudate can escape from the canister.

Removing the exudate canister from the negative pressure unit

- Press the blue unlock key **2** on the negative pressure unit.
- Slightly tilt the exudate canister and remove it.
- Properly dispose of the exudate canister. Observe local regulations.

IMPORTANT: Ensure correct disposal of the exudate canister, observing the locally applicable regulations.

Special notes

Keep out of the reach of children.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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