HARTMANN

Vivano® Med

Abdominal Kit

Wundverbandkit für die Unterdruck-Wundtherapie
Wound dressing kit for negative pressure wound therapy
Kit de pansements pour la thérapie par pression négative
Wondverbandkit voor wondtherapie met negatieve druk
Kit di medicazione per terapia con pressione negativa
Kit de apósitos para la terapia de presión negativa en heridas
Kit de pensos para a terapia de feridas por pressão negativa
Kit de curativos para a terapia de feridas por pressão negativa
Κιτ επιθέματος τραυμάτων για τη θεραπεία τραυμάτων με αρνητική πίεση
Sada krytí na rány pro podtlakovou terapii ran
Súprava krytia na rany na podtlakovú terapiu
Zestaw opatrunkowy do podciśnieniowej terapii ran
Kötszerkészlet negatívnyomás-terápiás sebkezeléshez
Набор раневой повязки для терапии ран отрицательным давлением
Превързочен комплект за терапия на рани с негативно налягане
Komplet potrošnog materijala za terapiju rana negativnim tlakom
Set de pansamente pentru terapia plăgilor cu presiune negativă
Komplet oblog za terapijo z negativnim tlakom
用於負壓傷口治療的傷口敷料套裝

DE Deutsch	3
EN English	11
FR Français	18
NL Nederlands	26
IT Italiano.....	34
ES Español.....	42
PT Português (PT)	50
BR Português (BR).....	58
GR Ελληνικά.....	66
CZ Česky	74
SK Slovensky	81
PL Polski.....	88
HU Magyar	96
RU Русский	103
BG Български.....	111
HR Hrvatski.....	119
RO Română.....	126
SI Slovenščina.....	134
HK 中文.....	141

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Stand der Information: 2022-04-12

Intended use

Foam dressing kit for negative pressure wound therapy intended to transfer negative pressure to the wound bed and to transport exudates from the wound. The integrated Organ Protection Layer is intended to protect the foam dressing against adhesion to organs and the organs against adhesion to the abdominal wall, additionally facilitating flow of wound exudate.

Individual components of VivanoMed Abdominal Kit are only intended for use in conjunction with the Vivano System by PAUL HARTMANN AG or the ATMOS S042 NPWT negative pressure wound therapy system.

VivanoMed Abdominal Kit is meant only for use in humans.

VivanoMed Abdominal Kit may only be used in acute care surgical environments enabling aseptic working conditions. The abdominal dressing will most often be applied in the operating theatre.

IMPORTANT: VivanoMed Abdominal Kit may only be used by a physician or a qualified person, as per the law of your country, in accordance with the physician's instructions.

Kit components identification

VivanoMed Foam

Hydrofilm
Transparent film
dressing

VivanoTec Port

VivanoMed Abdominal Organ
Protection Layer

Indications

VivanoMed Abdominal Kit is indicated for the temporary treatment of the open abdomen where closure of the abdominal wall is not possible or repeated access is required. This involves abdominal wounds with exposed visceral tissue, including abdominal wounds which can occur in the treatment of an abdominal compartment syndrome.

Contraindications

Contraindications for the use of the Vivano System:

- malignant tumour wounds
- non-enteric/unexplored fistulas
- untreated osteomyelitis
- necrotic tissue

IMPORTANT: Always protect vital structures with the Abdominal Organ Protection Layer. VivanoMed Foam must not be in direct contact with exposed nerves, anastomosis points, blood vessels or organs.

Using VivanoMed Abdominal Kit is not authorized in patients diagnosed with metastatic disease and bleeding diathesis.

NOTE: For more information on a particular contraindication, please refer to the **Warnings** and **Precautions** sections of this document.

Warnings

Please pay attention to the following warnings related to the use of VivanoMed Abdominal Kit.

Bleeding

NOTE: The Vivano System was not developed for the prevention or stopping of bleeding.

IMPORTANT: In the case of blood occurring suddenly or more frequently on the dressing, in the tubes or in the exudate container switch the negative pressure wound therapy unit off immediately, undertake haemostatic measures and inform the physician.

NOTE: Regardless of using the negative pressure wound therapy, certain medical conditions favor occurrence of bleeding complications.

The following circumstances increase the risk of possibly fatal bleeding, if not controlled with appropriate care:

- surgical sutures and/or anastomoses
- non-sutured haemostatic agents, e.g. spray wound sealant or bone wax
- trauma
- irradiation
- inadequate haemostasis
- wound infection
- treatment with anticoagulants or coagulation inhibitors
- protruding bone fragments or sharp edges

Patients with increased risk of bleeding complications should be monitored with additional level of care, under the responsibility of the physician.

IMPORTANT: In patients with diagnosed acute bleeding, coagulation disorders or being treated with anticoagulants, the **800 ml canister** should not be used for the exudate collection. Instead, a 300 ml canister should be utilized. Such a practice enables more frequent monitoring of the patient by the healthcare professionals, thereby reducing the potential risk of excessive blood loss.

IMPORTANT: When using non-sutured haemostatic agents, additional protective measures should be implemented, to prevent them from accidental displacement. Suitability for negative pressure wound therapy should be assessed by a physician on individual basis.

Monitoring of intra-abdominal pressure (IAP)

Due to possible risk of intra-abdominal hypertension (IAH) or abdominal compartment syndrome (ACS) occurrence, IAP of the patient should be closely monitored. Frequency of IAP inspections should be adapted to patient's overall condition, in accordance with the medical instructions of the physician.

IMPORTANT: In case of IAH or ACS occurrence, note the patient's IAP and disconnect the patient from the VivanoTec/VivanoTec Pro unit. If IAH or ACS persists without the negative pressure applied, discontinue the Vivano negative pressure wound therapy and address the medical condition according to the instructions of the physician. If IAH or ACS improves without the negative pressure applied, consider continuation of the Vivano negative pressure wound therapy with additional protective measures, or consider choosing an alternative approach. The time interval for safe interruption of the negative pressure wound therapy shall be based upon the evaluation of the physician, dependent on the overall wound and medical status of the patient. For more information regarding interruption of the negative pressure wound therapy please refer to the **Disconnection from the VivanoTec/VivanoTec Pro unit** Section of this document.

Malignant tumour wounds

Negative pressure wound therapy for malignant tumour wounds is contraindicated, as it is linked to the risk of enhanced tumour formation by proliferation support effect. However, it is considered as legitimate in a palliative context. For patients in the end-of-life stage for whom a complete cure is no longer the aim, the improvement of their quality of life by controlling the three most disabling elements: the odour, exudate, and pain associated with changing the dressings, outweighs the risk of accelerating the spread of tumours.

Non-enteric/unexplored fistulas

Application of the wound dressing on non-enteric or unexplored fistulas is contraindicated, as it may damage the intestinal structures and/or organs.

Untreated osteomyelitis

Application of the wound dressing in wounds with untreated osteomyelitis is contraindicated, as it may result in the spreading of the infection.

Necrotic tissue

Application of the wound dressing on necrotic tissue is contraindicated, as it may lead to local spreading of the infection.

Application of VivanoMed Foam on nerves, anastomosis points, blood vessels or organs

VivanoMed Foam must not be applied directly over exposed nerves, anastomosis points, blood vessels or abdominal organs as it may result in deterioration of the underlying structures.

Special precautions

Please pay attention to the following precautions.

Infected wounds

Wound dressings should be changed at regular intervals, according to instructions listed in the **Dressing change** section of this document. Infected abdominal wounds must be monitored more frequently, and may require wound dressings to be changed more often.

NOTE: For more information on abdominal wound monitoring in the context of negative pressure wound therapy, please refer to the **Monitoring** section of this document.

Typical signs of wound infection are redness, swelling, itching, increased warmth of the wound itself or in the neighbouring area, foul odour, etc.

Infected wounds may trigger systemic infection, manifesting by high fever, headache, dizziness, nausea, vomiting, diarrhoea, disorientation, erythroderma etc.

The consequences of systemic infection can be fatal.

IMPORTANT: If there is any suspicion of either local or systemic infection, contact the physician and consult if negative pressure wound therapy should be stopped, or an alternative therapy should be considered.

Blood vessels and organs

Blood vessels and organs should be adequately protected by means of fascias, tissues or other types of protective layers placed above them.

IMPORTANT: Special precautions need to be taken, when dealing with infected, weakened, irradiated or sutured blood vessels or body organs.

Bone fragments or sharp edges

Protruding bone fragments and sharp edges should be removed, or covered accordingly, before using VivanoMed Abdominal Kit, as they may damage blood vessels or abdominal organs and cause bleeding.

NOTE: For more information related to bleeding in the context of negative pressure wound therapy, please refer to the **Bleeding** section of this document.

IMPORTANT: Due to possible shifting in the relative position of tissues, vessels or abdominal organs within the wound, probability of contact with bone fragments, or sharp edges may be increased.

Enteric fistulas

In case of treatment of wounds containing explored enteric fistulas, additional level of precautions need to be implemented, if negative pressure wound therapy is to be applied. Presence of enteric fistula in close proximity to the wound increases the risk of wound contamination and/or infection. In order to mitigate the risk associated with the potential contact of intestinal content with the wound, enteric fistula needs to be surgically separated, following local guidelines or established surgical practices.

Spinal cord injuries with the development of autonomic hyperreflexia

Discontinue the negative pressure wound therapy if the patient has spinal cord injuries with the development of autonomic hyperreflexia.

Magnetic resonance imaging (MRI)

This device is not considered as MRI safe and must not be used in close proximity to a MRI unit.

Defibrillation

The VivanoTec/VivanoTec Pro unit must be disconnected if resuscitation of the patient using a defibrillator is necessary.

Hyperbaric oxygen therapy (HBO)

The VivanoTec/VivanoTec Pro unit must be disconnected for patients undergoing hyperbaric oxygen therapy, as its use is considered a potential fire hazard.

General precautions

Please pay attention to the following precautions.

Damaged, expired or contaminated product

Do not use any of the components of the kit, in case of damage, expiration or any suspicion of contamination. It may cause overall decline in the therapeutic efficiency, wound contamination and/or infection.

IMPORTANT: Before opening the product packaging, inspect it carefully for signs of damage. In particular, examine the integrity of the packaging surfaces and the entire seal area completely under suitable lighting conditions (white light or daylight). Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Single use only

All disposable components of the VivanoMed Abdominal Kit are meant for single use only. Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Re-sterilization

Components of VivanoMed Abdominal Kit are provided sterile. Do not re-sterilize any of the components as it may cause overall decline in the therapeutic efficiency of the kit and potentially lead to wound contamination and/or infection.

Product disposal

To minimize the risk of potential infection hazards or environmental pollution, disposable components of VivanoMed Abdominal Kit should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Safety measures for infection prevention

Implement and apply adequate personal protection and institutional infection control measures, when handling the components of VivanoMed Abdominal Kit (e.g. the use of sterile gloves, masks, gowns etc.).

IMPORTANT: Before and after using the stopper on the connector, it must be cleaned and disinfected.

Patient population

No general restrictions for the usage of VivanoMed Abdominal Kit in a different patient population (e.g. adults and/or children) is given. However, the VivanoTec/VivanoTec Pro unit has not been evaluated for use in paediatrics.

IMPORTANT: Before prescribing its use on a child, weight and height along with the overall health condition must be medically evaluated first.

Patient's health status

Weight and general condition of the patient should be taken into consideration during any application of negative pressure wound therapy.

Dressing size

The size of the dressing must be adapted to the size of the wound being treated with the negative pressure wound therapy. The Abdominal Organ Protection Layer must cover all exposed viscera and completely separate the viscera from contact with the abdominal wall.

Improper wound dressing size may either cause maceration and disintegration of the peri-wound tissue, or wound margins drying out and inefficient transfer of exudate.

NOTE: For more information on complications related to excessive coverage of the intact skin, please refer to the **Dressing on intact skin** section of this document.

Improper Abdominal Organ Protection Layer size may damage abdominal tissue/organs due to either high pressure applied to the organs, or dressing foam adhering to the adjacent structures.

IMPORTANT: In order to provide optimal conditions for negative pressure wound therapy, the film dressing should cover about 5 cm of intact skin around the wound.

Dressing placement

Use only dressings and Abdominal Organ Protection Layers coming directly from sterile packages.

Make sure that the Abdominal Organ Protection Layer completely separates the visceral tissue from the foam, as their prolonged contact may lead to fistula development and organ adhesion.

Do not force the placement of the foam, as it may lead to direct tissue damage or subsequent wound healing delay or even local necrosis, due to an elevated level of compression.

IMPORTANT: Always record the number of foams used for each wound.

The number of layers of film in the dressing can be adapted to each medical condition. Placement of multiple layers of film increases the risk of tissue maceration and consequently tissue irritation.

NOTE: Application of more than one Abdominal Organ Protection Layer to the abdominal wound, significantly impairs performance of the therapy and increases the risk of wound aggravation.

IMPORTANT: In case of tissue irritation, due to the use of multiple layers of the film, discontinue the Vivano negative pressure wound therapy.

Dressing removal

IMPORTANT: Always record the number of foams removed from the wound in order to guarantee the removal of all foams introduced.

Foam and/or Abdominal Organ Protection Layer left in the wound for a longer period of time than indicated in the **Dressing change** section, may cause growth of granulation tissue into the aforementioned components. This may increase the difficulty of dressing change and may promote wound infection, occurrence of bleeding, organ adhesion, damage or rupture of visceral tissue.

Dressing changes may possibly lead to disruption of the new granulation tissue, which may result in bleeding.

IMPORTANT: Implement additional protective measures, when changing the dressing of patients with identified increased risk of bleeding.

NOTE: For more information related to bleeding in the context of negative pressure wound therapy, please refer to the **Bleeding** section of this document.

Disconnection from the VivanoTec/VivanoTec Pro unit

The decision on how long the patient can be disconnected from the VivanoTec/VivanoTec Pro unit is a clinical assessment that has to be made by the physician.

The time interval for safe interruption of the therapy is strongly depended on the overall patient and wound status, as well as the exudate composition and amount of exudate extracted per time unit.

Long interruption might lead to exudate retention and local maceration effects, as well as to a blocked wound dressing due to coagulation effects within the foam matrix. Lack of the effective barrier between the wound and the non-sterile environment increases the risk of infection, as well as promotes occurrence of ACS, due to potential increase in IAP values.

IMPORTANT: Do not leave the dressing with the VivanoTec/VivanoTec Pro unit being switched off for prolonged time periods. In case of leaving the dressing for a longer time period, it is recommended for a physician to perform an evaluation of the wound condition along with the overall health status of the patient. According to the physician's evaluation, either rinsing of the wound along with a dressing change, or switching to an alternative therapy is recommended.

Pressure settings

IMPORTANT: Pressure applied to the abdominal wound should be adapted according to the medical evaluation.

Pressure settings below 50 mmHg may potentially lead to exudate retention and decreased therapeutic efficiency.

High pressure settings may increase the risk for micro-trauma, haematoma and bleeding, local hyperfusion, tissue damage or fistula formation.

Correct pressure setting for the Vivano negative pressure wound therapy must be decided by the physician and should be based upon the exudate output, values of measured IAP, overall patient status as well as recommendations from therapeutic guidelines.

NOTE: Intermittent pressure is not recommended for the application in the treatment of abdominal wounds.

Dressing on intact skin

Dressing on intact skin should cover about 5 cm around the wound. Prolonged or repeated dressing on bigger areas might result in tissue irritation.

IMPORTANT: In case of tissue irritation, discontinue the Vivano negative pressure wound therapy. Application of the wound dressing on intact skin may create wrinkles on the dressing surface. Formation of wrinkles significantly increases the risk of dressing leakage, and as a consequence occurrence of infection.

IMPORTANT: Special care should be taken when applying the wound dressing to the peri-wound area's fragile skin.

Dressings in the vicinity of the vagus nerve

Dressings in the vicinity of the vagus nerve should be performed under medical supervision, as its stimulation might cause bradycardia.

Allergies

Application of the Vivano negative pressure wound therapy is not recommended if the patient is allergic to any component of VivanoMed Abdominal Kit.

IMPORTANT: Hydrofilm contains an acrylic adhesive coating, which might present a risk of adverse reactions in patients who are allergic, or hypersensitive to acrylic adhesives.

Dressing application

Adequate personal protection and institutional infection control measures must be implemented, when applying the dressing.


Always apply the components of VivanoMed Abdominal Kit in the following order.

Wound preparation and application of Abdominal Organ Protection Layer

1. Before applying the dressing for the first time, and after each dressing change, the wound must be thoroughly cleaned and conditioned in accordance with the physician's instructions.
2. Open the sterile packaging of VivanoMed Abdominal Kit by opening the peel package.
3. Remove and open the packaging of the VivanoMed Abdominal Organ Protection Layer.
4. Unfold the VivanoMed Abdominal Organ Protection Layer and lay it out, so that application pockets are pointing towards you.


IMPORTANT: The smooth side should always face the abdominal opening. The side with the application pockets indicates the rough side of the layer which must point to the user.

5. Carefully apply the VivanoMed Abdominal Organ Protection Layer with the smooth side directly over the open abdominal cavity.

NOTE: The surplus margin of the VivanoMed Abdominal Organ Protection Layer may be either introduced into the abdominal wound or cut to match the size of the abdominal wound.

NOTE: Size and shape of the Abdominal Organ Protection Layer may be additionally adjusted by cutting, if it needs to be applied around tubes, drains or the falciform ligament.


IMPORTANT: When cutting the VivanoMed Abdominal Organ Protection Layer make sure that:

- application pockets are not entirely removed, as they facilitate the discrimination between the smooth and rough side of the Abdominal Organ Protection Layer
 - none of the residues fall into the wound
 - the size of the layer is sufficient to entirely cover internal organs
6. Insert your hands (or appropriate instruments) into the application pockets of the Abdominal Organ Protection Layer, in order to facilitate its placement between abdominal wall and internal organs.

IMPORTANT: Pay attention when using instruments instead of hands, to not break the application pockets e.g. due to excessive leverage.


7. Start placement of the VivanoMed Abdominal Organ Protection Layer by pushing it between the abdominal wall and the internal organs. Systematically work around the wound opening until all the internal organs are completely covered. Avoid any reinforced pushing as it may damage visceral tissue and sensitive structures.

Application of the dressing (Foam with Hydrofilm)

IMPORTANT: Before applying the dressing make sure that the VivanoMed Abdominal Organ Protection Layer is in place.

8. Match the size and shape of VivanoMed Foam to the wound opening by cutting with sterile scissors or a scalpel. The perforations may be used as a guidance to determine size and shape.


IMPORTANT: In order to minimize the risk of foam pieces falling into the wound, do not cut the foam in the wound vicinity.

NOTE: Make sure that the size of VivanoMed Foam matches the size of the VivanoMed Abdominal Organ Protection Layer covering the abdominal wound, and that it can be brought into contact with all remaining wound margins to provide tight wound closure.

9. Gently introduce the foam material into the wound.

NOTE: If necessary, multiple foam pieces may be applied on top of one another.

IMPORTANT: Make sure that no small foam pieces detach and contaminate the wound, when foam is being introduced into the wound.


IMPORTANT: Always record the number of dressing components used for each wound.

10. In order to seal the wound filled with VivanoMed Foam, the Hydrofilm film dressing should cover about 5 cm of intact skin around the wound. If necessary, the film dressing can be cut to shape.

IMPORTANT: Make sure that all release liners are removed from the film dressing, as leaving them might result in lower breathability of the film dressing.

IMPORTANT: Do not use any moisturizing/lipid transferring skin disinfection solution, before application of the film, as it may impair its adhesion to skin.


11. Apply the film dressing in accordance with the printed instructions and mould it to shape in order to achieve an airtight cover.

IMPORTANT: Do not stretch the film dressing during application as it may cause skin blisters and/or shear damage in attached skin areas.


IMPORTANT: Make sure you do not confuse the order of usage of the components of VivanoMed Abdominal Kit. Do not use the Abdominal Organ Protection Layer for wound sealing purpose.

12. Select a point in the centre of the film dressing, where VivanoTec Port is to be applied and cut a hole of approx. 2 x 4 cm. In doing so, make sure that the foam material is not perforated.


13. Apply VivanoTec Port so that the centre of the port is precisely above the prepared hole in the film dressing.

IMPORTANT: Make sure that the port holes are perfectly aligned with the prepared hole in the film dressing to guarantee optimal functionality of the system.

14. Connect the drainage tube to the canister tube of the Vivano System and commence the therapy in accordance with the user manual.

IMPORTANT: Improper connection of the canister to the Vivano System may result in impaired negative pressure performance, which may prolong the overall wound healing process.

IMPORTANT: Check that dressing collapses upon application of negative pressure.

Monitoring

IMPORTANT: The monitoring frequency must be adapted according to the overall patient's health status and the treated wound condition, evaluated by the physician.

When treating the wound with negative pressure using VivanoMed Abdominal Kit it should be periodically verified, if the nominal exudate flow/volume corresponds with the actual flow/volume. Additionally, the patient's health status and the wound condition should be closely monitored. Health status evaluation should include frequent IAP measurements to mitigate the risk of ACS occurrence, as well as cardiac output assessment. The wound condition should be monitored via checking of the wound dressing for signs of leakage, and the wound margins, as well as the exudate for signs of infection.

IMPORTANT: If there are any signs of infection, the physician must be notified immediately.

Dressing change

Dressing change must occur at regular time intervals, usually taking place between 48 to 72 hours, in order to provide optimal conditions for the wound healing process, as well as mitigate potential risk associated with prolonged presence of the dressing components within the abdominal cavity. In case of infected wounds, the dressing must be changed at intervals of 12 to 24 hours (or even more frequently if necessary).

IMPORTANT: It is a prerogative of the physician to decide, what frequency of wound monitoring and dressing change should be applied in a particular case. It is based upon the initial abdominal wound complexity, medical evaluation of the wound site, as well as the overall patient's health status. The wound dressing change interval should be increased accordingly, if the likelihood of organ protection layer adhesion seems high.

NOTE: For more information regarding monitoring of infected wounds, please refer to the **Infected wounds** section of this document.

Dressing change should be performed in the following order:

IMPORTANT: Ensure any bleeding has been stopped prior to dressing change. The physician must decide, if negative pressure wound therapy may be safely continued or an alternative therapy should be applied.

IMPORTANT: Each time the dressing is changed, all components of the dressing must be replaced with the components from the sterile package of VivanoMed Abdominal Kit.

1. Each time the dressing is changed the wound must be thoroughly cleansed and conditioned in accordance with the physician's instructions and applicable hygiene standards.
2. Inspect the wound and make sure that all pieces of VivanoMed Foam and the Abdominal Organ Protection Layer have been removed.
3. Apply the new dressing according to the instructions listed in the **Dressing application** section of this document.

NOTE: In case of adherence of the dressing to the wound, saline solution can be added to the foam material to enable dressing removal. After 15 to 30 minutes, the dressing should be carefully removed from the wound.

If the patient complains of pain while the dressing is being changed, premedication or the administration of a local anaesthetic should be considered.

Special notes

Keep out of the reach of children.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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