**HARTMANN**

# DermaPlast®

Ampoules & petites plaies  
Blaren & kleine wondjes  
Blasen & kleine Wunden  
Vesciche & piccole ferite  
Blisters & minor wounds

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## EN Instructions for use

### Product description

The hydrocolloid plaster is a non-sterile, single-use, self-adhesive plaster with a hydrocolloid formulation and occlusive properties (moist wound treatment). The product is meant for the application in humans by lay users.

### Composition

The plaster consists of one layer of adhesive hydrocolloid material (containing carboxymethyl-cellulose). The silicone paper protects the coated side of the plaster.

### Properties

Whilst traditional plasters work via dry wound healing, which allows the covered tissue to dehydrate and dry out, the hydrocolloid plaster is based on the mode of action of moist wound healing (occlusive properties).

Self-adhesive hydrocolloid plasters contain a layer which absorbs exudate and expands into a gel-like mass, thus providing cushioning for the treated damaged area and helping to relieve pressure and friction.

### Intended purpose

The hydrocolloid plaster is a single-use self-adhesive plaster for the protection and treatment of blisters and minor injuries, and for the absorption of exudate.

### Contraindications

The hydrocolloid plaster must not be used on untreated infected blisters or infected minor injuries, major injuries, highly bleeding or exudative wounds, for pressure ulcers or any other use apart from those mentioned in the intended purpose section.

### Side effects

Side effects may occur in patients who are allergic or hypersensitive to any of the product's components.

### Product disposal

The disposable components and used plasters can be disposed of normally with household waste.

### Mode of application

- Wash hands and ensure that the affected area is clean, dry and oil-free before applying the plaster.
- Cut the plaster to the required size.
- Warm the plaster with the hands before applying.
- Remove the first protective paper and apply the exposed part of the plaster, to the affected skin area.
- Then gently pull the second protective paper and apply the rest of the plaster. Gently mould the plaster edges to ensure optimum adhesion.
- The application maximum of the hydrocolloid plaster is up to 4 days per individual product. Should you wish to remove the plaster early or at the maximum time of application, first wet it briefly with warm water.

### Special precautions

- For single use only. Not reusable.
- Only apply the plaster if the hydrocolloid pad completely covers the blister or minor wound.
- Consult a doctor before applying the plaster on a potentially infected blister, minor wound, or if the blister or minor injury

shows signs of infection while plasters are in place (e.g. redness, overheating, swelling or pain).

- In cases of diabetes mellitus always consult a doctor before using plasters.
- Do not burst a sealed blister as this could risk an infection.
- Do not leave the product on the affected area for more than 4 days.
- Store plasters in a cool (not in the fridge), dry place, out of the reach of children.

### Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Date of revision of the text: 2021-12-02

 PAUL HARTMANN AG  
Paul-Hartmann-Straße 12  
89522 HEIDENHEIM, GERMANY  
www.hartmann.info

 Distributor

AT – PAUL HARTMANN Ges.m.b.H.  
2355 Wiener Neudorf  
BE – N.V. PAUL HARTMANN S.A.  
1480 Sint-Renelde  
CH – IVF HARTMANN AG · 8212 Neuhausen  
FR – Lab. PAUL HARTMANN S.à.r.l.  
Châtenois · 67607 Sélestat CEDEX  
GB – PAUL HARTMANN Ltd.  
Heywood/Greater Manchester OL10 2TT  
IT – PAUL HARTMANN S.p.A. · 37139 Verona  
NL – PAUL HARTMANN B.V. · 6546 BC Nijmegen

  
0 1 2 3

069 83777 (91022)