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Produktbeschreibung

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EN Instructions for use

Product description

Zetivut Plus Silicone Border is a sterile self-adhesive absorbent dressing with a silicone interface. The dressing comprises a perforated silicone wound contact layer, a semi-permeable polyurethane backing film and a green absorbent pad between these two layers. Acrylic adhesive is used to bind these layers together. The backing film is permeable to air but waterproof, which allows the patient to shower. The silicone layer facilitates the application of the dressing and promotes non-adherence to the wound. Above the silicone layer, the product features a white, hydrophilic non-woven fabric (viscose and polypropylene). The absorbent pad comprises cellulose and superabsorbent polyacrylate (SAP) for the absorption and retention of exudate and is wrapped in a hydrophilic cellulose tissue to distribute the liquid. It is covered with a hydrophobic green, polypropylene non-woven fabric and indicates the side facing away from the wound.

Properties and mode of action
Zetivut Plus Silicone Border absorbs and retains up to high levels of exudate in the absorbent pad. The wound dressing is able to redistribute pressure. The dressing pad and the dressing is able to redistribute pressure. It features a micro-adherent silicone interface and adhesive borders, so no additional materials are required to secure the dressing. The backing film is bacteria and waterproof. The dressing changes can be carried out atraumatically and almost painlessly.

Intended purpose

Zetivut Plus Silicone Border is a sterile single-use self-adhesive absorbent dressing with a silicone interface for long-term treatment of acute and chronic wounds, with moderate to high levels of exudate. It is used on adults by healthcare professionals in clinical or homecare environments and can be combined with cleansing solutions and local antiseptics, primary dressings and additional fixation. It can be used in context of compression therapy.

Indications
Zetivut Plus Silicone Border is suitable for the treatment of a wide range of wounds with moderate to high levels of exudate: acute wounds (burns, post-op, trauma wounds) and chronic wounds (venous ulcers, pressure ulcers, diabetic foot ulcers, tumour wounds).

Mode of application

- Clean the wound with a physiological solution or as instructed by a clinician.
- Select a suitable dressing size and shape for the wound, ensuring that the green absorbent pad exceeds the wound margins by approx. 1 to 2 cm. Do not cut the green absorbent pad. The transparent border can be cut to adapt the dressing to specific body areas.
- Before applying the dressing, after rinsing the wound, make sure that the surrounding skin is dry by using swabs to ensure that the dressing adheres.
- Remove all protective sheets. Place the adhesive side of the dressing onto the wound, so that the green side faces away from the wound.
- Press down the adhesive edges firmly and avoid wrinkles in the border to ensure the dressing stays securely in place.
- No additional fixation is required. If indicated, compression bandages can be applied as instructed by a clinician.
- The dressing can remain on the wound for up to 7 days, depending on the condition of the wound and the established clinical protocol.
- The dressing must be changed if clinically indicated or when exudate reaches the rim of the green absorbent pad.

Contraindications

Do not use Zetivut Plus Silicone Border in case of hypersensitivity to any of its components.

Special precautions

- Do not cut the green absorbent pad. The transparent border can be cut to adapt the dressing to specific body areas.
- Before treating a wound with healing impairment, a clinician has to be consulted to establish a clinical diagnosis. Treatment with Zetivut Plus Silicone Border cannot replace a causal treatment of the wound-healing impairment. If there are clinical signs of infection, the infection needs to be controlled with appropriate treatment before this dressing can be used.
- In all cases, follow the established clinical protocol.
- In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary, on these population groups this dressing should be used with caution and following a clinician's recommendation.
- It is not recommended to stand on the dressing in the case of wounds on the feet.
- Before opening the product packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.
- Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Side effects

Increased exudate, maceration, increased pain, critical colonization, adherence of the dressing to the wound, itching, burning, oedema, multiple ulceration, skin tears, blisters, erythema. If any of these side effects occur, consult a clinician.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices), if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with household waste.

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ES – Laboratorios HARTMANN S.A. - 08302 Mataró
IT – PAUL HARTMANN S.p.A. - 37139 Verona