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Produktbeschreibung

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EN Instructions for use

Product description

Omnistrip is a sterile wound closure strip for the closure of minor wounds and surgical incisions. It is intended to be used by healthcare professionals and lay users.

Composition

Omnistrip is made of discreet coloured, non-woven polyamide fabric that is coated with a polycarlate adhesive.

Intended purpose

Omnistrip wound closure strips are intended for primary closure of minor wounds healing by first intention (cuts and lacerations), surgical incisions and for support and stabilisation of surgical wounds. The device can be used on humans on breached and intact skin by healthcare professionals and lay-users.

Contraindications

Omnistrip is not suitable for the closure of wounds healing by secondary intention. It is also not suitable for infected, heavily soiled or bleeding wounds.

Special precautions
In the absence of evidence of any harmful effects and in view of the partly limited experience with the use of this product on sensitive population groups such as infants, children, pregnant or nursing women the product should be used with caution following a clinician's recommendation.

Side effects

Adverse reactions may occur in patients who are allergic or hypersensitive to any of its components.

Mode of application (see pictograms 1-4)

The skin has to be dry and clean before application.

- Peel off the middle section of the backing paper of the wound closure strips.
- Hold the edges of Omnistrip using the outer backing papers and place on the wound at right angles to the incision and press into place. When doing so, make sure that the wound closure strips are not subjected to stretching.
- Remove the external protective papers once at a time.
- Press the wound closure strips down again.

When removing Omnistrip, make sure to pull each end towards the centre of the wound and then lift up carefully.

A careful removal of the wound closure strips is advised especially on elderly or vulnerable skin.

Special instructions

It is also possible to apply the strips singly. Crosswise application is recommended on joints or other parts of the body where the edges of the incision are subject to stretching and tension. The device can remain in place for up to 7 days depending on clinical assessment. It can immediately be replaced by a new device up to 3 times.

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Before reusing the product's packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of Omnistrip should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the product together with regular (hospital) waste.

Incident reporting

For a patient/user (third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices), if during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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SE – HARTMANN ScandiCare AB - 33421 Anderstorp