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### Produktbeschreibung

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## 03 Instructions for use

### Product description

Zetivit Plus is a sterile superabsorbent dressing used on superficial, moderately to highly exuding wounds. The absorbent core absorbs and binds exudate.

### Composition

Zetivit Plus is a combined superabsorbent dressing which consists of four layers of different materials. On the wound side, the product features a soft, white, hydrophilic non-woven. The inner dressing core consists of soft cellulose fluff blended with liquid-retaining polyacrylate polymers. This absorbent core is enclosed in a thin nonwoven fabric that evenly distributes the liquid. On the side facing away from the wound, the product features a green layer of polypropylene nonwoven, which is water-repellent but permeable to air and allows for gas exchange.

### Properties and mode of action

Zetivit Plus rapidly absorbs exudate and binds it within the absorbent core. Exudate removal eliminates inhibitory factors from the wound, e.g. proteases. The increased absorption capacity of Zetivit Plus reduces the necessary frequency of dressing changes. This promotes wound rest and provides additional protection against contamination. Apart from its absorbent quality, Zetivit Plus also has a padding effect.

### Intended purpose

Single-use sterile superabsorbent dressing for long-term treatment of injured skin, acute and chronic, with moderate to high levels of exudate. It is used on adults only, by healthcare professionals in clinical or healthcare environments. It can be combined with local antiseptics, primary and secondary dressings and used under compression bandages.

### Indications

Zetivit Plus is suitable for the treatment of superficial, moderately to highly exuding wounds: acute wounds (traumatic wounds, post-operative wounds, lymphatic wounds) and chronic wounds (decubitus/pressure sores, venous, arterial or mixed leg ulcers, tumour wounds). Zetivit Plus can be used under compression bandages.

### Mode of application

- Select Zetivit Plus to match the wound size so that the dressing extends at least 1–2 centimetres beyond the wound margins.
- Place the white side of Zetivit Plus on the wound so that the green special nonwoven layer faces away from the wound.
- Secure the dressing e.g. with adhesive tape, conforming bandages or, if necessary, compression bandages.
- One dressing can remain on the wound for up to 7 days, depending on the condition of the wound. Change the dressing if clinically indicated or when exudate reaches the edges of the dressing or is noticeable through the green nonwoven top layer.
- In case of infected wounds, a combination with e.g. Atramaun Ag is possible.
- With decreasing amounts of exudate, the use of a suitable hydroactive wound dressing (e.g. HydroAct) is recommended.

### Contraindications

Do not use Zetivit Plus on dry wounds or on exposed bones, muscles or tendons. Do not use Zetivit Plus in case of hypersensitivity to any of its components.

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- U slučaju infekcijskih rana moguće je kombinirati i, primjerice, proizvodom Atrauman AG.