

Hydrofilm®



Spec.-No.: D 7.4021a

Department: SM-SPE

Date: 2007-10-08

1. General Product Description

Hydrofilm® is a self-adhesive film dressings made of semi-permeable, waterproof polyurethane film, transparent; protect against bacterial penetration; with hypoallergenic adhesive, can be removed painlessly without leaving any residues.

Hydrofilm® carries the CE mark according to EU directive 93/42/ EEC for medical devices. The product is classified as a class IIa medical device.

A conformity assessment has been performed for Hydrofilm® and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

The safe use of Hydrofilm® therefore, is ensured if the product is used in line with the intended purpose.

The semi permeable polyurethane film allows the free passage of water vapour and gases.

The hypoallergenic adhesive layer ensures good adhesion, even on conical areas of the body.

2. Application / Indication

Hydrofilm® as a protective layer over the wound to reduce the risk of secondary infection and protect against physical damage on dry, primary healing wounds and nearly healed epithelial wounds; to secure catheters or cannulae; ideal as a secondary dressing, e.g. to cover gel-forming calcium alginate dressings.

3. Contraindication

Hydrofilm® should not be used on clinically infected, bleeding or highly exudative wounds.

4. Presentations

Hydrofilm® sterile, individually sealed

Ref. No.	Dimensions	Number peel pouch per folding box
685 755	6 cm x 7 cm	P10
685 756	6 cm x 7 cm	P100
685 757	10 cm x 12,5 cm	P10
685 758	10 cm x 12,5 cm	P100
685 759	10 cm x 15 cm	P10
685 760	10 cm x 15 cm	P50
685 761	15 cm x 20 cm	P10
685 762	15 cm x 20 cm	P50
685 763	10 cm x 25 cm	P25
685 764	12 cm x 25 cm	P25
685 765	20 cm x 30 cm	P10

Hydrofilm®



Spec.-No.: D 7.4021a

Department: SM-SPE

Date: 2007-10-08

5. Product Characteristics

Composition:

- backing: semi-permeable, transparent polyurethane film
- adhesive: acrylic adhesive, hypoallergenic, colophonium-free
- handling aid two-piece supporting film (based on polyester)
- Moisture vapour transition rate (MVTR) according to DIN EN 13726-2: > 1.100 g/m² / 24h (typical value around 1.300 g/m² / 24h)
- Adaptiveness according to DIN EN 13726-4: 1,2 N/cm
- persistent deformation: 1 %

6. Product Requirements

- sterilization ethylenoxide, using a validated procedure (e.g.: ISO 11135)
- adhesive power > 4 N / in

7. Labelling

Lot-No. with 9-digit code:

e.g.:

7
year

XXX
for internal
purposes

12
week of production

XXX
for internal
purposes only

Expiry date:

e.g.:



2012
year

09
month

Shelf life: 5 years

Hydrofilm®



Spec.-No.: D 7.4021a

Department: SM-SPE

Date: 2007-10-08

8. Packaging

Separate sealed into peelpack out of sterile paper, in folding box; Transport carton acc. DIN, sealed with adhesive tapes; packed onto europallet.

Date: 2007-10-08

A. Bolz

- SM-PDW -