

# Foliodrape® Protect Stockinets

## General Product Description/Intended Purpose

Foliodrape® Protect Stockinets consist of a two layer material, a nonwoven material laminated with Polyethylene-film (PE), and are available in different sizes and for different applications. They are single-use sterile patient drapes and covers for short-term use for prevention of infection on humans in clinical environment..

Foliodrape® patient drapes and covers are intended to be used as sterile draping of the patient in operating theatres and all other surgical environments (e.g. day surgery, ambulance, doctor's surgery).

The main task of patient drapes and covers consists of the prevention of pathogen entry into the surgical wound. Consequently, the risks of infections and blood-borne diseases are reduced. The Foliodrape® products are placed on the market sterile and are in Medical Device Class I.

## Application/Indication

Draping is the perioperative process of covering a patient and surrounding areas in the operating theater with sterile drapes and covers to create and maintain a barrier, which protects the wound from the entry of microorganisms. The microorganisms can originate from the patient's skin, the OR, or the healthcare personnel. The aim of draping is to prevent the passage of pathogens between non-sterile and sterile areas. The different layers of hydrophilic nonwoven material laminated with Polyethylene-film act together as a fluid and bacteria barrier and minimize the transmission of micro-organisms. Depending on the procedure for which they are used, the Foliodrape® products usage time will vary from few minutes to several hours of continuous usage.

## Reference Numbers

Name of Product	Reference Number	Size of Product
Foliodrape® Protect Stockinet	258 980/3	35 x 120 cm
	258 981/3	35 x 80 cm
	258 982/3	25 x 80 cm

## Residual Risks, Contra-Indications and Undesirable Side Effects, Warnings

No contraindications known.

## Sterile Device

Do not use when sterile packaging or the product are damaged or were unintentionally opened before use.

## Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

## Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of Foliodrape® products should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. All infected waste is recommended to be disposed of by incineration.

Non-contaminated material such as waste paper, packaging and non-sterile but not biologically contaminated materials can be disposed separately. Separate disposal containers should be available where waste is created, so that staff can sort the waste as it is being discarded.

## Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.



# Foliodrape® Protect Stockinets

## Product Performance Characteristics

The Foliodrape® range is developed under high-quality standards to ensure safety and infection prevention. All products are free of chlorinated substances such as CMR/PBT/vPvB substances (carcinogenic, mutagenic and toxic to reproduction, etc.). There are no added flame-retardants (like polybrominated bi-phenyls, etc.).

## Material Characteristics

Composition	Details
<b>Nonwoven Material:</b>	Hydrophilic PP-spunbond, 40 g/m <sup>2</sup> , color: aqua
<b>Film:</b>	Pore-free polyethylene-film, 30 g/m <sup>2</sup> , color: aqua
<b>Total Laminate:</b>	Thermobonded, 70 g/m <sup>2</sup> , latex-free
<b>Sterilization:</b>	With ethylene oxide gas according EN 556-1 and EN 11135-1

1.  1. Hydrophilic PP spunbond
2.  2. PE-Film (impermeable for moisture and bacteria)

Foliodrape® PROTECT drapes corresponds to all requirements of European standard EN 13795-1:2019 for the highest performance level "High performance, critical area". For more detailed information, please visit page 5.

## Product Characteristics


Product Name	REF number	Size of drape	Additional Information
<b>Foliodrape® Protect Stockinet</b>	258 980/3	35 x 120 cm	Wrapped in crepe
	258 981/3	35 x 80 cm	Wrapped in crepe
	258 982/3	25 x 80 cm	Wrapped in crepe

## Labelling

The Foliodrape® labeling describes the intended purpose of the device as supported by clinical evaluation. Safety information supplied by the legal manufacturer on the labeling are based on the risk management documentation and elaborately present all warnings and precautions according the new requirements of the standards concerning the product.

### Lot-No. with 8-Digit Code

e.g.:

	0 year	12 Week of production	XXXXX for internal purposes only
-------------------------------------------------------------------------------------	-----------	-----------------------------	-------------------------------------

### Manufacturing Date

e.g.:  2020  
year      04  
month      07  
day


### Use-by-Date

e.g.:  2025  
year      04  
month      07  
day


**Shelf Life: 5 years**


# Foliodrape® Protect Stockinets


Medical Device **MD**

Unique Device Identification (UDI)  **UDI**

Single Sterile Barrier System  Single sterile barrier system

Do not use if the product sterile barrier system or its packaging is compromised  Do not use if package is damaged


Single use device 


Keep dry 


Sterilized using ethylene oxide **STERILE EO**


Legal manufacturer  PAUL HARTMANN AG  
Paul-Hartmann-Straße 12  
89522 HEIDENHEIM, GERMANY


European Conformity  0 1 2 3

Dispose of your waste in an appropriate manner 

Recycling symbols – Corrugated cardboard / Polyethylene Low-Density  PAP PE

Producer contributes to recycling of packaging 

Global trade identification number  **GTIN**

Indication arrow where to open the packaging 

Name and intended purpose of the product **Stockinette, Krepp**  
**Stockinet, crepe**

Reference / catalogue number **REF**

Pictogram with product look and dimensions 

## Storage and transport conditions

Transport and store in conditions for climatic regions I-IVa according to WHO Technical Report 953, 2009.

# Foliodrape® Protect Stockinets

## Logistic Data

Name of product	REF number	Peel Pouch size [mm]	Products per Peel Pouch	Peel Pouches per dispenser	Dispenser size	EAN of dispenser	Dispensers per carton <sup>1</sup>
<b>Foliodrape® Protect Stockinet</b>	258 980/3	Size: L 265 x 365	1	20	Size: L 275 x 374 x 234	4049500046060	2
	258 981/3	Size: M 130 x 290	1	30	Size: M 140 x 300 x 370	4049500046008	4
	258 982/3	Size: M 130 x 290	1	25	Size: M 140 x 300 x 370	4049500046039	4

Latest Date of Revision: 2020-09-07

# Foliodrape® Protect Stockinets

## Foliodrape® Protect stockinet, base material, testing results according EN 13795-1:2019

Property	Test method	Unit	Product requirement (PR) according EN 13795 -1:2019 for level "High Performance"	Test results <sup>1</sup>
Resistance to microbial penetration – Wet	EN ISO 22610	l <sub>8</sub>	Less critical product area: Not required Critical product area: 6.0	6.0 (no penetration)
Cleanliness – microbial	EN ISO 11737-1	CFU/ 100 cm <sup>2</sup>	Less critical product area: ≤ 300 Critical product area: ≤ 300	Not required <sup>3</sup>
Particle release	EN ISO 9073-10	log <sub>10</sub> (lint count)	Less critical product area: ≤ 4,0 Critical product area: ≤ 4,0	M <sub>d</sub> = 1.1 U <sub>q</sub> = 1.5
Resistance to liquid penetration <sup>4</sup>	EN 20811	cm (H <sub>2</sub> O)	Less critical product area: ≥ 10 Critical product area: ≥ 100	M <sub>d</sub> = 196 L <sub>q</sub> = 172
Bursting strength - Dry <sup>5</sup>	EN ISO 13938-1	kPa	Less critical product area: ≥ 40 Critical product area: ≥ 40	M <sub>d</sub> = 173 L <sub>q</sub> = 173
Bursting strength - Wet <sup>5</sup>	EN ISO 13938-1	kPa	Less critical product area: Not required Critical product area: ≥ 40	M <sub>d</sub> = 186 L <sub>q</sub> = 184
Tensile strength – Dry	EN 29073-3	N/50mm	Less critical product area: ≥ 20 Critical product area: ≥ 20	Longitudinal direction: M <sub>d</sub> = 93,9 L <sub>q</sub> = 89,8 Lateral direction : M <sub>d</sub> = 67,4 L <sub>q</sub> = 64,4
Tensile strength – Wet	EN 29073-3	N/50mm	Less critical product area: Not required Critical product area: ≥ 20	Longitudinal direction: M <sub>d</sub> = 96,6 L <sub>q</sub> = 96,4 Lateral direction: M <sub>d</sub> = 65,1 L <sub>q</sub> = 63,7

<sup>1</sup> Example REF T1615198 – Lot 699601007, Test report from independent Testing-Institute is available on request. Abbreviations: L<sub>q</sub> = lower quartile value, U<sub>q</sub> = upper quartile value and M<sub>d</sub> = median value. Product conforms when: L<sub>q</sub> ≥ PRmin resp. U<sub>q</sub> ≤ PRmax.

<sup>2</sup> Not required when the product fulfill the criteria for high performance in critical product area for resistance to microbial penetration wet (see Annex B7 in EN 13795:2019)

<sup>3</sup> see clause A1.2 in EN 13795-1:2019, Foliodrape® articles are sterile according to EN 556

<sup>4</sup> measured without use of supporting tissue (according to EN 20811),

<sup>5</sup> values after membrane correction (according to EN ISO 13938-1)