Veroval® compact

BPU 22

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Apreciado cliente:

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Upper arm blood pressure monitor

Dear customer

We are pleased that you have chosen to buy a blood pressure monitor from HARTMANN. The Veroyal® BPU 22 is a quality product for fully automatic blood pressure measurement on the upper arm of adults and is suitable for clinical and domestic use. Requiring no preliminary settings, this blood pressure monitor easily and automatically inflates for simple, quick and reliable measurement of the systolic and diastolic blood pressure and the pulse rate. The monitor also indicates if it detects an irregular heartbeat.

We wish you all the best for your health.

Please read these instructions for use carefully before using the device the first time because correctly measuring your blood pressure depends on appropriate handling of the device. These instructions for use are designed to teach you each step for self-measurement of blood pressure using the Veroyal® BPU 22. You will be given important and useful hints to ensure that you obtain reliable results for your personal blood pressure profile. Use this device in accordance with the information provided in the instructions for use. Please store these instructions carefully and ensure other users can access them. Check that the outside of the package is undamaged and that the contents are complete.

- Package contents:**
- Blood pressure monitor
 - Universal upper arm cuff
 - 4x 1.5 V AA batteries
 - Storage bag
 - Instructions for use with warranty certificate

1. Device and display description

Blood pressure monitor

- 1 Extra large LCD display
- 2 START/STOP button
- 3 Memory button for User 1
- 4 Memory button for User 2
- 5 Cuff connector socket
- 6 Battery compartment

- Cuff**
- 7 Cuff with application instructions
 - 8 Size indicator to correctly fit cuff
 - 9 Cuff plug
 - 10 Cuff tube

Display

- 11 Systolic blood pressure
- 12 Diastolic blood pressure
- 13 Pulse rate
- 14 Flashes when the device is measuring and the pulse rate is detected
- 15 Irregular heartbeat
- 16 Display to check the fit of the cuff
- 17 Battery symbol
- 18 Traffic light system for your values
- 19 Average value (A), morning (AM), evening (PM) / memory allocation number
- 20 User memory
- 21 Date and time display

2. Important information

Signs and symbols

The following symbols are used in the operating instructions, on the packaging and on the type plate of the device and accessories:

	0123	Labelling in accordance with Medical Devices Directive 93/42/EEC
		Manufacturer
		Authorised representative in the European Community
		Batch code
		Catalogue number
		Serial number
		Follow instructions for use
		Caution (Please note)
		Temperature limit
		Humidity limitation
		Protection against electric shock (type BF)
		Direct current
		Protected against solid objects with a diameter of ≥12.5 mm and against vertically falling drops of water
		Symbol for the marking of electrical and electronic equipment
		Dispose of packaging in an environmentally responsible way
		Cardboard recycling code

Important instructions for use

- Only use the device to take blood pressure measurements on the upper arm. Do not place the cuff on other areas of the body.
- Use only the supplied cuff or an original replacement cuff. The measured values may otherwise be incorrect.
- Only use the device on persons whose upper arm circumference is within the range indicated for the device.
- If measured values are questionable, repeat the measurement.

- Do not leave the device unattended near small children or persons who cannot operate it themselves. There is a risk of strangulation if the cuff tube becomes wrapped around the neck. Swallowing small parts that have detached from the device can also result in suffocation.
- Do not under any circumstances carry out blood pressure measurements on newborns, babies or toddlers.
- Please do not place the cuff over a wound as this may cause further injuries.
- Do not place the cuff on persons who have had a mastectomy.
- Please note that the pressure build-up in the cuff can lead to temporary disruption of medical devices being used at the same time on the same arm.
- Do not use the blood pressure monitor together with a high frequency surgical device.
- If an intravenous treatment is being carried out or a venous catheter is inserted in the arm, blood pressure measurements can lead to injuries. Never use the cuff on the arm to which these conditions apply.
- The function of the arm on which the cuff is placed may become impaired during inflation.
- If you are carrying out the measurement on another person, please ensure that using the blood pressure monitor does not have a lasting negative effect on the blood circulation.
- Taking measurements too often within a short time period as well as maintaining continuous cuff pressure can disrupt the blood circulation and cause injuries. Please take a break between measurements and do not bend or fold the air tube. If the monitor malfunctions, remove the cuff from the arm.
- Do not use the blood pressure monitor on patients with pre-eclampsia during pregnancy.

Important instructions for self-measurements

- Even minor changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, environmental factors) can lead to fluctuations in blood pressure. This explains why deviating values are often measured at the doctor or pharmacist.
- Measurement results depend fundamentally on the measurement location and the position of the patient (sitting, standing, reclining). They are also influenced by physical activity and physiological preconditions of the patient, for example. For comparable values, carry out the measurement in the same location and the same position.
- Diseases of the cardiovascular system can lead to erroneous readings or lower the accuracy of the measurement. The accuracy may also be affected if you have very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or tremors.

Consult your doctor before measuring your own blood pressure if ...

- you are pregnant. Blood pressure may vary during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because elevated values may affect foetal development. Check with your doctor in any case to find out if and when you should measure your own blood pressure, especially if you suffer from pre-eclampsia.
- you have diabetes, liver function disorders or narrowing of blood vessels (e.g. arteriosclerosis, peripheral arterial occlusive disease) because deviating measurements may occur in such cases.
- you have certain blood disorders (e.g. haemophilia) or severe circulatory disorders or are taking blood-thinning drugs.
- you wear a cardiac pacemaker because it may lead to deviating measured values. The blood pressure monitor itself has no impact on a cardiac pacemaker. Please note that the pulse rate displayed is not suitable for checking the frequency of cardiac pacemakers.
- you tend to bruise easily and/or are sensitive to pressure pain.
- you suffer from severe irregular heartbeat or cardiac arrhythmias. Because of the oscillometric measurement method that is used, in some cases incorrect measured values are obtained or the blood pressure cannot be measured.
- If this symbol appears often, it may indicate heart rhythm disorders (arrhythmias). If this is the case, please consult your doctor. Severe heart rhythm disorders may sometimes produce false measurements or lower the measurement accuracy. Please discuss with your doctor whether taking your own blood pressure readings is suitable for you.
- Your self-measured values are for your information only – they are not a substitute for a medical examination! Discuss your measured values with your doctor and under no circumstances make your own medical decisions based on these measurements (e.g. medications and their dosage)!
- Measuring your own blood pressure is not a substitute for medical treatment! Do not interpret your measured values by yourself and do not use them for self-prescribed treatment. Take measurements based on the instructions of your doctor and trust his or her diagnosis. Take medications as prescribed by your doctor and never change the dose on your own. Discuss the appropriate time to measure your blood pressure with your doctor.

Handling batteries

- Observe the polarity labels plus (+) and minus (-).
- Only use high-quality batteries (see specification in Chapter 11 'Technical data'). If you use low-quality batteries, we cannot guarantee the specified number of measurements.
- Never mix old and new batteries or batteries from different manufacturers.
- Remove empty batteries immediately.
- Replace batteries if the battery symbol remains illuminated.
- Always replace all the batteries at the same time.
- If the device remains unused for a longer period of time, you should remove the batteries to prevent possible leakage.

Batteries

Choking hazard

Small children could swallow batteries and suffocate on them. Keep batteries out of the reach of children!

Risk of explosion

- Do not throw batteries into a fire.
- Batteries must not be charged or short-circuited.
- If a battery has leaked, wear protective gloves and clean the battery compartment with a dry cloth. If liquid from a battery cell comes into contact with skin or eyes, clean the affected area with water and seek medical attention if necessary.
- Protect batteries from excessive heat.
- Do not disassemble, open or crush batteries.

Safety information concerning the device

- This blood pressure device is not waterproof!
- This blood pressure device is made of high-quality electronic precision components. The accuracy of the measured values and the lifetime of the device depend on careful handling.
- Protect the device from strong shocks, impacts or vibrations and do not let it drop on the floor.
- Do not excessively bend or kink the cuff or the air tube.
- Never open the device. The device must not be modified, dismantled or repaired by the user. Repairs may only be carried out by authorised specialists.
- Never inflate the cuff when it is not properly applied to the upper arm.
- Use the device only with the approved upper arm cuff. The device may otherwise sustain external or internal damage.
- The cuff tube may only be removed from the unit by pulling on the corresponding connector plug. Never pull on the tube itself!
- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight because this may cause it to malfunction.
- Keep the packaging, batteries and device out of the reach of children.
- Please comply with the storage and operating conditions defined in Chapter 11 'Technical data'. Storage or use outside the specified temperature and humidity range can affect the accuracy of the measurement or the function of the device.
- Do not use the device near strong electromagnetic fields and keep away from radio equipment or mobile phones. Portable and mobile high-frequency and communication devices, such as telephones and mobile phones, can impair the functionality of this electronic medical device.

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.

Instructions for the measurement function check

Every Veroyal® BPU has been carefully tested by HARTMANN for metrological accuracy and has been developed with a view to a long service life. We recommend carrying out a metrological check every 2 years for **devices in professional use**, for example, in pharmacies, medical practices or hospitals. Please also observe the national regulations specified by the legislator. Metrological checks may only be carried out by the competent authorities or authorised maintenance providers against cost reimbursement.

Disposal information

- To protect our environment, empty batteries must not be disposed of in household waste. Please comply with the relevant waste disposal regulations or use public collection points.
- This product is subject to the European Directive 2012/19/EU on Used Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please learn about your local regulations on the proper disposal of electrical and electronic products. Proper disposal protects the environment and human health.

3. Blood pressure information

To determine your blood pressure you need to measure two values:

- **SYS** – Systolic (upper) blood pressure, which is produced when the heart contracts and pumps blood into the blood vessels.
- **DIA** – Diastolic (lower) blood pressure, which is the value measured when the heart muscle is dilated and fills again with blood.
- Blood pressure values are displayed in mmHg.

So that you can understand your results more easily, a traffic-light system is available on the left side of the Veroyal® BPU 22 that directly indicates the result, making it easier to categorise the measured value. The World Health Organisation (WHO) and the European Society of Cardiology (ESC) & European Society of Hypertension (ESH) have developed the following summary for classifying blood pressure values:

Results indicator	Evaluation	Systolic pressure	Diastolic pressure
Red	Grade 3 hypertension	above 179 mmHg	and/or above 109 mmHg
Orange	Grade 2 hypertension	160 – 179 mmHg	and/or 100 – 109 mmHg
Yellow	Grade 1 hypertension	140 – 159 mmHg	and/or 90 – 99 mmHg
Green	High normal	130 – 139 mmHg	and/or 85 – 89 mmHg
Green	Normal	120 – 129 mmHg	and/or 80 – 84 mmHg
Green	Optimal	up to 119 mmHg	and up to 79 mmHg

Classification of Office Blood Pressure and definition of Hypertension Grade
(Source: 2018 ESC/ESH Guidelines)

4. Preparing the measurement

Inserting/changing the batteries

- Open the battery cover on the underside of the device. Insert the batteries (see Chapter 11 'Technical data'). Ensure correct polarity (+ and -) when inserting batteries. Close the battery lid. flashes on the display. Set date and time as described below.
- If the 'Change battery' symbol is permanently illuminated, blood pressure can no longer be measured and you must replace all the batteries.

Setting the date and time

- Be sure to set the date and time correctly. This is the only way to save your measured values correctly with date and time for subsequent retrieval.

- To switch to setting mode, reinsert the batteries or hold the START/ STOP button down for 5 seconds. Proceed as follows:

- Hours**
- The hour format flashes on the display.
- Select the desired hour format using the Save buttons / and confirm using the START / STOP button .

- Date**
- The display shows the year (a), the month (b) and the day (c).
- Depending on the display, use the Save buttons / to select the year, the month or the day and confirm using the START / STOP button .

- If the 12-hour format is selected, the month is displayed before the day.

- Time**
- The hour (d) and the minutes (e) flash consecutively on the display.
- Depending on the display, use the Save buttons / to select the current hour or minutes and confirm using the START / STOP button .

The device switches off automatically once all data has been set.

5. Measuring your blood pressure

Application of the cuff

- Before applying the cuff, insert the cuff socket into the connector socket on the left side of the device.
- Do not mechanically constrict, compress or bend the sleeve tube.
- The blood pressure must be measured on a bare upper arm. If the cuff is completely open, thread the end of the cuff through the metal ring to form a loop (see Fig. 1). Position the cuff on the upper arm so that the lower edge is 2-3 cm above the crook of the arm and above the artery (see Fig. 2). The tube points to the centre of the palm.


Fig. 1


Fig. 2

- Now bend your arm slightly, grasp the free end of the cuff, wrap it firmly round your arm and close the hook-and-loop fastener.

- The cuff should fit securely but not tightly. You should be able to insert two fingers between the arm and the cuff. Make sure the tube is not bent or damaged.

- Important: It is essential that the cuff is correctly applied to obtain correct measurements. The cuff included with the device is a universal cuff for upper arm circumferences of 22 cm to 42 cm. The white arrow must point to an area within the sizing scale. If it is outside the range, it is not possible to guarantee a correct measurement result.

- This innovative Veroyal® device with Comfort Air Technology ensures comfortable measurement. The individual inflation pressure is determined independently for each measurement and depends on the corresponding systolic blood pressure value. This allows for more comfortable upper arm measurements.

Carrying out the measurement

- Rest for approx. 5 minutes before measurement.
- You should take your blood pressure in a quiet place, in a relaxed and comfortable seated position.
- Measurement can be taken on the right or left arm. We recommend you perform the measurement on the left upper arm. Over the long term, the arm giving higher readings should be used for blood pressure monitoring. However, if there is a very clear difference between readings on either arm, you need to check with your doctor which arm you should use for the measurement.
- Always measure on the same arm and place your forearm in a relaxed position on a support.
- We recommend that you measure your blood pressure while sitting with your back supported by the back of the chair. Place both feet flat on the floor next to each other. Do not cross your legs. Relax your forearm and hand with your palm facing upwards on a support and make sure the cuff is at the height of your heart.
- Always take measurements at the same time of day. Only regular measurements taken at the same time each day over a longer period of time give a meaningful assessment of blood pressure values.
- Do not measure when you have a strong urge to urinate. A full bladder can increase blood pressure by about 10 mmHg.
- Do not take your blood pressure after taking a bath or playing sport.
- Do not eat, drink or exercise for at least 30 minutes before the measurement.
- Please wait for at least one minute between two measurements.
- Start a measurement only after applying the cuff. Press the START/ STOP button . The appearance of all display segments followed by the time and date indicates that the device is carrying out its self-test and is ready for use.
- Check the display segments for completion (see Chapter 1).
- After about 0.5 seconds the cuff will automatically inflate. If this inflation pressure is insufficient or if the measurement is interrupted, the device continues to pump in increments of 40 mmHg until a high enough pressure is reached. During inflation, the results indicator on the left also increases at the same time.

- If the cuff has been applied securely enough on the arm, the cuff symbol appears on the display. If the cuff symbol does not appear in the display, the cuff has not been fitted securely enough and the error message 'E3' appears on the display after a few seconds.

Important:

- As the air pressure in the cuff decreases, the heart symbol and the dropping pressure in the cuff are shown on the display.
- After the measurement is complete, the systolic and diastolic blood pressure as well as the pulse rate (see Fig. 3) appear simultaneously.

Fig. 3

- In addition to the measured values, the time, date, the associated user memory or as well as the associated memory numbers (e.g.) appear. The measured value is automatically stored in the displayed user memory. As long as the result of the measurement is displayed, you can assign values to the corresponding user memory by pressing the or button. If you do not assign them, the measured value is automatically stored in the displayed user memory. You can use the results indicator on the left of the display to classify your measurement result (see table in Chapter 3 'Blood pressure information').
- To switch off the device, press the START/STOP button , otherwise the device will automatically switch off after 1 minute.

If you wish to stop the measurement for any reason, simply press the START/STOP button . The cuff inflation or measurement is interrupted and the pressure is automatically released.

- If this symbol is shown in the display, the unit has detected an irregular heartbeat during the measurement. However, the measurement may also have been disrupted by body movement or speaking. It is best to repeat the measurement. If you see this symbol regularly when measuring your blood pressure, we recommend that your doctor checks your heart rhythm.

6. Memory function

User memory

- The Veroyal® BPU 22 stores up to 100 measurements in each user memory. Once all memory slots are occupied, the oldest value is deleted. Memory recall is activated by pressing the or button when the device is switched off. For values saved in the first user memory, press the button, for the second user memory press the button.

Average values

- After selecting the particular user memory, the corresponding symbol or and an A appear on the display. The average value of all the data stored for the particular user memory is displayed (see Fig. 4).


Fig. 4

- By pressing the button again (or button if you are in user memory 2), the average values for all morning measurements 'PM' (5 am to 9 am) taken on the last 7 days appear.

- By pressing the button again (or button if you are in user memory 2) the average values for all evening measurements 'PM' (6.00 pm to 8.00 pm) taken on the last 7 days appear.

Individual measured values

- By pressing the button again (or button, if you are in user memory 2), all memory values can be retrieved one after the other, starting with the most recent measured value.

- If an irregular heartbeat was detected during measurement, that information is also stored and displayed when the measured value is retrieved from the device's memory along with the systolic and diastolic blood pressure readings, pulse rate, time, date and year.
- You can cancel the memory recall at any time by pressing the START/STOP button . Otherwise, the device will switch off automatically after a few seconds.
- Even if the power supply is interrupted, e.g. when changing the batteries, the stored values will still be available.

Deleting saved readings

All data saved for each user can be separately deleted for user memory and user memory . To do so, press the button of the corresponding user memory (or). The average value will then appear on the display.

Press and hold the user memory button for 5 seconds: 'CL 00' appears on the display. All data in the selected user memory are now deleted. If you release the button ahead of time, no data will be deleted.

7. Explanation of error messages

Error that has occurred	Possible causes	Remedy
Monitor cannot be switched on	Batteries are missing, inserted incorrectly or flat.	Check batteries and insert four identical new batteries if necessary.
Cuff will not inflate	Cuff connector plug is incorrectly positioned in the connector socket of the monitor. Wrong cuff type connected.	Check the connection between the cuff connector plug and socket. Verify that only an approved Veroyal® cuff and connector are being used.
E1	Pulse could not be detected correctly.	Check whether the cuff has been applied correctly. Do not talk or move during the measuring procedure.
E2	The pulse could not be detected.	Check whether the cuff has been positioned correctly. Do not talk or move during the measuring procedure.
E3	The cuff was applied too firmly or too loosely.	Apply the cuff so that two fingers fit between the cuff and the upper arm. Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.
E4	There was an error during the measurement.	Please rest for at least 1 minute and take the measurement again.
E5	Inflation pressure is higher than 300 mmHg.	Please rest for at least 1 minute and take the measurement again.
E6	There is a system error.	Contact customer service if this error message appears.
	The batteries are almost empty.	Replace batteries.
Implausible measured values	Implausible measured values often occur due to inappropriate handling of the device or mistakes during the measuring procedure.	Please observe Chapter 5 'Measuring your blood pressure' as well as the safety information. Then repeat the measurement.

Switch the device off if an error symbol appears. Check for all possible causes and note the instructions for self-measurement from Chapter 2 'Important information'. Relax for a minute and then take the measurement again.

8. Device care

- Only ever use a soft, damp cloth to clean the device. Please do not use thinner, alcohol, detergents or solvents.
- The cuff can be cleaned carefully using a slightly damp cloth and mild soap solution. Do not completely immerse the cuff in water.
- It is recommended to clean and disinfect the cuff regularly or after each use, especially when used by several users, to prevent infection. The cuff should be disinfected, particularly on the inside, by wiping with a disinfectant. Use a disinfectant that is compatible with the cuff materials, e.g. 75 % ethanol or isopropyl alcohol. To protect from external influences, keep the device, the cuff and these instructions in the storage bag.
- If you store a device and cuff, do not place any heavy objects on the device and the cuff. Remove the batteries.

9. Warranty conditions

- We offer a 3-year warranty for this premium quality blood pressure monitor from the day of purchase and in accordance with the conditions listed below.
- Claims must be made during the warranty period. The date of purchase must be documented by an appropriately completed and stamped warranty certificate or proof of purchase.
- Within the warranty period, HARTMANN will replace or repair free of charge any faulty device components caused by material faults or manufacturing errors. This does not extend the warranty period.
- The device is intended for the purpose described in this instruction manual only.
- The warranty does not cover damage resulting from improper use or unauthorised interference. Accessory parts that are subject to wear and tear (batteries, cuffs, etc.) are excluded from the warranty. Claims for compensation are limited to the value of the merchandise; compensation for consequential damage or injury is expressly excluded.
- In warranty cases, please send the device with cuff together with the fully completed and stamped warranty certificate or purchase receipt directly or via your dealer to the Customer Services department for your country.

10. Contact details for customer queries

- HK PAUL HARTMANN Asia-Pacific Ltd.
Suite 3102-3103, 31/F Manhattan Place
23 Wang Tai Road, Kowloon Bay,
Kowloon, Hong Kong
www.veroyal.info
Tel: +852 2953 7100
Fax: +852 2796 7610
service@hk.hartmann.info
- ZA HARTMANN South Africa
Epsom Avenue, Northriding,
2169 Johannesburg
www.veroyal.info
Tel: +27 860 4278 6266
phzaelpdesk@hartmann.info

11. Technical data

Model	Veroyal® BPU 22
Type	GCE606
Measuring method	oscillometric
Display range	0 – 300 mmHg
Measuring range	Systolic (SYS): 50–280 mmHg Diastolic (DIA): 30–200 mmHg Pulse: 40–199 pulse beats per minute Displaying of correct values outside the measuring range cannot be guaranteed.
Display unit	1 mmHg
Technical measuring accuracy	Cuff pressure: ±3 mmHg Pulse: ±5 % of the displayed pulse rate
Clinical measuring accuracy	Complies with the requirements of DIN EN 1060-4; Kortokoff validation method: Phase I (SYS), Phase V (DIA)

Operating mode	Continuous mode
Nominal voltage	DC 6V
Power supply	4x 1.5 V alkaline manganese (AA/LR06) batteries
Expected service life	20,000 measurements
Battery capacity	Approx. 1,000 measurements
Protection against electric shock	Medical-electrical equipment with internal power supply (only when batteries are used); applied part: type BF
Protection against harmful ingress of water or solid materials	IP21 (not protected against moisture)
Inflation pressure	Comfort Air Technology: Individual inflation pressure depends on systolic blood pressure
Automatic switch-off function	1 minute after measurement/otherwise 30 sec.
Cuff	Veroyal® cuff for BPU 22, cuff for arm circumference of 22 to 42 cm Model number: 30601
Memory capacity	2 x 100 measurements with average of all measurements and morning/evening average of the last 7 days
Operating conditions	Ambient temperature: +10 °C to +40 °C Relative humidity: <90 %, non-condensing Air pressure: 800 hPa – 1050 hPa
Storage/transport conditions	Ambient temperature: –20 °C to +55 °C Relative humidity: <90 %, non-condensing
Serial number	inside the battery compartment
Reference to standards	IEC 60601-1; IEC 60601-1-2 (in accordance with CISPR11, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8)
Weight	Approx. 220 g (without batteries)
Dimensions	Approx. 134(L) x 48(W) x 91(H) mm

Legal requirements and guidelines

- The Veroyal BPU 22 conforms to the European regulations, which are based on the European Medical Device Directive 93