DE

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AT — PAUL HARTMANN Ges.m.b.H. · 2355 Wiener Neudorf

Instructions for use

Product description
Regocell Silicone Border is a sterile self-adhesive absorbent dressing with a silicone interface. The dressing comprises a perforated silicone wound contact layer, a semipermeable polyurethane backing film and a green absorbent pad in between those two layers. Acrylic adhesive is used to bond the dressing to the skin. The backing film is permeable to air but waterproof, which allows the patient to shower. The silicone layer facilitates the application of the dressing and prevents non-adherence to the wound. Above the silicone layer, the product features a white, hydrophilic non-woven fabric (viscose and polyester). The absorbent pad comprises cellulose and superabsorbent polyacrylate (SAP) for the absorption and retention of exudate and is wrapped in a hydrophilic cellulose tissue to distribute the liquid. It is covered with a hydrophobic green, polypropylene non-woven fabric to indicate the side facing away from the wound.

Properties and mode of action
Regocell Silicone Border absorbs and retains up to high levels of exudate in the absorbent pad and the dressing is able to redistribute pressure. It features a micro-adherent silicone interface and adhesive border, so no additional material is required to secure the dressing. The backing film is bacteria and showerproof. The dressing changes can be carried out atraumatically and almost painlessly.

Intended purpose
Regocell Silicone Border is a sterile single-use soft adhesive absorbent dressing with a silicone interface for the management of exuding wounds. It can be used as part of a prophylactic therapy for preventing skin damage such as pressure injuries. It can be used in clinical or homecare environments and can be combined with cleaning solutions and local antiseptics, primary dressings and additional fixation. It can be used in context of compression therapy. Regocell Silicone Border can be used by healthcare professionals and lay persons under the supervision of healthcare professionals.

Indications
Regocell Silicone Border is specifically for the management of a wide variety of exuding wounds: acute wounds (burns, surgical wounds, traumatic wounds) and chronic wounds (venous ulcers, diabetic foot ulcers, tumour wounds). It can also be used on intact skin as part of pressure injury prevention protocol.

Duration of use
Regocell Silicone Border can remain on the wound for up to 7 days, depending on the condition of the wound and the clinical decision/protocol. It is safe to continue using Regocell Silicone Border for as long as clinically indicated.

Mode of application
- Clean the wound according to your protocol. Make sure that the surrounding skin is dry and clean to ensure that the dressing sticks.
- Select the dressing size and shape, so that the green absorbent pad exceeds the wound edges by approx. 1 to 2 cm.
- If used for prevention, assure that the dressing pad covers the area at risk of pressure injury/ tissue damage. Where applicable, consider a dressing shape optimized for the specific anatomical site.
- Do not cut the green absorbent pad. The transparent border can be cut to adapt the dressing to specific body areas.
- Remove the protective sheets according to the application instructions. Apply the adherent side of the dressing, so that the green pad is on the wound.
- Press down the adhesive edges firmly and avoid wrinkles in the border to ensure the dressing adheres securely. Be careful when pressing down.
- No additional fixation is required. If indicated, compression bandages can be applied as instructed by a clinician.
- If used for prevention, the area at risk of pressure injury should be regularly inspected according to clinical practice. For the purpose of wound assessment, the dressing may be gently lifted and re-applied after.
- The dressing must be changed if clinically indicated or when the absorbent pad/wound margin of the green absorbent pad.

Contraindications
Do not use Regocell Silicone Border in case of hypersensitivity to any of its components.

Special precautions
- Do not cut the green absorbent pad. The transparent border can be cut to adapt the dressing to specific body areas.
- Before beating a wound that is not healing as expected, consult a clinician. Treatment with Regocell Silicone Border cannot replace a causal treatment of the wound-healing impairment. If the wound is infected, treatment of the infection needs to be controlled with appropriate treatment before this dressing can be used.
- In all cases, follow the instructions of the clinician.
- In the absence of available data supporting the use of this medical device on sensitive population groups, such as pregnant or nursing women, it should be used with caution and following a clinician's prescription.
- It is not recommended to stand on the dressing.
- The use of Regocell Silicone Border in the case of prophylactic pressure injury prevention therapy does not preclude the need to continue to develop and follow a comprehensive pressure injury prevention protocol, e.g. nutrition, hydration, repositioning, skin care and mobility.
- Before opening the product packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or moisture, do not use and contact your supplier.
- Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Side-effects
If any of the following side-effects occur, consult a clinician: increased exudate, maceration, increased pain, irritation, adherence of the dressing to the wound, itching, burning, oedemata, multiple ulceration, skin tears, blisters, erythema.

Incident reporting
If a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices), if, during the use of this device or as a result of its use, a serious incident occurs, report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal
To minimize the risk of potential infections hazards, or environmental pollution, the components of the medical device should follow disposal procedures according to applicable medical laws, rules, regulations and infection prevention standards. Dispose of the medical device with household waste.

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