

MoliCare® Premium Soft, incontinence slip

General Product Description/Intended Purpose

Intended Purpose:

MoliCare® Premium Soft products are non-sterile absorbent slips for incontinence care for incontinent persons, used in home care environment as well as healthcare facilities for the usage of bladder weakness, urinary and faecal incontinence with several severities, applied by lay person as well as healthcare professionals. The products are for single use only but a new device can be applied indefinitely.

Application/Indication

For several severities of urinary and faecal incontinence.



Reference Numbers*

Packaging		XS	S	M	L	XL
	Hip-/waist girth (cm)	40–60	60–90	90–120	120–150	150–175
P30	extra	169 248	169 448	169 648	169 848	
	super		169 450	169 650	169 850	
P26	super plus			169 671 / 169 308		
P24	super plus				169 871 / 169 408	
P14	extra					169 948
	super					169 950
P10	extra		169 498	169 098	169 198	
	super		169 698	169 298	169 398	
P2	extra			169 273	169 373	
	super			169 275	169 375	

*the full assortment may not be available in all countries

Single use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare® Premium Soft should follow disposal procedures according to applicable and

MoliCare® Premium Soft, incontinence slip

local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product performance characteristics

Characteristics	unit	Size	extra	super	super plus
total weight:*	[g]	XS	63		
		S	74	86	
		M	93	110	128
		L	110	126	143
		XL	127	161	

* $\leq \pm 10\%$ tolerance

Material characteristics

Part	Material
back sheet:	polyethylene, nonwoven-film-laminate, white
absorbent core (3-layer):	1. layer: kraft pulp* 2. layer: kraft pulp* - with superabsorbent polymer 3. layer: kraft pulp, „Curled Fiber“ *bleached elementary chlorine free
adhesive:	hotmelt adhesive
top sheet:	polypropylene-nonwoven, hydrophilic, white
cuff:	polypropylene-nonwoven, hydrophobic, white
elastics:	elastics for leg gather and standing cuffs
wetness indicator:	water soluble ink
side panel	polypropylene-nonwoven, hydrophobic, white
fixation tapes	3-layered with polypropylene-nonwoven, adhesive, polypropylene-strips (hooks)

MoliCare® Premium Soft, incontinence slip

Product characteristics

Dimensions	unit	Size	extra	super	super plus
total length:*	[mm]	XS	610		
		S	700	700	
		M	815	815	815
		L	980	980	980
		XL	980	980	
total width abdominal side:*	[mm]	XS	410		
		S	490	490	
		M	650	650	650
		L	780	780	770
		XL	930	930	
total width back side*	[mm]	XS	410		
		S	490	490	
		M	650	650	650
		L	780	780	770
		XL	930	930	
total width crotch:*	[mm]	XS	220		
		S	240	240	
		M	320	320	320
		L	320	320	320
		XL	320	320	
absorbent core length:*	[mm]	XS	520		
		S	540	540	
		M	660	660	650
		L	780	780	780
		XL	780	780	
absorbent core width abdominal side:*	[mm]	XS	190		
		S	190	190	
		M	200	200	200
		L	200	200	200
		XL	240	240	
absorbent core width back side:*	[mm]	XS	190		
		S	190	190	
		M	245	245	245
		L	265	265	255
		XL	350	350	
absorbent core width crotch:*	[mm]	XS	135		
		S	135	135	
		M	170	170	170
		L	170	170	170
		XL	180	180	

*≤ ± 10% tolerance

MoliCare® Premium Soft, incontinence slip

Labelling

9-digit LOT-code
e.g.: 900102822



9
year 201X


001
internal key

02
internal key

82
internal key

2
cumulated
check sum

Date of manufacture


e.g.: 

2019
year

08
month

07
day

Use-by-Date

e.g.: 

2024
year

08
month

01
day

Shelf life: 5 years

Medical Device



Keep dry



Do not reuse



CE-marking



Latest date of revision: 2020-05-28