

MoliCare® Comfort, incontinence slip

General Product Description/Intended Purpose

Intended Purpose:

MoliCare® Comfort products are non-sterile absorbent slips for incontinence care for incontinent persons, used in home care environment as well as healthcare facilities for the usage of bladder weakness, urinary and faecal incontinence with several severities, applied by lay person as well as healthcare professionals. The products are for single use only but a new device can be applied indefinitely.

Application/Indication

For several severities of urinary and faecal incontinence.



Reference Numbers*

Packaging		S	M
	Hip-/waist girth (cm)	60-90	90-120
P30	super	169 183	169 283

*the full assortment may not be available in all countries

Single use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare® Comfort should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

MoliCare® Comfort, incontinence slip

Product performance characteristics

Characteristics	unit	Size	super
		S	78
		M	108
		L	124

* $\leq \pm 10\%$ tolerance

Material characteristics

Part	Material
back sheet:	polyethylene
absorbent core (3-layer):	1. layer: kraft pulp* 2. layer: kraft pulp* - with superabsorbent polymer 3. layer: kraft pulp, „Curled Fiber“ *bleached elementary chlorine free
adhesive:	hotmelt adhesive
top sheet:	polypropylene-nonwoven, hydrophilic, white
cuff:	polypropylene-nonwoven, hydrophobic, white
elastics:	elastics for leg gather and standing cuffs
wetness indicator:	water soluble ink
side panel	polypropylene-nonwoven, hydrophobic, white
fixation tapes	polypropylene, adhesive

Product characteristics

Dimensions	unit	Size	super
total length:*	[mm]	S	700
		M	815
total width abdominal side:*	[mm]	S	490
		M	650
total width back side*	[mm]	S	490
		M	650
total width crotch:*	[mm]	S	215
		M	270
absorbent core length:*	[mm]	S	540
		M	650
absorbent core width abdominal side:*	[mm]	S	190
		M	300
absorbent core width back side:*	[mm]	S	190
		M	320
absorbent core width crotch:*	[mm]	S	135
		M	175

* $\leq \pm 10\%$ tolerance

MoliCare® Comfort, incontinence slip

Labelling

9-digit LOT-code
e.g.: 900102822



9
year 201X


001
internal key

02
internal key

82
internal key

2
cumulated
check sum

Date of manufacture


e.g.: 

2019
year

08
month

07
day

Use-by-Date

e.g.: 

2024
year

08
month

01
day

Shelf life: 5 years

Medical Device



Keep dry



Do not reuse



CE-marking



Latest date of revision: 2020-05-28