

MoliCare® Fixpants long leg, fixation pants

General Product Description/Intended Purpose

Intended Purpose:

MoliCare® Fixpants long leg are non-invasive, non-sterile fixation systems for safe and hygienic securing of all types of incontinence pads for long term contact with the body used by lay persons and/or health care professionals for humans on intact skin.

MoliCare® Fixpants long leg are light, soft and permeable to air, made of an extra-strong material with especially tight structure, full run resist.

These unisex fixation pants are reusable and can be used with MoliCare® Premium Form, MoliCare® Form and other types of incontinence pads as well as pads for women in childbed and are suitable for both immobile/ bedridden and mobile/active patients/residents (women and men).


Thanks to their high stretchability, the change and the correct positioning of the pad is facilitated. The size identification is given by the woven-in colour coding system.

The products are suitable for usage of bladder weakness and urinary/faecal incontinence with several severities.

MoliCare® Fixpants long leg can be used in a home care environment and professional healthcare facilities.

Application/Indication

Several severities of urinary and faecal incontinence



Reference Numbers*

Packaging	P5	P50
Size M	947763	947741
Size L	947764	947742
Size XL	947765	947743
Size XXL	947766	947744

*the full assortment may not be available in all countries

MoliCare[®] Fixpants long leg, fixation pants

Reusable Device

The product can be washed up to 25 times according to the washing instructions on the bag. Not following the washing instructions might have an impact on the elastic fibre elasticity and product fit/ fixation performance.
Information available on request.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare[®] Fixpants long leg should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product performance characteristics

Characteristics	unit	M	L	XL	XXL
total weight:*	[g]	11,8	13,5	14,0	14,3

* $\pm 1,0$ g tolerance

Material characteristics

Material	Quantity [%]
Polyamide	96
Elastane	4

Product characteristics

Dimensions	M	L	XL	XXL
waist 200g load: [mm]:*	260	270	280	290
seam height 1200 g load: ** [mm]:	450	510	510	530
waist 1200g load*** [mm]:	510	550	590	630
leg 1000g load:*** [mm]:	365	385	415	435

* ± 40 mm tolerance

** ± 70 mm tolerance

*** ± 50 mm tolerance

MoliCare® Fixpants long leg, fixation pants

Labelling

8-digit LOT-code

e.g.: 20089822



XX

XX

XX

XX


month of production

day of production

shift number

machine number

Manufacturing date


e.g.: 

2020
year

01
month

08
day

Expiry date

e.g.: 

2025
year

01
month

08
day

Shelf life: 5 years

Medical Device



Keep dry



CE-marking



Wash at or below 60 °C



Do not bleach



Tumble drying (low temperature)



Do not iron



Do not dry clean



Latest date of revision: 2020-02-03