

Technical Data Sheet

Foliodrape® Protect, Nephroscopy Set (Art. No. 938 752)



SpecNo.:	D 6.1004d
Deparment:	CMO-DOE
Date:	2016-01-25

1. General Product Description

Foliodrape® Protect Nephroscopy Set for surgical draping in nephroscopic surgery.

The Foliodrape® Protect Nephroscopy Set carries the CE-mark according to EU directive 93/42/EWG for medical devices and is classified as class I sterile medical devices. A conformity assessment has been performed for Foliodrape® Protect and it has been shown to be in compliance with all applicable requirements of the above mentioned directive. The safe use of Foliodrape® Protect Nephroscopy Set, therefore, is ensured if the product is used in line with the intended purpose.

2. Application / Indication

Foliodrape® Protect Nephroscopy Set is a disposable product for sterile draping of the patient and the theatre equipment.

3. <u>Presentations</u>

Art. No. 938 752	Foliodrape® Protect Nephroscopy Set	
	sterile, 10 Sets per dispenser, 1 dispenser per transport carton	
	Dimensions of dispenser: 545 x 374 x 234 mm (L x W x H)	
	Dimensions of transport carton: 590 x 390 x 250 mm (L x W x H)	

4. Product Characteristics

Content:	1 reinforced table cover 140 x 190 cm	
	1 nephroscopy drape sheet 300 x 170 cm with fenestration (15 x 20 cm, self-adhesive) and fluid collection pouch (with integrated sieve)	
300x170cm	2 cellulose towels 30 x 33 cm	
Sterilization:	With ethylene oxide gas according to EN 556-1 and EN 11135-1	
Self adhesives:	Skin-friendly synthetic rubber, latex free and colophony free	



Technical Data Sheet

Foliodrape® Protect, Nephroscopy Set (Art. No. 938 752)



SpecNo.:	D 6.1004d
Deparment:	CMO-DOE
Date:	2016-01-25

5. Product Requirements

Foliodrape® Protect corresponds to all requirements of European standard for surgical drapes EN 13795 for the highest performance level "High performance, critical area".

For Details see Technical Data sheet of the base material (D 6.1005a) and back table cover (D 6.1005d).

6. Labelling

Lot-No. with 9-digit code:

e.g.: 6023XXXXX

LOT	6	023	XXXXX
	year of production	serial production order	for internal use

Expiry date:

e.g. 2021 – 01

year month

Shelf life: 5 years

The product will be sterile during this period as long as the packaging is not damaged and the storage conditions are kept.

7. Packaging

Foliodrape® Protect Nephroscopy Set is sealed into sterile packaging (peeling pouch) according to EN ISO 11607 and EN 868. The peel pouch is packed into folding box and transport carton, sealed with adhesive tapes and stored on euro-pallet.

Date: 2016-01-25

PAUL HARTMANN AG
Development OTM & Exam. Gloves (CMO-DOE)

i. A.

Evelyn Hartmann

Dr. Evelyn Hartmann