

		<h1>Technical Data Sheet</h1>
--	--	-------------------------------

Foliodrape® Protect, Arthroscopy Drape (Art. No. 938 820)		Spec.-No.:	D 6.1008a
		Department:	CMO-DOE
		Date:	2015-01-21

1. General Product Description

Foliodrape® Protect Arthroscopy Drape, for knee-joint arthroscopic procedures.

Foliodrape® Protect Arthroscopy Drape carries the CE-mark according to EU directive 93/42/EWG for medical devices and is classified as class I sterile medical devices.

A conformity assessment has been performed for *Foliodrape® Protect* and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

The safe use and effectiveness of *Foliodrape® Protect Arthroscopy Drape*, therefore, is ensured if the product is used in line with the intended purpose.


2. Application / Indication

Foliodrape® Protect Arthroscopy Drape is a disposable product for sterile draping of the patient and the theatre equipment.

3. Presentations

Art. No. 938 820	<i>Foliodrape® Protect Arthroscopy Drape</i>
	sterile, 6 drapes per dispenser, 1 dispenser per transport carton
	Dimensions of dispenser: 545 x 374 x 234 mm (L x W x H)
	Dimensions of transport carton: 590 x 390 x 250 mm (L x W x H)

4. Product Characteristics

Content: 		1 arthroscopy drape 320 x 245 cm with elastic fenestration, fluid collection pouch and velcro fastener 1 drainage tube (silicone) and clamp
Sterilization:	With ethylene oxide gas according to EN 556-1 and EN 11135-1	
Self adhesives:	Skin-friendly synthetic rubber, latex free and colophony free	

		Technical Data Sheet
--	--	----------------------

Foliodrape® Protect, Arthroscopy Drape (Art. No. 938 820)		Spec.-No.:	D 6.1008a
		Department:	CMO-DOE
		Date:	2015-01-21

5. Product Requirements

Foliodrape® Protect corresponds to all requirements of European standard for surgical drapes EN 13795, for the highest performance level "High performance, critical area".

For Details see Technical Data sheet of the base material (D 6.1005a).


6. Labelling

Lot-No. with 9-digit code:

e.g.: 4023XXXXX

LOT	4	023	XXXXX
	year of production	serial production order	for internal use

Expiry date:

e.g.  2019 – 01
year month

Shelf life: 5 years

The product will be sterile during this period as long as the packaging is not damaged and the storage conditions are kept.

7. Packaging

Foliodrape® Protect Arthroscopy Drape is wrapped in crepe paper and sealed into sterile packaging (peeling pouch) according to EN ISO 11607 and EN 868. The peel pouch is packed into folding box and transport carton, sealed with adhesive tapes and stored on euro-pallet.

Date: 2015-01-21

PAUL HARTMANN AG
Development OTM & Exam. Gloves (CMO-DOE)

i. A.



Dr. Evelyn Hartmann