
		Technical Data Sheet		
Foliodrape® COMFORT Cardiovascular Set I (Art. No. 257 128)			Spec.-no.:	D 6.5339a
			Department:	CMO-DOE
			Date:	2013-01-31

1. General Product Description

Foliodrape® COMFORT Cardiovascular Set I for the sterile draping in cross-leg bypass surgery.

The Foliodrape® COMFORT Cardiovascular Set I carries the CE-mark according to EU directive 93/42/EEG for medical devices and is classified as class I medical devices.

A conformity assessment has been performed for Foliodrape® COMFORT and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

The safe use of Foliodrape® COMFORT Cardiovascular Set I, therefore is ensured if the product is used in line with the intended purpose.

2. Application / Indication



Foliodrape® COMFORT Cardiovascular Set I is a disposable product for sterile draping of the patient and the theatre equipment.

3. Presentations:

Art. No. 257 128	Foliodrape® Comfort Cardiovascular Set I
	sterile, 3 Sets per dispenser, 1 dispenser per transport carton
	Dimensions of dispenser: 565 x 370 x 295 mm (L x W x H)
	Dimensions of transport carton: 585 x 385 x 315 mm (L x W x H)

4. Product characteristics

Content:	1 table cover 140 x 190 cm
	1 thorax drape sheet 260 x 315 cm
	1 slit drape sheet 200 x 260 cm with perineal cover, self-adhesive
	1 drape sheet 150 x 175 cm (2-layer)
	2 foot covers 35 x 25 cm
	1 mayo stand cover 80 x 145 cm
	2 adhesive tape 10 x 50 cm
	4 cellulose towels 33 x 33 cm
Sterilization:	With ethylene oxide gas

		<h1>Technical Data Sheet</h1>		
Foliodrape® COMFORT Cardiovascular Set I (Art. No. 257 128)			Spec.-no.: Department: Date:	D 6.5339a CMO-DOE 2013-01-31

5. Product requirements:

Foliodrape® COMFORT corresponds to all requirements of European standard for surgical drapes EN 13795, for the highest performance level "High performance, critical area".

For details see Technical Data sheet of the base material (D 6.5030), table cover (D 6.5410a), and the mayo stand cover (D 6.5403).

6. Labelling Lot no. with 9-digit code:

e.g.: 9XXX38YYY

LOT	9	XXX	38	YYY
	year of production	for internal purpose only	week of production	for internal purposes only

Expiry date:

e.g. ⌚ 2014 – 01
year month

Shelf life: 5 years

The product will be sterile during this period as long as the packaging is not damaged and the storage conditions are kept.

7. Packaging

Foliodrape® COMFORT Cardiovascular Set I is sealed into sterile packaging (peeling pouch) according to EN ISO 11607 and EN 868. The peel pouch is packed into folding box and transport carton, sealed with adhesive tapes and stored on euro-pallet.

Date: 2013-01-31

PAUL HARTMANN AG
CMO-DOE Development OTM & Exam Gloves

i. A.



Dr. Thomas Hoffmann