

Peha-crepp® E



Spec.-No.: D 2.5410
 Department: Dev-KOB
 Date: 2017-04-27

1. General Product Description

Peha-crepp® E is a highly-elastic conforming bandage acc. DIN 61 634-FB, consisting of a softly special fabric with a high proportion of natural fibres, elasticity approx. 160% after 24h standard atmosphere; quick and easy to apply without reverse turns to secure the bandage; reliable non-slip dressing retention; may be sterilized; for single use only; white.

Peha-crepp® E carries the CE mark according to EU directive 93/42/EEC for medical devices. The product is classified as a class I medical device.

A conformity assessment has been performed for **Peha-crepp® E** and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

The safe use of **Peha-crepp® E**, therefore, is ensured if the product is used in line with the intended purpose.

2. Application / Indication

Dressing retention, especially on joints and conical or rounded parts of the body.

3. Presentations

Dimensions	Unit	Units per Case
individually wrapped in PP-foil		
4cm x 4m (stretched length)	OP 1 bandage	20/160 OP
6cm x 4m (stretched length)	OP 1 bandage	20/200 OP
8cm x 4m (stretched length)	OP 1 bandage	20/200 OP
10cm x 4m (stretched length)	OP 1 bandage	20/120 OP
12cm x 4m (stretched length)	OP 1 bandage	20/80 OP
in cardboard box		
4cm x 4m (stretched length)	AP 20 bandages	8 AP
6cm x 4m (stretched length)	AP 20 bandages	10 AP
8cm x 4m (stretched length)	AP 20 bandages	10 AP
10cm x 4m (stretched length)	AP 20 bandages	6 AP
12cm x 4m (stretched length)	AP 20 bandages	4 AP
in compact packaging		
6cm x 4m (stretched length)	AP 100 bandages	8 AP
8cm x 4m (stretched length)	AP 100 bandages	6 AP
10cm x 4m (stretched length)	AP 100 bandages	6 AP

Peha-crepp® E



Spec.-No.: D 2.5410
Department: Dev-KOB
Date: 2017-04-27

4. Product Characteristics

Material Composition: 70% Viscose
30% Polyamide

Construction: knitted fabric
elastic bandage with fast edges

5. Product Requirements

Extensibility: approx. 160% after 24 hours
approx. 125% directly

Weight per unit area: approx. 27g/m² stretched

6. Labelling

Lot-No. with 8-digit code:

e.g.:

5
year

50
week of production

XXXXX
for internal purposes only

7. Packaging

Peha-crepp® E individually wrapped in PP-foil and packed in folding box;

or:

Peha-crepp® E bandages packed in cardboard box;

Transport carton acc. DIN, sealed with adhesive tapes, packed onto europallet.

Date: 2017-04-27



- Department Development (Dev-KOB) -