

Peha-Haft color®

General Product Description/Intended Purpose

Peha-Haft color® is a cohesive conforming bandage for retention of all kinds of primary and secondary dressings, splints, casts, catheters and cannulas (etc.) particularly on joints as well as on conical and round parts of the body. It is used for the stabilization of muscles, tendons, ligaments, joints with low immobilization effect as well. It is manufactured from 40% cotton, 31% viscose and 29% polyamide. It is for single-use and can be used by instructed lay persons and health professionals on intact skin in clinical and home environment.

Application/Indication

Peha-Haft color® is a cohesive conforming bandage that may be used in medical interventions when acute and chronic wounds need to be treated short or long term or when slight immobilization is needed in combination with medical device. It is suitable for fixation of:

- primary and secondary dressings
- cannulae, catheters, tubes and similar medical devices
- splints, casts and similar medical devices

Reference Numbers

REF	Name	Size	Bandages per folding box	Individually wrapped
9324602	Peha-haft color	6cmx20m	1	No
9324612	Peha-haft color	8cmx20m	1	No
9324622	Peha-haft color	10cmx20m	1	No
9324732	Peha-haft color	6cmx20m	1	No
9324742	Peha-haft color	8cmx20m	1	No
9324752	Peha-haft color	10cmx20m	1	No

Residual Risks, Contra-Indications and Undesirable Side Effects, Warnings

Not known.

Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of Peha-Haft color® should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Peha-Haft color®

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.


Product Characteristics

Property	Requirements
Width	2,5 cm; 4,0 cm; 6,0 cm; 8,0 cm; 10,0 cm; 12,0 cm
Length, stretched	4m; 20m
Elongation, MD	80 %
Cohesive force	45 ± 25 cN/cm
Unwinding force	15-125 cN/cm
Weight per area (stretched)	39,5 g/m ²


Labelling

Lot-No. with 9-Digit Code


e.g.:

	9	999	01	00	1
explanation	[9 = 2019]	[999-000]			
definition	year	serial production order			for internal purposes onl

Manufacturing Date

e.g.:		2015	04	07
		year	month	day

Expiry Date

e.g.:		2018	04	07
		year	month	day

Shelf Life: 3 years

Medical Device



Latest Date of Revision: 2020-01-08