

## Hydrofilm I.V.



Spec.-No.: D7.4019b

Department: SM-PDW

Date: 29.11.2010

### 1. General Product Description

Hydrofilm I.V. is a sterile self-adhesive dressing for the fixation of peripheral intra venous cannulae.

Hydrofilm I.V. is made of a transparent semi-permeable, waterproof polyurethane film coated with a hypoallergenic acrylate adhesive.

Hydrofilm I.V. carries the CE-mark according to EU directive 93/42/EEC for medical devices.

The product is classified as a class I, rule 4 (1) medical device.

A conformity assessment has been performed for Hydrofilm I.V. control and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

The safe use and effectiveness of Hydrofilm I.V., therefore, is ensured if the product is used in the line with the intended purpose.

### 2. Application / Indication

To fix peripheral intra venous cannulae and for the sterile treatment of the exit site.

### 3. Presentations

Sterile, individually sealed.

Ref. No.	Dimensions	Number peel pouches per folding box
685740/0	7 cm x 9 cm;	P50

### 4. Product Characteristics

Composition:

Backing: transparent polyurethane film

Adhesives: hydroallergenic acrylate adhesive

Liner: silicone paper

Handling aid: two-pieces supporting film (based on polyester)

## Hydrofilm I.V.



Spec.-No.: D7.4019b  
Department: SM-PDW  
Date: 29.11.2010

### 5. Product Requirements

Hydrofilm I.V. creates a strong and reliable fixation of the cannula.

Easy handling due to a four-step application system.

Moisture vapor transmission rate (MVTR) according to DIN EN 13726-2 > 1100 g/m<sup>2</sup>/24h

ETO sterilized

### 6. Labelling

Lot-No. with 9-digit code:  
e.g.:

X	XXX	XX	XXX
year	order number	week	for internal purposes only


Expiry date:  
e.g.: 2015 1  
year month

Shelf life: 5 years

### 7. Packaging

Hydrofilm I.V. control is individually sealed in peel pouches; packed into folding boxes; transport carton acc. DIN; sealed with adhesive tapes, packed onto europallet.

Date: 2010-11-29



-Sybille Hörmann/SM-PDW-