

MoliCare® Premium MEN PANTS

General Product Description/Intended Purpose

Intended Purpose:

MoliCare® Premium MEN PANTS are pull-up pants for incontinence care for men.

The products are suitable for usage in case of bladder weakness with several severities, urinary and fecal incontinence with several severities and for usage of toilet training.

The products are intended for mobile to semi-mobile patients and restless or dement patients. The products are for single use only and should be regular changed as needed.

The products can be used in home care environment and professional healthcare facilities. The application can be performed by lay person and/or health care professional.

Application/Indication

MoliCare® Premium MEN PANTS are suitable for usage in case of bladder weakness with several severities, urinary and fecal incontinence with several severities and for usage of toilet training.



Reference Numbers*

Packaging		M	L
	Hip-/waist girth (cm)	80–120	100–150
P8	5 drops / green	915817	-
	7 drops / blue	915827	-
P7	5 drops / green	-	915818
	7 drops / blue	-	915828
P1	5 drops / green	915802	915803
	7 drops / blue	915807	915808

*the full assortment may not be available in all countries

MoliCare® Premium MEN PANTS

Single use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare® Premium MEN PANTS should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product performance characteristics

MoliCare® Premium MEN PANTS 5 drops / green:

Characteristics	unit	M	L
total weight:*	[g]	55	62

* ≤ ± 10% tolerance

MoliCare® Premium MEN PANTS 7 drops / blue:

Characteristics	unit	M	L
total weight:*	[g]	65	71

* ≤ ± 10% tolerance

Material characteristics

Part	Material
chassis	polypropylene non-woven, hydrophobic, breathable, grey
back sheet:	polyethylene film with polypropylene nonwoven layer, breathable, grey
absorbent core (3-layer):	1. layer: kraft pulp* 2. layer: kraft pulp* - with superabsorbent polymer 3. layer: kraft pulp*, „Curled Fiber“ *bleached elementary chlorine free
adhesive:	hotmelt adhesive
top sheet:	polypropylene - nonwoven, hydrophilic, white
cuff:	polypropylene - nonwoven, hydrophobic, white
elastics:	polyurethane based elastics for leg gathers, hip area, waist and standing cuffs
product information:	printing on inside of waistband (rear side)

MoliCare® Premium MEN PANTS


Product characteristics

Dimensions	M	L
total length*:	750	850
total width*:	660	810
total width crotch*:	245	245
absorbent core length*:	470	520
absorbent core width*:	180	180
absorbent core width crotch*:	110	110


* ≤ ± 10% tolerance

Labelling


9-digit LOT-code

e.g.:  9 001 02 82 2
year 201X internal key internal key internal key cumulated
check sum

Manufacturing date

e.g.:  2019 08 07
year month day

Expiry date

e.g.:  2024 08 01
year month day

Shelf life: 5 years

Medical Device



Keep dry



Do not reuse



CE-marking



Latest date of revision: 2019-10-04