

Peha-Haft®

General Product Description/Intended Purpose

Peha-Haft® is a cohesive conforming bandage for retention of all kinds of primary and secondary dressings, splints, casts, catheters and cannulas (etc.) particularly on joints as well as on conical and round parts of the body. It is used for the stabilization of muscles, tendons, ligaments, joints with low immobilization effect as well. It is manufactured from 37% cotton, 43% viscose and 20% polyamide. It is for single-use and can be used by instructed lay persons and health professionals on intact skin in clinical and home environment.

Application/Indication

Peha-Haft® is a cohesive conforming bandage that may be used in medical interventions when acute and chronic wounds need to be treated short or long term or when slight immobilization is needed in combination with medical device. It is suitable for fixation of:

- primary and secondary dressings
- cannulae, catheters, tubes and similar medical devices
- splints, casts and similar medical devices

Reference Numbers

REF	Name	Size	Bandages per folding box	Individually wrapped
9324414	Peha-haft	4cmx4m	1	No
9324426	Peha-haft	6cmx4m	1	No
9324435	Peha-haft	8cmx4m	1	No
9324445	Peha-haft	10cmx4m	1	No
9324454	Peha-haft	12cmx4m	1	No
9324462	Peha-haft	4cmx20m	1	No
9324472	Peha-haft	6cmx20m	1	No
9324482	Peha-haft	8cmx20m	1	No
9324492	Peha-haft	10cmx20m	1	No
9324502	Peha-haft	12cmx20m	1	No
9324372	Peha-haft	6cmx20m	1	No
9324382	Peha-haft	8cmx20m	1	No
9324392	Peha-haft	10cmx20m	1	No
9324414	Peha-haft	4cmx4m	1	No
9324454	Peha-haft	12cmx4m	1	No

Peha-Haft®

Residual Risks, Contra-Indications and Undesirable Side Effects, Warnings

Not known.

Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of Peha-Haft® should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product Characteristics

Property	Requirements
Width	2,5 cm; 4,0 cm; 6,0 cm; 8,0 cm; 10,0 cm; 12,0 cm
Length, stretched	4 m; 20 m
Elongation, MD	80 %
Cohesive force	45 cN/cm
Unwinding force	15-125 cN/cm
Weight per area (stretched)	39,5 g/m ²

Labelling

Lot-No. with 9-Digit Code

e.g.:

	9	999	01	00	1
explanation	[9 = 2019]	[999-000]			
definition	year	serial production order			for internal purposes onl

Manufacturing Date

e.g.:		2015	04	07
		year	month	day

Peha-Haft®

Expiry Date

e.g.:  2018 04 07
year month day

Shelf Life: 3 years

Medical Device

MD

Technical Data Sheet

Peha-Haft®

Latest Date of Revision: 2020-08-01