

# MoliCare® MEN PAD


## General Product Description/Intended Purpose

### Intended Purpose

MoliCare® MEN PAD are single use, non-sterile unisex (except 2 drops) devices with fluff insert for short long-term management of light urine & fecal incontinence in adults used in genital area used by lay persons and health professionals in clinical and/or home environment for more than 30 days.

### Application/Indication

MoliCare® MEN PAD products are for absorption of urine and fecals due to incontinence.



### Reference Numbers\*

Products	Version	Reference Number	Packaging
MoliCare® MEN PAD	MoliCare® MEN PAD 2 drops (light grey)	168600	14
	MoliCare® MEN PAD 3 drop (green)	168603	14
	MoliCare® MEN PAD 4 drop (blue)	168705/168707	14
	MoliCare® MEN PAD 5 drops (orange)	168801/168807	14
	MoliCare® MEN PAD 2 drops (light grey)	168 503	2
	MoliCare® MEN PAD 3 drop (green)	168 604	2
	MoliCare® MEN PAD 4 drop (blue)	168 501	2
	MoliCare® MEN PAD 5 drops (orange)	168 609	2

\* may not be available in all countries

### Single use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

### Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare®MEN PAD should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

### Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

# MoliCare® MEN PAD

## Product performance characteristics

MoliCare®MEN PAD

Characteristics	2 drops	3 drops	4 drops	5 drops
Total weight* [g]	13.9	15.0	22.1	34.4
Absorbent core weight* [g]	11.3	12.0	18.0	30.3

\* $\leq \pm 10\%$  tolerance

## Material characteristics

Component	Material
Topsheet	Polypropylene – nonwoven, hydrophilic,
Absorbent core (3 Layer)	1. Layer kraft pulp* 2. Layer kraft pulp* + superabsorbent polymer 3. Layer:kraft pulp "curled fiber" (except 3 drops)
Backsheet	Polyethylene, nonwoven-film-laminate, white
Elastics	polyurethane foam
Cuffs	Polypropylene – nonwoven, hydrophobic, white (except 2 drops)
Adhesives	Hotmelts
Backing	Silicon Paper

\*bleached elementary chlorine free

## Product characteristics

MoliCare®MEN PAD

Dimensions in mm	2 drops	3 drops	4 drops	5 drops
Total length	260	260	345	345
Total width	145	180	180	180
Length of absorbent core	240	215	307	307
Width of absorbent core	110	140	152	152

\* $\leq \pm 10\%$  tolerance

# MoliCare® MEN PAD

## Labelling


### 9-digit LOT-code

e.g.: 900102822




9	001	02	82	2
year 201X	internal key	internal key	internal key	cumulated check sum

### Manufacturing date

e.g.: 

2019	10	10
year	month	day

### Expiry date

e.g.: 

2022	10	10
year	month	day

Shelf life 5 years

### Medical Device



### Unique Device Identification (UDI)



### Keep dry



### Do not reuse



### CE-marking



**Latest date of revision:** 2020-01-09