

MoliCare® Pad (P30)

General Product Description/Intended Purpose


Intended Purpose

MoliCare® Pads are single use, non-sterile unisex devices with fluff insert for short and long-term management of light urine & fecal incontinence in adults used in genital area used by lay persons and health professionals in clinical and/or home environment for more than 30 days.

Application/Indication

MoliCare® Pad products are for absorption of urine and fecals due to incontinence.

Application technique



Reference Numbers*

Products	Version	Reference Number	Packaging
MoliCare® Pad	MoliCare® Pad 2 drops mini P30 (green)	166806	30
	MoliCare® Pad 3 drops midi P30 (turquoise)	166807	30
	MoliCare® Pad 4 drops maxi P30 (purple)	166808	30

** may not be available in all countries

Single use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare®Pad should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

MoliCare® Pad (P30)

Product performance characteristics

Characteristics	MoliCare®Pad 2 drops mini	MoliCare®Pad 3 drops midi	MoliCare®Pad 4 drops maxi
Total weight* [g]	12,7	18,7	33,6
Absorbent core weight* [g]	10,5	15,6	29,2

*≤ ± 10% tolerance

Material characteristics

Component	Material
Topsheet	Polypropylene – nonwoven, hydrophilic,
Absorbent core (3 Layer)	1. Layer kraft pulp* 2. Layer kraft pulp* + superabsorbent polymer 3. Layer:kraft pulp "curled fiber"
Backsheet	Polyethylene, nonwoven-film-laminate, white
Elastics	polyurethane foam
Adhesives	Hotmelts
Backing	Silicon paper

*bleached elementary chlorine free

Product characteristics

Dimensions in mm	MoliCare®Pad 2 drops mini	MoliCare®Pad 3 drops midi	MoliCare®Pad 4 drops maxi
Total length	270	335	430
Total width	100	132	162
Length of absorbent core	240	305	395
Width of absorbent core	75	112	138

*≤ ± 10% tolerance

Labelling

9-digit LOT-code

e.g.: 900102822



9
year 201X


001
internal key

02
internal key

82
internal key

2
cumulated
check sum

Manufacturing date


e.g.: 

2019
year

10
month

10
day

Expiry date

e.g.: 

2022
year

10
month

10
day

Technical Data Sheet

MoliCare® Pad (P30)

Shelf life 5 years

Medical Device



Unique Device Identification (UDI)



Keep dry



Do not reuse



CE-marking



Latest date of revision: 2020-01-09