

MoliCare® Premium Bed Mat Textile

General Product Description/Intended Purpose

MoliCare® Premium Bed Mat Textile are reuse, non-invasive, non-sterile bed mats for bladder weakness and urinary incontinence care with several severities for humans used by lay persons or healthcare professionals.

Application/Indication

The products are suitable for usage of bladder weakness and urinary incontinence with several severities.

Reference Numbers*

Packaging	7 drops	7 drops with wings
P10	155 011	155 013
P1	155 010	155 012

* the full assortment may not be available in all countries

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare® Premium Bed Mat Textile should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Material Characteristics

Part	Material
backsheet	laminated of polyester and polyurethane film
absorbent core	mix of polyester and viscose
topsheet	polyester
seams	polyester
wings*	mix of polyester and cotton

* 7 drops with wings

Product Characteristics

Characteristics	unit	7 drops	7 drops with wings
total weight:*	[g]	435	491
total length*	[cm]	85	75
total width*	[cm]	90	185

* $\leq \pm 10\%$ tolerance

MoliCare® Premium Bed Mat Textile

Labelling

Lot-No. with 8 Digit Code

e.g.: 900102822



2020

02

07

Year

month

day

Manufacturing Date

e.g.:



2020

02

07

year

month

day

Expiry Date

e.g.:



2025

02

01

year

month

day

Medical Device



Keep dry



CE-marking



Distributor



Distributor

Latest Date of Revision: 2019-11-18