

MoliCare® Bed Mat Eco

General Product Description/Intended Purpose

MoliCare® Bed Mat Eco single use, non-invasive, non-sterile bed mats for bladder weakness and urinary incontinence care with several severities for humans used by lay persons or healthcare professionals.

Application/Indication

The products are suitable for usage of bladder weakness and urinary incontinence with several severities.

Reference Numbers

Packaging	40x60	60x90
P300	5 drops	159 001
P100	5 drops	159 002
	7 drops	159 003
	9 drops	159 005
P50	7 drops	159 004
	9 drops	159 006
P5	9 drops	159 007

* the full assortment may not be available in all countries

Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare® Bed Mat Eco should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Material Characteristics

Part	Material
backsheet	polyethylene, green polyethylene, blue*
absorbent core	recycled cellulose layers
adhesive	hotmelt adhesive
topsheet	polypropylene-nonwoven, hydrophilic, white

* MoliCare Bed Mat Eco 9 drops

MoliCare[®] Bed Mat Eco

Product Characteristics

Characteristics	unit	5 drops 40x60	5 drops 60x90	7 drops 40x60	7 drops 60x90	9 drops 40x60	9 drops 60x90
total weight:*	[g]	27,2	64,1	47,2	113,3	73,8	178,3
total length**	[cm]	60	90	60	90	60	90
total width**	[cm]	40	60	40	60	40	60

* ≤ ± 10 % tolerance

** ≤ - 3 % tolerance on the surface area

Labelling

Lot-No. with 9-Digit Code

e.g.: 900102822



9	001	02	82	2
Year 201X	Internal key	Internal key	Internal key	cumulated check sum

Manufacturing Date

e.g.:  2020 year 02 month 07 day

Expiry Date

e.g.:  2025 year 02 month 01 day

Medical Device



Keep dry



Do not reuse



CE-marking



Distributor



Latest Date of Revision: 2019-11-18