

MoliCare® Premium Bed Mat *(French Tender Assortment)*

General Product Description/Intended Purpose

MoliCare® Premium Bed Mat *(French Tender Assortment)* are single use, non-invasive, non-sterile bed mats for bladder weakness and urinary incontinence care with several severities for humans used by lay persons or healthcare professionals.

Application/Indication

The products are suitable for usage of bladder weakness and urinary incontinence with several severities.

Reference Numbers*

Packaging		40x60	60x60	60x90
P60	5 drops		161 080	
P30	7 drops	161 081 161 082		161 085
P60	7 drops		161 083 161 084	
P30	8 drops	161 086	161 087	161 088 161 089
P25	9 drops			161 090

*the full assortment may not be available in all countries

Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare® Premium Bed Mat should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Material Characteristics

Part	Material
backsheet	polyethylene, green
absorbent core	kraft pulp**
adhesive	hotmelt adhesive
Intermediate layer	tissue, white*
topsheet	polypropylene-nonwoven, hydrophilic, white

*except 40x60 articles

**bleached elementary chlorine free

MoliCare® Premium Bed Mat *(French Tender Assortment)*

Product Characteristics

Characteristics	unit	5 drops 60x60	7 drops 40x60	7 drops 60x60	7 drops 60x90	8 drops 40x60	8 drops 60x60	8 drops 60x90	9 drops 60x90
total weight:*	[g]	38,7	26,3	47,5	81,0	32,5	60,7	94,2	170,8
total length**	[cm]	60	60	60	60	60	60	60	60
total width**	[cm]	60	40	60	90	40	60	90	90

* ≤ ± 10 % tolerance

** ≤ - 3 % tolerance on the surface area

Labelling

Lot-No. with 9-Digit Code

e.g.: 900102822



9

Year 201X

001

Internal key

02

Internal key


82

Internal key

2

cumulated
check sum

Manufacturing Date


e.g.: 

2020
year

02
month

07
day

Expiry Date

e.g.: 

2025
year

02
month

01
day

Medical Device



Keep dry



Do not reuse



CE-marking



Distributor



Distributor

Latest Date of Revision: 2019-11-18