

Foliodrape® General Surgery Sets

General Product Description/Intended Purpose

Foliodrape® General surgery sets sterile are single-use products for short-term use and consist of a variety of products like patient drapes, equipment covers, fixation and collection accessories, commodity products (e.g. towels); included into sterile packages. The Foliodrape® standard sets are intended to be used in different fields of applications/disciplines. Foliodrape® products are to prevent the passage of pathogens between non-sterile and sterile areas. The Polyethylene-film or different layers of hydrophilic nonwoven material laminated with Polyethylene-film act together as a fluid and bacteria barrier and minimize the transmission of micro-organisms. Foliodrape® General surgery sets are placed onto the market sterile and are in Medical Device Class I.

Application/Indication

The Foliodrape® standard sets are used in different fields of applications/disciplines, such as general surgery, gynecology, orthopedics, accident surgery, ear, nose and throat (Otorhinolaryngology or ENT)/oral and maxillofacial surgery, ophthalmology, urology, and cardiothoracic/cardiovascular and neurosurgery. Moreover, universal sets can be used for any type of surgery. Foliodrape® products are for short term use* and prevention of infection on humans in clinical environment. Depending on the procedure for which they are used, the Foliodrape® products usage time will vary from few minutes to several hours of continuous usage.

* Definitions according to Reg. (EU) 2017/745 (MDR)

Reference Numbers

Name of Product	Reference Number
Foliodrape® Protect Plus Abdominal / Spinal set	938 734/3
Foliodrape® Protect Plus Abdomino-Perineal set I, crepe	938 911/3
Foliodrape® Protect Plus Abdomino-Perineal set II	936 807/0
Foliodrape® Protect Plus Abdomino-Perineal set III	936 808/0
Foliodrape® Protect Plus Abdomino-Perineal set IV, crepe	936 805/0
Foliodrape® Protect Plus Abdomino-Perineal set V	936 804/0
Foliodrape® Protect Plus Abdomino-Perineal set VI	938 714/3
Foliodrape® Protect Plus Abdomino-Perineal set VII	936 528/0
Foliodrape® Protect Plus Abdomino-Perineal set VII, crepe	936 279/0
Foliodrape® Protect Basic set I, crepe	938 913/3
Foliodrape® Protect Basic set II, crepe	938 925/3
Foliodrape® Protect Basic set III	938 728/2
Foliodrape® Protect Fenestrated drape set I, crepe	936 761/0
Foliodrape® Protect Fenestrated drape set II	936 762/0
Foliodrape® Protect Plus Laparoscopy set I	938 709/4
Foliodrape® Protect Plus Laparoscopy set II, leggings	938 713/3
Foliodrape® Protect Minor surgery set I	938 722/3
Foliodrape® Protect Minor surgery set I, crepe	938 914/2
Foliodrape® Protect Universal split sheet set I	936 735/0
Foliodrape® Protect Plus Universal split sheet set I, crepe	938 932/3
Foliodrape® Protect Plus Universal split sheet set II	938 708/3
Foliodrape® Protect Plus Universal split sheet set III	938 739/3
Foliodrape® Protect Universal / Paediatrics set I	938 754/2
Foliodrape® Protect Plus Universal / Paediatrics set I, crepe	936 876/0

Foliodrape® General Surgery Sets

Foliodrape® Protect Plus Universal / Paediatrics set II	938 743/4
Foliodrape® Protect Universal set I	938 721/3
Foliodrape® Protect Universal set II, crepe	938 900/3
Foliodrape® Protect Universal set III, crepe	277 011/5
Foliodrape® Protect Universal set IV	938 706/4
Foliodrape® Protect Plus Universal set I, crepe	938 930/2
Foliodrape® Protect Plus Universal set II, crepe	938 901/3
Foliodrape® Protect Plus Universal set III	938 712/3
Foliodrape® Protect Plus Universal set IV	938 707/3
Foliodrape® Protect Plus Universal set V	936 800/0
Foliodrape® Protect Plus Universal set VI	938 744/3
Foliodrape® Protect Plus Universal set VII, extra reinforced, crepe	938 931/2

Residual Risks, Contra-Indications and Undesirable Side Effects, Warnings

No contraindications known.

Preparatory Treatment

Sterile surgical drapes are used during surgery to prevent contact with unprepared surfaces and to maintain the sterility of environmental surfaces, equipment and the patient's surroundings.

Skin and mucosis disinfection prior drape application: Prepare the surface with antiseptic solution, the solution should remain in contact with the skin as per the instruction on the disinfectant label. Do not make contact with any other area of the patient. Lay the drapes on the patient as per their intended purpose recommended by the product name. Leave uncovered only the operative field and those areas necessary for the maintenance of anesthesia.

Sterile Device

Do not use when sterile packaging or the product are damaged or were unintentionally opened before use.

Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of Foliodrape® products should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. All infected waste is recommended to be disposed of by incineration.

Non-contaminated material such as waste paper, packaging and non-sterile but not biologically contaminated materials can be disposed separately. Separate disposal containers should be available where waste is created, so that staff can sort the waste as it is being discarded.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Foliodrape® General Surgery Sets

Product Performance Characteristics

The Foliodrape® range is developed under high-quality standards to ensure safety and infection prevention. All products are free of chlorinated substances such as CMR/PBT/vPvB substances (carcinogenic, mutagenic and toxic to reproduction, etc.). There are no added flame-retardants (like polybrominated bi-phenyls, etc.).


Material Characteristics

For Performance Indicators, visit individual Technical Data Sheets for components listed below.

Protect	Foliodrape® Protect drapes
Protect Plus	Foliodrape® Protect Plus drapes
Drapes With Viscose	Foliodrape® Protect Plus Extra Reinforced
Back table cover	Foliodrape® Back table cover reinforced
Instrument table cover	Foliodrape® Instrument table covers
Adhesive bag	Foliodrape® Adhesive Bags
Cellulose Towels	Foliodrape® Cellulose Towels
OR tape adhesive	Foliodrape® OR Tape
Hook and loop fastener	Foliodrape® Hook and loop fastener
Leggings	Foliodrape® Leggings
Pouches/Bags on drapes	Foliodrape® Collection Bags


Foliodrape® PROTECT and PROTECT PLUS drapes correspond to all requirements of European standard EN 13795-1:2019 for the highest performance level "High performance, critical area". Sets are sterilized with ethylene oxide gas according EN 556-1 and EN 11135-1.

Product Characteristics

Reference Number	Contents Picture	Content Description	Count
938 734/3		Back table cover (Wrapping) 190 x 150 cm	1
		Instrument table cover 80 x 145 cm	1
		Cellulose Towels 33 x 30 cm	4
		OR tape adhesive 10 x 50 cm	1
		Protect Plus Abdominal / Spinal drape 260 x 320 cm 9 x 22 cm fenestration, aperature adhesive, hook and loop fastener	1


Foliodrape® General Surgery Sets

938 911/3




Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	2
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Hook and loop fastener	1
2 x 23 cm	
Protect Plus Abdomino-perineal drape	1
250/280 x 295 cm	
26 x 30 cm fenestrated incision foil, 12 x 15 cm fenestration with 2,5 cm wide adhesive, bilateral suction bags, and leggings	
Protect Surgical drape	1
90 x 75 cm	

936 807/0



Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Surgical drape, adhesive	1
90 x 75 cm	
86 cm adhesive	
Protect Plus Abdomino-perineal drape	1
260/280 x 305 cm	
32 x 28 cm fenestrated incision foil, 12 x 15 cm fenestration, with bilateral suction bags, with leggings	


936 808/0



Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Surgical drape, adhesive	1
90 x 75 cm	
86 cm adhesive	
Protect Plus Anesthetic drape, adhesive	1
300 x 170 cm	
80 cm adhesive	
Protect Plus Abdomino-perineal split drape	1
200 x 200 cm	
32 x 104 cm split, 15 x 20 cm fenestration, with bilateral suction bags. with leggings	

Foliodrape® General Surgery Sets

936 805/0



Back table cover (Wrapping) 1
190 x 150 cm

Cellulose Towels 2
33 x 30 cm


OR tape adhesive 1
10 x 50 cm

Hook and loop fastener 1
2 x 23 cm

Protect Plus Surgical drape, adhesive 1
100 x 90 cm
96 cm adhesive

Protect Plus Gynecology / Abdomino-Perineal drape 1
260/280 x 320 cm
30 x 40 cm fenestrated incision foil, 12 x 15 cm fenestration, with bilateral suction bags, with leggings

936 804/0



Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 2
33 x 30 cm

OR tape adhesive 1
10 x 50 cm


Hook and loop fastener 1
2 x 23 cm

Protect Surgical drape, adhesive 1
90 x 75 cm
86 cm adhesive

Protect Plus Gynecology / Abdomino-Perineal drape 1
260/280 x 320 cm
30 x 40 cm fenestrated incision foil, 12 x 15 cm fenestration, with bilateral suction bags, with leggings

Foliodrape® General Surgery Sets

938 714/3



Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 4
33 x 30 cm


OR tape adhesive 1
10 x 50 cm

Hook and loop fastener 2
2 x 23 cm

Protect Surgical drape 1
90 x 75 cm

Protect Plus Abdomino-Perineal drape 1
250/260 x 290 cm
32 x 28 fenestrated incision foil, 12 x 15 cm fenestration, with leggings

936 528/0



Back table cover (Wrapping) 1
190 x 150 cm

OR tape adhesive 2
10 x 50 cm

Back table cover (Folded) 1
190 x 150 cm

Protect Plus Surgical drape, adhesive 2
100 x 90 cm
96 cm adhesive

Protect Surgical drape, adhesive 1
75 x 75 cm
71 cm adhesive


Protect Plus Surgical drape, adhesive 1
250 x 170 cm
90 cm adhesive

Protect Patient underlay 1
100 x 120 cm

Protect Legging 2
75 x 140 cm

Foliodrape® General Surgery Sets

936 279/0



Back table cover (Wrapping) 1

190 x 150 cm

OR tape adhesive 2

10 x 50 cm

Back table cover (Folded) 1

190 x 150 cm

Protect Plus Surgical drape, adhesive 2

100 x 90 cm

96 cm adhesive

Protect Surgical drape, adhesive 1

75 x 75 cm

71 cm adhesive

Protect Plus Surgical drape, adhesive 1

250 x 170 cm

90 cm adhesive


Protect Patient underlay 1

100 x 120 cm

Protect Legging 2

75 x 140 cm

938 913/3



Back table cover (Wrapping) 1

190 x 150 cm

Cellulose Towels 2


33 x 30 cm

Protect Surgical drape, adhesive 4

75 x 75 cm

71 cm adhesive

938 925/3



Back table cover reinforced (Wrapping) 1

190 x 150 cm

Cellulose Towels 2

33 x 30 cm

OR tape adhesive 1


10 x 50 cm

Protect Surgical drape, adhesive 4

90 x 75 cm

86 cm adhesive

938 728/2



Protect Back table cover (Wrapping) 1

100 x 150 cm

Protect Surgical drape, adhesive 1

175 x 170 cm

90 cm adhesive

Protect Surgical drape, adhesive 1

175 x 170 cm

90 cm adhesive


Protect Surgical drape, adhesive 2

90 x 75 cm

86 cm adhesive

Foliodrape® General Surgery Sets

936 761/0




Back table cover (Wrapping) 1
190 x 150 cm

Cellulose Towels 2
33 x 30 cm

Protect Fenestrated drape, adhesive 1
170 x 250 cm
15 x 12 cm fenestration with aperture adhesive

936 762/0




Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 2
33 x 30 cm

Protect Fenestrated drape, adhesive 1
170 x 250 cm
15 x 12 cm fenestration with aperture adhesive

938 709/4



Back table cover (Wrapping) 1
190 x 150 cm


Instrument table cover 1
80 x 145 cm

Cellulose Towels 4
33 x 30 cm

OR tape adhesive 1
10 x 50 cm

Protect Plus Abdominal drape 1
250/180 x 320 cm
32 x 28 cm fenestrated incision foil, with bilateral suction bags

938 713/3



Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 4
33 x 30 cm


OR tape adhesive 1
10 x 50 cm

Hook and loop fastener 1
2 x 23 cm

Protect Plus Abdominal drape 1
250/260 x 280 cm
32 x 28 cm fenestrated incision foil, with bilateral suction bags, with leggings

Foliodrape® General Surgery Sets

938 722/3




Back table cover (Wrapping) 1
190 x 150 cm

Cellulose Towels 2
33 x 30 cm

Protect Fenestrated drape, adhesive 1
200 x 240 cm
Ø14,5 cm fenestration with aperture adhesive

938 914/2




Back table cover (Wrapping) 1
190 x 150 cm

Cellulose Towels 1
33 x 30 cm

Protect Fenestrated drape, adhesive 1
200 x 240 cm
Ø14,5 cm fenestration with aperture adhesive

936 735/0



Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 4
33 x 30 cm


OR tape adhesive 1
10 x 50 cm

Protect Split sheet, adhesive 1
200 x 260 cm
20 x 115 cm split with adhesive

Protect Surgical drape, adhesive 1
75 x 90 cm
71 cm adhesive

Protect Surgical drape, adhesive 1
240 x 150 cm
90 cm adhesive

938 932/3



Back table cover (Wrapping) 1
190 x 150 cm

Cellulose Towels 2
33 x 30 cm


OR tape adhesive 2
10 x 50 cm

Protect Plus Split sheet, adhesive 1
200 x 260 cm
6,5 x 95 cm split with adhesive

Protect Plus Surgical drape, adhesive 1
250 x 170 cm
90 cm adhesive


Foliodrape® General Surgery Sets

938 708/3




Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Plus Split sheet, adhesive	1
245 x 150 cm	
6,5 x 65 cm split with adhesive	
Protect Plus Surgical drape, adhesive	1
200 x 245 cm	
90 cm adhesive	

938 739/3



Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Plus Split sheet, adhesive	1
200 x 260 cm	
15 x 95 cm split with adhesive	
Protect Plus Surgical drape, adhesive	1
250 x 170 cm	
90 cm adhesive	


938 754/2



Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Surgical drape, adhesive	1
175 x 200 cm	
90 cm adhesive	


Foliodrape® General Surgery Sets

936 876/0




Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Plus Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Plus Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Plus Surgical drape, adhesive	1
175 x 200 cm	
90 cm adhesive	

938 743/4



Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Plus Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Plus Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Plus Surgical drape, adhesive	1
175 x 200 cm	
90 cm adhesive	


938 721/3



Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Surgical drape, adhesive	1
175 x 170 cm	
90 cm adhesive	


Foliodrape® General Surgery Sets

938 900/3




Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	2
10 x 50 cm	
Protect Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Surgical drape, adhesive	1
175 x 170 cm	
90 cm adhesive	

277 011/5



Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	2
33 x 30 cm	
Hook and loop fastener	1
2 x 23 cm	
Protect Surgical drape, adhesive	2
100 x 75 cm	
100 cm adhesive	
Protect Surgical drape, adhesive	1
250 x 170 cm	
90 cm adhesive	
Protect Surgical drape, adhesive	1
200 x 170 cm	
75 cm adhesive	


938 706/4



Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Surgical drape, adhesive	1
175 x 200 cm	
90 cm adhesive	


Foliodrape® General Surgery Sets

938 930/2




Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	2
10 x 50 cm	
Protect Plus Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Plus Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Plus Surgical drape, adhesive	1
175 x 170 cm	
90 cm adhesive	

938 901/3



Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	2
10 x 50 cm	
Protect Plus Surgical drape, adhesive	2
100 x 90 cm	
96 cm adhesive	
Protect Plus Surgical drape, adhesive	1
250 x 170 cm	
90 cm adhesive	
Protect Plus Surgical drape, adhesive	1
200 x 170 cm	
90 cm adhesive	


938 712/3



Back table cover (Wrapping)	1
190 x 150 cm	
Instrument table cover	1
80 x 145 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Plus Surgical drape, adhesive	2
90 x 75 cm	
86 cm adhesive	
Protect Plus Surgical drape, adhesive	1
240 x 150 cm	
90 cm adhesive	
Protect Plus Surgical drape, adhesive	1
200 x 200 cm	
90 cm adhesive	

Foliodrape® General Surgery Sets

938 707/3



Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 4
33 x 30 cm


OR tape adhesive 1
10 x 50 cm

Protect Plus Surgical drape, adhesive 2
100 x 90 cm
96 cm adhesive

Protect Plus Surgical drape, adhesive 1
300 x 170 cm
100 cm adhesive

Protect Plus Surgical drape, adhesive 1
200 x 200 cm
90 cm adhesive

936 800/0



Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 4
33 x 30 cm


OR tape adhesive 1
10 x 50 cm

Protect Plus Surgical drape, adhesive 2
100 x 90 cm
96 cm adhesive

Protect Plus Surgical drape, adhesive 1
300 x 170 cm
100 cm adhesive

Protect Plus Surgical drape, adhesive 1
200 x 200 cm
90 cm adhesive

938 744/3



Back table cover (Wrapping) 1
190 x 150 cm

Instrument table cover 1
80 x 145 cm

Cellulose Towels 4
33 x 30 cm

OR tape adhesive 1
10 x 50 cm


Protect Plus Surgical drape, adhesive 2
120 x 100 cm
116 cm adhesive

Protect Plus Surgical drape, adhesive 1
300 x 170 cm
100 cm adhesive

Protect Plus Surgical drape, adhesive 1
200 x 200 cm
90 cm adhesive

Foliodrape® General Surgery Sets

938 931/2



Back table cover (Wrapping)	1
190 x 150 cm	
Cellulose Towels	4
33 x 30 cm	
OR tape adhesive	1
10 x 50 cm	
Protect Plus Surgical drape, adhesive, extra reinforced	2
120 x 100 cm	
116 cm adhesive	
With Viscose	
Protect Plus Surgical drape, adhesive, extra reinforced	1
300 x 170 cm	
98 cm adhesive	
With Viscose	
Protect Plus Surgical drape, adhesive, extra reinforced	1
200 x 200 cm	
98 cm adhesive	
With Viscose	
Two-section adhesive bag	1
43 x 38 cm	

Labelling

The Foliodrape® labeling describes the intended purpose of the device as supported by clinical evaluation. Safety information supplied by the legal manufacturer on the labeling are based on the risk management documentation and elaborately present all warnings and precautions according the new requirements of the standards concerning the product.

Lot-No. with 8-Digit Code

e.g.:

LOT	0	12	XXXXX
	year	Week of production	for internal purposes only

Manufacturing Date

e.g.:		2020	04	07
		year	month	day

Use-by-Date

e.g.:		2025	04	07
		year	month	day

Shelf Life: 5 years

Medical Device **MD**

Unique Device Identification (UDI)


Single Sterile Barrier System  Single sterile barrier system

Foliodrape® General Surgery Sets

Do not use if the product sterile barrier system or its packaging is compromised



Single use device



Keep dry



Sterilized by Ethylene Oxide



Legal manufacturer

PAUL HARTMANN AG
Paul-Hartmann-Straße 12
89522 HEIDENHEIM, GERMANY

European Conformity



Dispose of your waste in an appropriate manner



Recycling symbols – Corrugated cardboard / Polyethylene Low-Density



Producer contributes to recycling of packaging



Global trade identification number



GTIN

Indication arrow where to open the packaging




Name and intended purpose of the product

Universal-Set III, Kreppeinschlag
Universal set III, crepe

Reference / catalogue number



Pictogram with product look and dimensions



Storage and transport conditions

Transport and store in conditions for climatic regions I-IVa according to WHO Technical Report 953, 2009.

Foliodrape® General Surgery Sets

Logistic Data

Name of product	REF number	Peel Pouch size [mm]	Products per Peel Pouch	Peel Pouches per dispenser	Dispenser size	EAN of dispenser	Dispensers per carton
Foliodrape® Protect Plus Abdominal / Spinal set	938 734/3	Size: XL 345 x 510	1	5	Size: XL 374 x 545 x 234	4052199522616	1
Foliodrape® Protect Plus Abdomino-Perineal set I, crepe	938 911/3	Size: XL 345 x 510	1	4	Size: XL 374 x 545 x 234	4049500032872	1
Foliodrape® Protect Plus Abdomino-Perineal set II	936 807/0	Size: XL 345 x 510	1	3	Size: XL 374 x 545 x 234	4052199521251	1
Foliodrape® Protect Plus Abdomino-Perineal set III	936 808/0	Size: XL 345 x 510	1	4	Size: XL 374 x 545 x 234	4052199523835	1
Foliodrape® Protect Plus Abdomino-Perineal set IV, crepe	936 805/0	Size: XL 345 x 510	1	3	Size: XL 374 x 545 x 234	4052199521336	1
Foliodrape® Protect Plus Abdomino-Perineal set V	936 804/0	Size: XL 345 x 510	1	3	Size: XL 374 x 545 x 234	4052199525600	1
Foliodrape® Protect Plus Abdomino-Perineal set VI	938 714/3	Size: XL 345 x 510	1	5	Size: XL 374 x 545 x 234	4052199213330	1
Foliodrape® Protect Plus Abdomino-Perineal set VII	936 528/0	Size: L 360 x 580	1	5	Size: XL 368 x 568 x 300	4052199517582	1
Foliodrape® Protect Plus Abdomino-Perineal set VII, crepe	936 279/0	Size: L 360 x 580	1	5	Size: XL 368 x 568 x 300	4052199517520	1
Foliodrape® Protect Basic set I, crepe	938 913/3	Size: M 266,5 x 360	1	10	Size: L 274 x 368 x 296	4049500043083	2
Foliodrape® Protect Basic set II, crepe	938 925/3	Size: M 266,5 x 360	1	8	Size: L 274 x 368 x 296	4052199516219	2
Foliodrape® Protect Basic set III	938 728/2	Size: L 360 x 580	1	10	Size: XL 368 x 568 x 300	4052199530567	1
Foliodrape® Protect Fenestrated drape set I, crepe	936 761/0	Size: L 360 x 580	1	10	Size: XL 368 x 568 x 300	4052199516042	1
Foliodrape® Protect Fenestrated drape set II	936 762/0	Size: L 360 x 580	1	10	Size: XL 368 x 568 x 300	4052199515960	1

Foliodrape® General Surgery Sets

Foliodrape® Protect Plus Laparoscopy set I	938 709/4	Size: XL 345 x 510	1	5	Size: XL 374 x 545 x 234	4052199020181	1
Foliodrape® Protect Plus Laparoscopy set II, leggings	938 713/3	Size: XL 345 x 510	1	5	Size: XL 374 x 545 x 234	4052199213309	1
Foliodrape® Protect Minor surgery set I	938 722/3	Size: XL 345 x 510	1	10	Size: XL 374 x 545 x 234	4052199523231	1
Foliodrape® Protect Minor surgery set I, crepe	938 914/2	Size: XL 345 x 510	1	8	Size: XL 374 x 545 x 234	4052199523606	1
Foliodrape® Protect Universal split sheet set I	936 735/0	Size: L 360 x 580	1	6	Size: XL 368 x 568 x 300	4052199516721	1
Foliodrape® Protect Plus Universal split sheet set I, crepe	938 932/3	Size: L 360 x 580	1	6	Size: XL 368 x 568 x 300	4052199516417	1
Foliodrape® Protect Plus Universal split sheet set II	938 708/3	Size: L 360 x 580	1	7	Size: XL 368 x 568 x 300	4052199516691	1
Foliodrape® Protect Plus Universal split sheet set III	938 739/3	Size: L 360 x 580	1	5	Size: XL 368 x 568 x 300	4052199516172	1
Foliodrape® Protect Universal / Paediatrics set I	938 754/2	Size: L 360 x 580	1	8	Size: XL 368 x 568 x 300	4052199517339	1
Foliodrape® Protect Plus Universal / Paediatrics set I, crepe	936 876/0	Size: L 360 x 580	1	7	Size: XL 368 x 568 x 300	4052199515762	1
Foliodrape® Protect Plus Universal / Paediatrics set II	938 743/4	Size: L 360 x 580	1	7	Size: XL 368 x 568 x 300	4052199510217	1
Foliodrape® Protect Universal set I	938 721/3	Size: L 360 x 580	1	9	Size: XL 368 x 568 x 300	4052199530482	1
Foliodrape® Protect Universal set II, crepe	938 900/3	Size: L 360 x 580	1	9	Size: XL 368 x 568 x 300	4052199517292	1
Foliodrape® Protect Universal set III, crepe	277 011/5	Size: L 360 x 580	1	9	Size: XL 368 x 568 x 300	4052199516264	1
Foliodrape® Protect Universal set IV	938 706/4	Size: L 360 x 580	1	8	Size: XL 368 x 568 x 300	4052199517377	1
Foliodrape® Protect Plus Universal set I, crepe	938 930/2	Size: L 360 x 580	1	8	Size: XL 368 x 568 x 300	4052199517179	1
Foliodrape® Protect Plus Universal set II, crepe	938 901/3	Size: L 360 x 580	1	7	Size: XL 368 x 568 x 300	4052199517254	1

Foliodrape® General Surgery Sets

Foliodrape® Protect Plus Universal set III	938 712/3	Size: L 360 x 580	1	7	Size: XL 368 x 568 x 300	4052199517773	1
Foliodrape® Protect Plus Universal set IV	938 707/3	Size: L 360 x 580	1	6	Size: XL 368 x 568 x 300	4052199517018	1
Foliodrape® Protect Plus Universal set V	936 800/0	Size: L 360 x 580	1	6	Size: XL 368 x 568 x 300	4052199500300	1
Foliodrape® Protect Plus Universal set VI	938 744/3	Size: L 360 x 580	1	6	Size: XL 368 x 568 x 300	4052199500195	1
Foliodrape® Protect Plus Universal set VII, extra reinforced, crepe	938 931/2	Size: L 360 x 580	1	6	Size: XL 368 x 568 x 300	4052199517131	1

Latest Date of Revision: 2020-09-04