

MoliCare® Premium Form (ADL)

General Product Description/Intended Purpose

Intended Purpose:

MoliCare® Premium Form products are absorbent pads for incontinence care. The products are, in combination with fixation systems, suitable for usage in case of bladder weakness with several severities, urinary and faecal incontinence with several severities.

The products are for single use only and should be regularly changed as needed.

MoliCare® Premium Form products can be used in home care environment and professional healthcare facilities. The application can be performed by lay person and/or health care professional.

Application/Indication

Several severities of urinary and faecal incontinence



Reference Numbers*

Packaging	P36
4 drops	168 031
5 drops	168 032

*the full assortment may not be available in all countries

Single use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request. Reuse of the product could lead to microbiological contamination and/or mild infections of the skin.

Product disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of MoliCare® Premium Form should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

MoliCare® Premium Form (ADL)

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product performance characteristics

Characteristics	unit	normal plus (4 drops)	extra (5 drops)
total weight:*	[g]	57,9	62,4

* $\leq \pm 10\%$ tolerance

Material characteristics

Part	Material
back sheet:	polyethylene, nonwoven-film-laminate, white
absorbent core	kraft pulp* superabsorbent polymer buffer treated non-woven, hydrophilic, white
adhesive:	hotmelt adhesive
top sheet:	polypropylene - nonwoven, hydrophilic, white
cuff:	polypropylene - nonwoven, hydrophobic, white
elastics:	elastics for leg gather and standing cuffs
wetness indicator:	water soluble ink

* bleached elementary chlorine free

Product characteristics

Dimensions	normal plus (4 drops)	extra (5 drops)
total length:*	620	620
total width abdominal side:*	280	280
total width back side*	295	295
total width crotch:*	235	235
absorbent core length:*	560	560
absorbent core width abdominal side:*	240	240
absorbent core width back side:*	240	240
absorbent core width crotch:*	160	160

* $\leq \pm 10\%$ tolerance

Labelling

9-digit LOT-code

e.g.: 900102822



9	001	02	82	2
year 201X	internal key	internal key	internal key	cumulated check sum

MoliCare® Premium Form (ADL)

Manufacturing date

e.g.:  2019 08 07
year month day

Expiry date

e.g.:  2024 08 01
year month day

Shelf life: 5 years

Medical Device



Keep dry



Do not reuse



CE-marking



Latest date of revision: 2020-08-03