

MediSet[®], sterile

Set de badigeon, 4783702

Intended Purpose

Intended Use:

The pre-operation sets are single-use, sterile sets made with various components, intended to prepare the care site or the surgical site before a care or a surgical act.

During a pre-operation procedure, such operations can be done:

- realization of the pre-operative treatment
- realization of deterision or antisepsis
- pre-interventional antisepsis
- final brush before the surgical act

The care is performed by a healthcare professional on all kinds of patient without restriction on age or gender.

The pre-operation set is used less than 60 minutes and there is no accumulated use of the set.

The components might be in contact with skin, mucosal membranes, injured skin, for a transient use.

Target user/group: Healthcare professionals

Medical device class of the set (based on component with the highest classification): II a

The set is supervised by TÜV – 0123

Application/Indication

The indication results from the necessity to prepare the care site or the surgical site before a care or a surgical act.

That could include:

- The realization of the pre-operation treatment
- The realization of deterision or antisepsis
- The pre-interventional antisepsis
- The final brush before the surgical act

Reference Number

4783702

Quantity	Component description	Included substances
5	Pagasling No. 4 egg-size 24x24 cm	See set label
1	Sponge Forcep atraumatic straight 19 cm	See set label

Instuction for use, Contra-Indications and Undesirable Side Effects, Warnings:



Consult instructions for use



Caution

For possible contraindications please follow the instructions/warning hints of the individual set items listed in the instructions for use enclosed within the set.

MediSet[®], sterile

Set de badigeon, 4783702

Sterile Device:



Sterilized using ethylene oxide



Do not use if package is damaged



Do not re-sterilize

Single Use Device:



Do not re-use

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.



MediSet[®], sterile

Set de badigeon, 4783702

year month day

Shelf Life: Maximum 3 years, or the shelf life of the component with the shortest shelf-life


Medical Devices in the set are indicated by following symbol



Unique Device Identification (UDI) for the set:



Double Sterile Barrier System



Packaging: Multivac in Soft Blister

Primary packaging /transport carton	Primary packaging /dispenser	transport carton/pallet	Primary packaging dimension	Dispenser dimension	Transport carton dimension	Pallet dimension
50	n.a.	15	345 x 210 x 35 mm	n.a.	767 x 377 x 344 mm	1200 x 800 x 1928 mm

Latest Date of Revision: 2020-11-12