

## MoliCare Premium Elastic



Spez.-Nr.: D 3.2109

Abteilung: CMO-ID

Datum: 2017-01-03

### 1. General Product Description

Modern all-in-one incontinence brief with a three-layer anatomically shaped absorbent core, quick dry distribution layer, soft anti-leak cuffs, elastic threads in the crotch area, textile back sheet, wetness indicator, partly breathable side panels in the hip area and reclosable hook & loop tapes. Thanks to a special designed absorbent core, MoliCare Premium Elastic offers maximum leakage protection through high absorbent and holding capacity. MoliCare Premium Elastic has a high retention, supports therefore skin health and provides a high wearing comfort.

Classification acc. EEC-Directive for non-sterile medical devices class I

### 2. Application / Indication

Suitable for usage in cases of moderate to severest urinary and faecal incontinence and also for restless or highly dependent patients as well as for mobile patients providing maximum security.

### 3. Presentations

	<b>sizes</b>			
	1 = small	2 = medium	3 = large	4 = extra large
MoliCare Premium Elastic 6 drops	•	•	•	•
MoliCare Premium Elastic 7 drops	•	•	•	•
MoliCare Premium Elastic 8 drops	•	•	•	•
MoliCare Premium Elastic 9 drops	•	•	•	•

## MoliCare Premium Elastic



Spez.-Nr.: D 3.2109

Abteilung: CMO-ID

Datum: 2017-01-03

4. Product Characteristics**material:**

- backsheet: nonwoven-film-laminate, white printed
- absorbent core:  
(3 layers)
  - 1. layer: 100 % kraft pulp, bleached elementary chlorine free
  - 2. layer: kraft pulp with superabsorbent polymer
  - 3. layer: cellulose, curled fibre
- adhesive: hotmelt adhesive
- nonwoven: polypropylene nonwoven, hydrophilic, white
- cuffs: polypropylene nonwoven, hydrophobic, white
- nonwoven side panels: polypropylene nonwoven, hydrophobic, partly breathable, white
- elastics: elastics for leg gathers and standing cuffs
- tape: nonwoven hook & loop tapes

MoliCare Premium Elastic




Spez.-Nr.: D 3.2109

Abteilung: CMO-ID

Datum: 2017-01-03

5. Product Drawing



Dimensions: [mm]		size 1 small	size 2 medium
total length:	A	815 ± 20	815 ± 20
total width abdominal side:	B	320 -20/+40	320 -20/+40
total width crotch:	F	320 ± 10	320 ± 10
absorbent core length:	C	660 ± 30	660 ± 30
absorbent core width abdominal side:	d	200 ± 10	200 ± 10
absorbent core width tape side:	D	245 ± 10	245 ± 10


MoliCare Premium Elastic



Spez.-Nr.: D 3.2109

Abteilung: CMO-ID

Datum: 2017-01-03



Dimensions: [mm]		size 3 large	size 4 extra large
total length:	A	980 ± 20	980 ± 20
total width abdominal side:	B	480 -20/+40	480 -20/+40
total width crotch:	F	320 ± 10	320 ± 10
absorbent core length:	C	780 ± 40	780 ± 40
absorbent core width abdominal side:	d	200 ± 10	240 ± 10
absorbent core width tape side:	D	265 ± 10	350 ± 10

MoliCare Premium Elastic



Spez.-Nr.: D 3.2109

Abteilung: CMO-ID

Datum: 2017-01-03

6. Product Characteristics

**MoliCare Premium Elastic 6 drops**

characteristics:		size 1 small	size 2 medium	size 3 large	size 4 extra large
total weight:	[g]	85*	96*	113*	125*
1. absorption time:**	[s]	max. 80	max. 80	max. 80	max. 80
2. absorption time:**	[s]	max. 110	max. 110	max. 110	max. 110
1. rewetting rate:**	[g]	max. 0,4	max. 0,4	max. 0,4	max. 0,4
2. rewetting rate:**	[g]	max. 5,0	max. 5,0	max. 5,0	max. 5,0

**MoliCare Premium Elastic 7 drops**

characteristics:		size 1 small	size 2 medium	size 3 large	size 4 extra large
total weight:	[g]	95*	100*	120*	148*
1. absorption time:**	[s]	max. 70	max. 70	max. 70	max. 120
2. absorption time:**	[s]	max. 100	max. 100	max. 100	max. 200
1. rewetting rate:**	[g]	max. 0,4	max. 0,4	max. 0,4	max. 0,4
2. rewetting rate:**	[g]	max. 5,0	max. 5,0	max. 5,0	max. 5,0

**MoliCare Premium Elastic 8 drops**

characteristics:		size 1 small	size 2 medium	size 3 large	size 4 extra large
total weight:	[g]	98*	129*	147*	160*
1. absorption time:**	[s]	max. 70	max. 120	max. 120	max. 120
2. absorption time:**	[s]	max. 100	max. 200	max. 200	max. 200
1. rewetting rate:**	[g]	max. 0,4	max. 0,4	max. 0,4	max. 0,4
2. rewetting rate:**	[g]	max. 5,0	max. 5,0	max. 5,0	max. 5,0

MoliCare Premium Elastic



Spez.-Nr.: D 3.2109

Abteilung: CMO-ID

Datum: 2017-01-03

MoliCare Premium Elastic 9 drops

characteristics:		size 1 small	size 2 medium	size 3 large	size 4 extra large
total weight:	[g]	112*	131*	157*	182*
1. absorption time:**	[s]	max. 120	max. 120	max. 120	max. 120
2. absorption time:**	[s]	max. 200	max. 200	max. 200	max. 200
1. rewetting rate:**	[g]	max. 0,4	max. 0,4	max. 0,4	max. 0,4
2. rewetting rate:**	[g]	max. 5,0	max. 5,0	max. 5,0	max. 5,0

\* tolerance  $\pm 10\%$  (equates to  $\pm 3$  s)

\*\* with 100 or 200 ml saline solution (depending on fluff and SAP content)

7. Labelling

9-digit LOT-code with production time arranged in front:


e.g.: 1221046 – **500102822**

Produktionszeitpunkt			LOT-Code					
Tag	Stunde	Minute						
***	**	**	-	X	XXX	XX	XX	X
				Jahr 201X / 202X	Interner Schlüssel	Interner Schlüssel	Interner Schlüssel	Kumulierte Prüfziffer

8. Packaging

MoliCare Premium Elastic packed into polyethylene-bag; into transport carton acc. DIN, sealed with adhesive tapes, packed onto euro pallet.

January 3<sup>rd</sup> 2017



PAUL HARTMANN AG

– Product Development Incontinence Management (CMO-ID) –