

MoliPants® soft Net pants



Spec.-No.: D 9.0398

Department: CMO-SPE

Date: 2013-12-03

1. General Product Description

The fixation pants MoliPants® soft guarantee a secure fitting of anatomically shaped incontinence pads. The dense structure of the material, which makes it particularly resistant to ladders, gives MoliPants® soft an underwear-like appearance with the wetness indicator still being visible. The soft cuffs, the skin friendly, air permeable material as well as the cross elasticised threads ensure a high wearing comfort. The colour coding of the size and the washing instruction on the inner side of the waistband facilitate the daily work of the caregiver. MoliPants® soft keep their stability of fit for a long time: they are washable up to 50 times at 60°C in a laundry net and can be tumble dried at low heat.

Classification for non sterile medical devices class I, acc. EEC-Directive

2. Application / Indication

MoliPants® soft ensure a safe and hygienic fixation of incontinence pads of all types. These unisex fixation pants can be used with MoliForm® and MoliNea® pads and are suitable for both bedridden and mobile patients / residents. Thanks to their high stretchability, the change and the correct positioning of the pad is facilitated. The size identification is given by the colour of the woven-in MoliPants® soft brand name

Available in six sizes – from small to extra extra extra large:

art.-no.	size	waist / hip girth of wearer / patient: [cm]	
947 795 947 790	Small	45 – 80	P5, OK 60 P25, OK 8
947 796 947 791	Medium	60 – 100	P5, OK 40 P25, OK 8
947 797 947 792	Large	80 – 120	P5, OK 40 P25, OK 8
947 798 947 793	X-large	100 – 160	P5, OK 40 P25, OK 8
947 799 947 794	XX-large	140 – 180	P5, OK 40 P25, OK 8
947 789	XXX-large	160 – 200	P25, OK8

MoliPants[®] soft Net pants



Spec.-No.: D 9.0398

Department: CMO-SPE

Date: 2013-12-03

3. Product Characteristics

characteristics:	Small	Medium	Large	X-large	XX-large	XXX-large
colour-code (waist):	yellow	blue	brown	green	grey	orange
total weight: [g]	19,5 ± 6,0	21,0 ± 6,0	22,5 ± 6,0	25,5 ± 6,0	27,0 ± 6,0	39,5 ± 6,0
product dimensions:						
waistband width, unstretched: [cm]	20,0 ± 4,0	22,5 ± 4,0	25,0 ± 4,0	27,0 ± 4,0	28,5 ± 4,0	35,0 ± 6,0
leg width, unstretched: [cm]	14,0 ± 4,0	16,0 ± 4,0	17,0 ± 4,0	17,5 ± 4,0	18,0 ± 4,0	21,0 ± 6,0

Tolerances are typically according to textile products of this type. / All values are subject to minor modifications.

4. Product Labelling

Writing on waistband (inside):

Washing instructions, machine code, size and 4-digit code:

e.g.:

XX	XX
calendar week	year

5. Packaging

MoliPants[®] soft are packed into folding boxes (except for OP 5 packaging). Transport cartons correspond to DIN, are sealed with adhesive tapes and packed onto euro-pallets.

MoliPants[®] soft OP 5 are packed into bags. Transport cartons correspond to DIN, are sealed with adhesive tapes and packed onto euro-pallets.

MoliPants[®] soft Net pants

Spec.-No.: D 9.0398

Department: CMO-SPE

Date: 2013-12-03

Lot-No. with 9-digit or 10-digit code on folding box / bag:

e.g.:

XX	XX	XX	XXX
year of packing	month of packing	day of packing	for internal purposes only

or

XXX	XX	X	XXXX
for internal purposes only	month of production	year of production	for internal purposes only

December 3rd 2013

PAUL HARTMANN AG
- Product Development CMO-ID -

i.A.

