

MoliForm® Premium Soft



Spec.-No.: D 3.2230

Department: SM-PDI

Date: 2011-03-30

1. General Product Description

Incontinence pad with Dry-Plus topsheet covering, fluid distribution layer, ultra absorbent core and wetness indicator; extra large absorbent pad, complete nonwoven covering, excellent fitting anti-leak-cuffs, anatomically shaped with soft elastics in the crotch area; comfortable textile backsheet. MoliForm® Premium Soft offers optimum leakage protection by ultimate fluid absorption and storage capacity. It features inner cuffs for complete security and confidence. Provides maximum fluid retention as well as enhanced skin and wearing comfort. Absorption capacities are colour-coded for ease of product selection; easy to apply with elasticised net panties (MoliPants®).

Classification acc. EEC-Directive for non sterile medical devices class I

2. Application / Indication

Suitable for any grade of incontinence, also for bed-ridden or mobile and active patients. Selected according to the individual user's level of incontinence. For most severe degrees of incontinence and secure overnight care MoliForm® Premium Soft Maxi / Super / Super Plus are ideal with especially high absorbency. Flexible incontinence system for individual requirements of the wearer.

3. Assortment

Three different sizes with 7 individual absorbent capacities:

size:	absorbent capacity:
MoliForm® Premium Soft Light	for light incontinence
MoliForm® Premium Soft Normal	for light incontinence
MoliForm® Premium Soft Plus	for moderate incontinence
MoliForm® Premium Soft Extra	for severe incontinence
MoliForm® Premium Soft Maxi	for severe incontinence
MoliForm® Premium Soft Super	for very severe incontinence
MoliForm® Premium Soft Super Plus	for very severe incontinence

MoliForm[®] Premium Soft

Spec.-No.: D 3.2230

Department: SM-PDI

Date: 2011-03-30

4. Product Characteristics

Material:

- Liquid protection (backsheet): white textile backsheet printed (polypropylene/polyethylene laminate) yellow, green, blue, pink or purple
- Absorbent core:
(3 layers)
 - 1. layer: cellulose, curled fibre
 - 2. layer: kraft pulp^{*)} with superabsorbent polymer
 - 3. layer: layer: 100 % kraft pulp^{*)}

^{*)} bleached elementary chlorine free
- Nonwoven: polypropylene nonwoven, hydrophilic, white
- Cuffs: polypropylene nonwoven, hydrophobic, white
- Adhesive: hotmelt adhesive
- Elastics: elastics for leg gathers and standing cuffs

MoliForm® Premium Soft




Spec.-No.: D 3.2230

Department: SM-PDI

Date: 2011-03-30

5. Product Drawing



Dimensions [mm]		Light	Normal	Plus	Extra	Maxi	Super	Super Plus
Total length:	A	525 ± 30	620 ± 30	620 ± 30	690 ± 30	690 ± 30	690 ± 30	690 ± 30
Total width abdominal side:	b	260 ± 10	280 ± 10	280 ± 10	300 ± 10	300 ± 10	300 ± 10	300 ± 10
Total width rear side:	B	280 ± 10	295 ± 10	295 ± 10	310 ± 10	310 ± 10	310 ± 10	310 ± 10
Total width crotch:	F	210 ± 10	235 ± 10	235 ± 10	245 ± 10	245 ± 10	245 ± 10	245 ± 10
Absorbent core length:	C	474 ± 20	560 ± 20	560 ± 20	630 ± 20	630 ± 20	630 ± 20	630 ± 20
Absorbent core width abdominal side:	d	215 ± 10	240 ± 10	240 ± 10	250 ± 10	250 ± 10	250 ± 10	250 ± 10
Absorbent core width rear side:	D	242 ± 10	255 ± 10	255 ± 10	260 ± 10	260 ± 10	260 ± 10	260 ± 10
Absorbent core width crotch	E	134 -4/+10	160 -4/+10	160 -4/+10	170 -4/+10	170 -4/+10	170 -4/+10	170 -4/+10

MoliForm® Premium Soft



Spec.-No.: D 3.2230

Department: SM-PDI

Date: 2011-03-30

6. Product Characteristics

Characteristics:		Light	Normal	Plus	Extra	Maxi	Super	Super Plus
Total weight ¹⁾	[g]	51	60	64	95	100	107	129
Weight Superabsorbent Polymer ¹⁾	[g]	5,5	6,6	11,4	12,5	16,4	23,5	20,4
Weight of fluff pulp ¹⁾	[g]	35,5	40,0	40,0	67,5	68,5	68,5	93,9
1 st Absorption time ²⁾	[s]	max. 40	max. 90	max. 90	max. 90	max. 140	max. 140	max. 140
2 nd Absorption time ²⁾	[s]	max. 80	max. 140	max. 140	max. 140	max. 190	max. 190	max. 190
1 st Rewet ²⁾	[g]	max. 0,5	max. 0,5	max. 0,5	max. 0,5	max. 0,5	max. 0,5	max. 0,5
2 nd Rewet ²⁾	[g]	max. 5,0	max. 5,0	max. 5,0	max. 5,0	max. 5,0	max. 5,0	max. 5,0

¹⁾ ± 10% tolerance (equals to ± 3s)

²⁾ with 50 ml for Light; with 100 ml saline solution for Normal/Plus/Extra; with 200 ml saline solution for Super/Super Plus

7. Labelling

Lot-No. with 9-digit code:

e.g.:

0
year

162
day of production

XXXXX
for internal purposes only

8. Packaging

MoliForm® Premium Soft packed into polyethylene bags; transport carton acc. DIN, sealed with adhesive tapes, stacked on euro-pallet.

March 30th 2011

PAUL HARTMANN AG
- Product Development SM-PDI -

i.A. 