

MoliMed® for men activ



Spec.-No.: D 3.2423a

Department: SM-PDI

Date: 2012-03-12

CHANGE CONTROL REPORT

Page /Edition	replaces Page/Edition	Reason for Change
1-3 from 2003-03-13	-----	New Technical Data Sheet
1-3 from 2003-09-09	1-3 from 2003-03-13	Revised version (updated)
1-3 from 2005-09-01	1-3 from 2003-09-09	Revised version (updated) / Relaunch 2005 with white cellulose
1-3 from 2006-07-28	1-3 from 2005-09-01	Revised version (updated) Presentation of an expiry date
1-3 from 2006-07-28	1-3 from 2006-07-28	Expiry date deleted
1-3 from 2012-03-12	1-3 from 2006-07-28	Revised product weight

Department	Date	Signature
Product Development SM-PDI	2012-03-12	<i>J. Eibes</i>

MoliMed® for men activ



Spec.-No.: D 3.2423a

Department: SM-PDI

Date: 2012-03-12

1. General Product Description

Incontinence pad specially designed for men's anatomy; with Dry Plus Distribution Layer and superabsorbent core; neutralizes odour; complete nonwoven covering; air active backsheet with textile comfort; provides optimum security by highest absorption and storing capacity and guarantees protection against re-wetting; particularly comfortable and discreet to wear; easily fixed in underwear with adhesive strip.

Classification acc. EEC-Directive for non sterile medical devices class I

2. Application / Indication

Suitable for active male patients in cases of slight urinary incontinence;

3. Presentations

MoliMed® for men activ

4. Product Characteristics**Material:****MoliMed® for men activ**

- | | |
|------------------------|---|
| - liquid protection: | Nonwoven-film-laminate, air active |
| - absorbent core: | 100 % kraft pulp, bleached elementary chlorine free, kraft pulp with superabsorbent polymer |
| - distributor edition: | Curled Fiber, white |
| - adhesive: | Hotmelt adhesive |
| - nonwoven: | Polypropylene - nonwoven, hydrophilic, white |
| - backing: | Silicon paper |

MoliMed® for men activ




Spec.-No.: D 3.2423a

Department: SM-PDI

Date: 2012-03-12

5. Product Drawing



Dimensions [mm]		MoliMed® for men activ
total length:	A	147 ± 10
total width of outstanding side:	B	115 ± 10
total width of closed side:	b	75 ± 10

MoliMed® for men activ



Spec.-No.: D 3.2423a

Department: SM-PDI

Date: 2012-03-12

6. Product Requirements

Characteristics:		MoliMed® for men activ
weight:	[g]	13,9 min. 9,9
weight of superabsorbent polymer:	[g]	2,7 min. 2,1
weight of pulp:	[g]	6,6 min 5,4
centrifugal test:	[g]	100 min. 40

7. Labelling

Lot-No. with 9-digit code:

e.g.:

2
year


066
day of production

XXXXX
for internal purposes only

8. Packaging

MoliMed® for men activ packed into polyethylene-bag; into transport carton acc. DIN, packed onto europallet.

March 12th 2012



- Product Development SM-PDI -