

MoliMed® Comfort Mini P 28
 MoliMed® Comfort Midi P 28
 MoliMed® Comfort Maxi P 28



Spec.-No.: D 3.2425c
 Department: SM-PDI
 Date: 2012-02-22

1. General Product Description

MoliMed® Comfort is an anatomically shaped pad, with Dry Plus distribution layer, thin absorbent pad with Superabsorber polymer, complete nonwoven covering, waterproof polyethylene backsheet and soft elastics in the crotch area.

MoliMed® Comfort pads provide optimum security by highest absorption and storing capacity and guarantee protection against re-wetting. MoliMed® Comfort can be inserted securely and discretely in the underwear by an adhesive strip.

Classification: non sterile medical devices class I, acc. EEC-Directive

2. Application / Indication

Suitable for use in cases of slight urinary incontinence, as protection following vaginal surgery or as a comfortable, leakproof maternity pad.

3. Presentations

MoliMed® Comfort Mini P 28	art.-no.: 168 100, 168 103
MoliMed® Comfort Midi P 28	art.-no.: 168 101, 168 104
MoliMed® Comfort Maxi P28	art.-no.: 168 102, 168 105

4. Product Characteristics


material:	
- liquid protection:	polyethylene film, white
- absorbent core: (two layered)	100% kraft pulp, bleached elementary chlorine free, kraft pulp with superabsorbent polymer
- distribution-layer:	cellulose, curled fibre, white
- adhesive:	hotmelt adhesive
- nonwoven:	polypropylene nonwoven, hydrophilic
- elastics:	polyurethane foam
- backing:	silicon paper

MoliMed® Comfort Mini P 28
 MoliMed® Comfort Midi P 28
 MoliMed® Comfort Maxi P 28



Spec.-No.: D 3.2425c
 Department: SM-PDI
 Date: 2012-02-22

5. Product Drawing:



dimensions [mm]		Mini	Midi	Maxi
total length:	A	265 ± 6	335 ± 8	430 ± 10
total width:	E	112 ± 6	132 ± 6	162 ± 6
absorbent core length:	B	235 ± 10	305 ± 10	395 ± 10
absorbent core width:	F	90 ± 8	112 ± 8	138 ± 8

MoliMed® Comfort Mini P 28
 MoliMed® Comfort Midi P 28
 MoliMed® Comfort Maxi P 28



Spec.-No.: D 3.2425c

Department: SM-PDI

Date: 2012-02-22

6. Product Characteristics

Characteristics:	Dim.	Mini	Midi	Maxi
total weight:	[g]	14*	20*	36*
superabsorbent polymer weight:	[g]	2,6*	4,0*	7,3*
cellulose weight	[g]	9,3*	13,6*	24,9*
1. absorption time**	[s]	max. 15	max. 45	max. 70
2. absorption time**	[s]	max. 40	max. 80	max. 160
1. rewet**	[g]	max. 0,5	max. 0,5	max. 0,5
2. rewet**	[g]	max. 1,5	max. 4	max. 10

* tolerance $\pm 10\%$ (equates to ± 3 s)

** with 25ml, 50ml or 100ml 0,9% saline solution (depending on cellulose and SAP content)

7. Labelling

Lot-No. with 9-digit code:

e.g.:

1
Year

066
day of the year

XXXXX
for internal purposes only

8. Packaging

MoliMed® Comfort are packed into polyethylene bags; transport cartons acc. DIN, packed onto europallets.

February 22nd 2012

PAUL HARTMANN AG
 - Product Development SM-PDI -

i.A.

